



Hoes, Universal Curettes Langer



supra/subgingival use

supra/subgingival use



GOLDMANN-
FOX
DI-11-101
GF 5

Premolars

Molars

buccal/lingual

175 mm, 7"



GOLDMANN-
FOX
DI-11-102
GF 6

Anteriors e

Premolars

Molars

labial/lingual

mesial/distal

175 mm, 7"



LANGER
DI-11-103
L 1/2

Premolars

Molars

universal

Lower jaw

170 mm, 6³/₄"



LANGER
DI-11-104
L 3/4

Premolars

Molars

universal

Upper jaw

170 mm, 6³/₄"



LANGER
DI-11-105
L 5/6
Anteriors

universal

175 mm, 7"



Universal Curettes



supra/subgingival use



GOLDMANN-
FOX
DI-11-106
GF 2
Anteriors
Premolars

universal

175 mm, 7"



GOLDMANN-
FOX
DI-11-107
GF 4
Premolars
Molars

universal

175 mm, 7"



COLUMBIA
DI-11-108
2R/2L
Anteriors

universal

175 mm, 7"



COLUMBIA
DI-11-109
4R/4L
Premolars
Molars

universal

170 mm, 6³/₄"



COLUMBIA
DI-11-110
13/14
Premolars
Molars

universal

175 mm, 7"



Universal Curettes, GRACEY Curettes



supra/subgingival use

subgingival use



McCALL
for infected tissue
DI-11-111
M 17S/18S

Premolars
Molars

165 mm, 6¹/₂"



DI-11-112
M 23 A
Anteriors

Premolars
Molars
universal

175 mm, 7"



GRACEY
DI-11-113
Fig. 1/2
Anteriors

universal

175 mm, 7"



GRACEY
DI-11-114
Fig. 3/4
Anteriors

Premolars

universal

170 mm, 6³/₄"



GRACEY Curettes



subgingival use



GRACEY
DI-11-115
Fig. 5/6

*Anteriors
Premolars*

universal

175 mm, 7"



GRACEY
DI-11-116
Fig. 7/8

*Premolars
Molars*

buccal/lingual

170 mm, 6³/₄"



GRACEY
DI-11-117
Fig. 9/10

*Premolars
Molars*

buccal/lingual

170 mm, 6³/₄"



GRACEY
DI-11-118
Fig. 11/12

*Premolars
Molars*

*mesio buccal/
mesio lingual*

175 mm, 7"



GRACEY
DI-11-119
Fig. 13/14

*Premolars
Molars*

*disto buccal/
disto lingual*

175 mm, 7"



After-GRACEY Curettes



*subgingival
use*

Shaft 3 mm longer than Standard-GRACEY narrow working tip
– for easier access to the deep root surfaces



After-GRACEY
DI-11-120
Fig. 1/2
Anteriors

universal

170 mm, 6³/₄"



After-GRACEY
DI-11-121
Fig. 3/4
Anteriors

Premolars

universal

170 mm, 6³/₄"



After-GRACEY
DI-11-122
Fig. 5/6
Anteriors

Premolars

universal

170 mm, 6³/₄"



After-GRACEY
DI-11-123
Fig. 7/8
Anteriors

Premolars

Molars

buccal/lingual

170 mm, 6³/₄"



After-GRACEY
DI-11-124
Fig. 11/12
Anteriors

Premolars

Molars

*mesio buccal/
mesio lingual*

170 mm, 6³/₄"



After-GRACEY
DI-11-125
Fig. 13/14
Anteriors

Premolars

Molars

*disto buccal/
disto lingual*

170 mm, 6³/₄"



After-GRACEY Currettes



subgingival
use

Shaft 3 mm longer than Standard-GRACEY, narrow working tip shortened
- for easier access to the deep, narrow root surfaces



Mini-GRACEY
DI-11-126
Fig. 1/2

Anteriors

universal

175 mm, 7"



Mini-GRACEY
DI-11-127
Fig. 3/4

Anteriors
Premolars

universal

170 mm, 6³/₄"



Mini-GRACEY
DI-11-128
Fig. 5/6

Anteriors
Premolars

universal

175 mm, 7"



Mini-GRACEY
DI-11-129
Fig. 7/8

Premolars
Molars
buccal/lingual

170 mm, 6³/₄"



Mini-GRACEY
DI-11-130
Fig. 11/12

Premolars
Molars
mesio buccal/
mesio lingual

175 mm, 7"



Mini-GRACEY
DI-11-131
Fig. 13/14

Premolars
Molars
disto buccal/
disto lingual

175 mm, 7"



Mini-GRACEY Curettes



subgingival
use

Shaft 3 mm longer than Standard-GRACEY, narrow working tip shortened
– for easier access to the deep, narrow root surfaces



Mini-GRACEY
DI-11-132
Fig. 1/2

Anteriors

universal

175 mm, 7"



Mini-GRACEY
DI-11-133
Fig. 3/4

Anteriors
Premolars

universal

170 mm, 6³/₄"



Mini-GRACEY
DI-11-134
Fig. 5/6

Anteriors
Premolars

universal

175 mm, 7"



Mini-GRACEY
DI-11-135
Fig. 7/8

Premolars
Molars
buccal/lingual

170 mm, 6³/₄"



Mini-GRACEY
DI-11-136
Fig. 11/12

Premolars
Molars
mesio buccal/
mesio lingual

175 mm, 7"



Mini-GRACEY
DI-11-137
Fig. 13/14

Premolars
Molars
disto buccal/
disto lingual


175 mm, 7"



Curettes With Interchangeable Working Tips



Fast and clear identification by colour-coded titanium handles



COLUMBIA
DI-11-138

for molars,
purple

length 165 mm, 6 $\frac{1}{2}$ "

McCALL
DI-11-139

for anteriors
and canines,
yellow

165 mm, 6 $\frac{1}{2}$ "

GRACEY
DI-11-140

for anteriors
and canines,
yellow

165 mm, 6 $\frac{1}{2}$ "

GRACEY
DI-11-141

for bicuspid,
green

160 mm, 6 $\frac{1}{4}$ "

GRACEY
DI-11-142

for molars,
purple

165 mm, 6 $\frac{1}{2}$ "

GRACEY
DI-11-143

for molars,
blue

165 mm, 6 $\frac{1}{2}$ "



Curettes, Subgingival



1/1

17,5 cm / 6 3/4"
Universal curette
8 mm, hollow handle



1/1
Columbia
DI-11-144
Fig. 13/14



"Incisors and lateral teeth;
universal"



1/1
Columbia
DI-11-145
Fig. 2R/2L



Incisors, universal



1/1
Columbia
DI-11-146
Fig. 4R/4L




Lateral teeth, universal



1/1
DI-11-147
Fig. GX4



"Incisors and lateral teeth;
universal"




1/1
Mc Call
DI-11-148
Fig. 17S/18S



universal
bifurcate root surfaces



Curettes, Subgingival




17,5 cm / 6 3/4"
Universal curette
8 mm, hollow handle



Langer

DI-11-149

Fig. 1/2



Langer

DI-11-150

Fig. 3/4



Langer

DI-11-151

Fig. 5/6



"Mandible, lateral teeth;
mesial/distal"




"Maxilla, lateral teeth;
mesial/distal"



"Maxilla and mandible,
incisors; universal"



GRACEY-Mini Currettes, Subgingival



1/4

Gracey-MINI
17,5 cm / 6 3/4"
8 mm, hollow handle



DI-11-152
Fig. 1/2

Incisors, universal
First shaft longer,
working end shorter



DI-11-153
Fig. 3/4

"Incisors, premolars;
universal"
First shaft longer,
working end shorter



DI-11-154
Fig. 5/6

"Incisors, premolars;
universal"
First shaft longer,
working end shorter



DI-11-155
Fig. 7/8

"Premolars, molars;
lingual/buccal"
First shaft longer,
working end shorter





GRACEY-Mini Currettes, Subgingival



DI-11-156
Fig. 9/10

"Premolars, molars;
lingual/buccal"
First shaft longer,
working end shorter




DI-11-157
Fig. 11/12

"Premolars, molars;
mesial, mesiolingual,
mesiobuccal"
First shaft longer,
working end shorter




DI-11-158
Fig. 13/14

"Premolars, molars;
distal"
First shaft longer,
working end shorter



DI-11-159
Fig. 15/16

"Premolars, molars;
mesial"
First shaft longer,
working end shorter




DI-11-160
Fig. 17/18

"Premolars, molars;
distal"
First shaft longer,
working end shorter

Gracey-MINI
17,5 cm / 6 3/4"
8 mm, hollow handle



GRACEY-Standard Currettes, Subgingival



DI-11-161
Fig. 1/2



"Incisivos; universal"
Incisives, universel
Incisivi, universali



DI-11-162
Fig. 3/4



"Incisors, premolars;
universal"



DI-11-163
Fig. 5/6



"Incisors, premolars;
universal"



DI-11-164
Fig. 7/8




"Premolars, molars;
lingual/buccal"

Gracey-STANDARD
17,5 cm / 6 3/4"
8 mm, hollow handle



GRACEY-Rigid Curettes, Subgingival



DI-11-165
Fig. 9/10



"Premolars, molars;
lingual/buccal"



DI-11-166
Fig. 11/12



"Premolars, molars;
mesial, mesiolingual,
mesiobuccal"



DI-11-167
Fig. 13/14




"Premolars, molars;
distal"




DI-11-168
Fig. 15/16



"Premolars, molars;
distal"



DI-11-169
Fig. 17/18




"Premolars, molars;
distal"

Gracey-STANDARD
17,5 cm / 6 3/4"
8 mm, hollow handle



GRACEY-Rigid Curettes, Subgingival



DI-11-170
Fig. 1/2



Incisors, universal
extra rigid, for stubborn
dental plaque



DI-11-171
Fig. 3/4



"Incisors, premolars;
universal"
extra rigid, for stubborn
dental plaque



DI-11-172
Fig. 5/6



"Incisors, premolars;
universal"
extra rigid, for stubborn
dental plaque



DI-11-173
Fig. 7/8




"Premolars, molars;
lingual/buccal"
extra rigid, for stubborn
dental plaque



Gracey-RIGID
17,5 cm / 6 3/4"
8 mm. hollow handle



GRACEY-Rigid Currettes, Subgingival



1/4

DI-11-174
Fig. 9/10



"Premolars, molars;
lingual/buccal"
extra rigid, for stubborn
dental plaque



1/4

DI-11-175
Fig. 11/12



"Premolars, molars;
mesial, mesiolingual,
mesiobuccal"
extra rigid, for stubborn
dental plaque



1/4

DI-11-176
Fig. 13/14



"Premolars, molars;
distal"
extra rigid, for stubborn
dental plaque




1/4

DI-11-177
Fig. 15/16




"Premolars, molars;
mesial"
extra rigid, for stubborn
dental plaque



1/4

DI-11-178
Fig. 17/18



"Premolars, molars;
distal"
extra rigid, for stubborn
dental plaque



1/4

Gracey-RIGID
17,5 cm / 6 3/4"
8 mm, hollow handle



GRACEY-Curettes with 10 mm-handle



1/4

Big - Gracey
17,5 cm / 6 3/4"
10 mm, hollow handle



DI-11-179
Fig. 1/2



Incisors, universal



DI-11-180
Fig. 3/4



"Incisors, premolars;
universal"



DI-11-181
Fig. 5/6



"Incisors, premolars;
universal"




DI-11-182
Fig. 7/8



"Premolars, molars;
lingual/buccal"



GRACEY-Curettes with 10 mm-handle



GR
1/4

DI-11-183
Fig. 9/10

"Premolars, molars;
lingual/buccal"

V
1/4

DI-11-184
Fig. 11/12

"Premolars, molars; mesial,
mesiolingual/mesiobuccal"

B
1/4

DI-11-185
Fig. 13/14

"Premolars, molars;
distal"

V
1/4

DI-11-186
Fig. 15/16

"Premolars, molars;
mesial"

B
1/4

DI-11-187
Fig. 17/18


"Premolars, molars;
distal"

1/4

Big - Gracey
17,5 cm / 6 3/4"
10 mm, hollow handle



Sinus Lift Currettes



1/4

Körner/Westermann




Fig. 1
DI-11-188
19 cm / 7 1/2"



Fig. 2
DI-11-189
20 cm / 8"




Fig. 3
DI-11-190
19 cm / 7 1/2"




1/4

Fig. 4
DI-11-191
17 cm / 6 3/4"




1/4

Fig. 5
DI-11-192
18 cm / 7"







Sinus Lift Currettes



DI-11-193



DI-11-194



DI-11-195