



Serving Humanity Since 1908

WELMED
INSTRUMENTS CO.

ENT Catalogue
Group 20, 21, 22, 23, 24 & 25



Manufacturers & Exporters of Quality Surgical Instruments

**Anaesthesia, Otolaryngology, Rhinology, Tonsillectomy,
Laryngo-Bronchoscopy, Oral Maxillo-Facial Surgery
and Cranio-Maxillo-Facial Surgery Instruments**



Otology 21

WELMED Instruments Company
Defence Road
Near Wapda Grid Colony
Sialkot 51310, Pakistan.

Phone : +92 52 3554782, 3558704
Fax : +92 52 3550296
Cell : +92 333 8680640

Email : info@welmedinst.com
welmedpk@gmail.com
wmm@brain.net.pk
welmed@cyber.net.pk

Website : www.welmedinst.com

About Us

Quality begins at **WELMED** Instruments Company where master craftsmen forge and manufacture each instrument to its exact specifications and measurements. Instruments must undergo a series of quality control inspections in order to be worthy of **WELMED** name and warranty. We offer a wide selection of medical instruments including Surgical and Dental instruments.

We have acquired a unique and unmatched understanding of manufacturing through a large experience, working since 1908. Years of experience & state of the art manufacturing facilities helped us to keep our prices low. Our Surgical & Dental instruments are created by expert craftsmen, providing unfaltering quality you can rely on.

We have launched a fully interactive website where you can get online and track your orders 24 hours a day, 7 days in a week, and 365 days a year. Ask for our free e-catalogue on CD-Rom to save time and enjoy hassle free ordering.

We provide our customers with most reliable customer service as your satisfaction is our top priority. Our customers always enjoy our wonderful customer service and quality tools at the best price.

Lastly we invite you to have a work experience with **WELMED** Instruments Company and we will prove ourselves as one of the most valuable source you have.

You can find us at

WELMED Instruments Company
Defence Road, Near Wapda Grid Colony
Sialkot 51310, Pakistan.

Phone : +92 52 3554782, 3558704
Fax : +92 52 3550296
Cell : +92 333 8680640

Emails :
Info@welmedinst.com
welmedpk@gmail.com
wmm@brain.net.pk
welmed@cyber.net.pk

Website : www.welmedinst.com



History

WELMED Instruments Company was founded in 1995. Our Great Grandfather **Mian Muhammad Fazal Elahi** was the founder of M/s M.F. Elahi & Company (E.P.) Pvt. Ltd, established in 1908. WELMED is a sister concern of M.F. Elahi & Co. It is our pride to be a member of M.F. Elahi Family.

WELMED Instruments Company is located in the heart of beautiful city, Sialkot, Pakistan. We are manufacturers and exporters of quality Surgical, Dental instruments and Scissors of all sorts in various categories. We are closely attached with highly skilled surgeons and other technicians of the surgical industry to make the new instruments development possible.

WELMED is in a continuous process of development and testing of new instruments for their precision and perfect functioning. We are also working at imported forging supply by foreign customers.

Certificates

Welmed Instruments Company has achieved the certificates for ISO-9001, ISO-13485, cGMP (FDA) and CE Marking.

Quality Control

All material use in the manufacturing of our instruments are tested from relevant laboratories. Our continuous quality control is our best warranty so that you built confidence in our products.

Deliveries

Deliveries are quick and in-time.

Mode Of Payment

Terms of Payment is ADVANCE or LETTER OF CREDIT (Irrevocable)

Quality Policy

WELMED Instruments Company aims to produce quality Surgical, Dental instruments and Scissors of all sorts at reasonable prices to become an exemplary manufacturer. Our total emphasize is to achieve customer's satisfaction by clearly understanding their requirements and following a quality management system which is a primary concern of all our work force. We do our best to produce instruments of finest quality.



Warranty

All WELMED Instruments are produced from finest stainless steel available to be free from any defect in material or workmanship. Any product that proves to be defective in material or workmanship will either be repaired or replaced without any charge. This warranty is void if the instrument is abused, altered, modified or is not used in the manner for which it is intended.

Otology 21



WELMED
INSTRUMENTS CO.

Allport Incus Hook	63	Gorney Septum	95
Attic Dissectr	67	Greven	27
Austin Duckbill Knife	100	Guilford-Wright	71
Austin Duckbill	45	Guilford-Wright Double Edge Knife	49
Austin Pick	52,53	Guilford-Wright Drum Elevator	101
Austin Sickle	45	Guilford-Wright Footplate Pick	71
Ball End Elevator	96	Guilford-Wright Stapes Pick	56
Barbara Pick	82	Guilford-Wright	30,48
Bellucci Ear	26	Hajek-Ballenger	97
Bellucci Pick	55,56	Halle Elevator	95
Beyer	25	Hartmann	33
Billeau Ear Loop	37	Hartman-Herzfeld	28
Boucheron	18	Hartmann	18,19,32
Boies Nasal	99	Haberman	96
Buck Ear Curette	34,35	Henner Elevator	95
Canal Elevatoi	86	Holmgreen	21
Canal Incision Knife	86	Hook	79
Citelli	25	Hough Crura	29
Cleft Palate	100	Hough Curette	38
Converse Skin Hook	64,65	Hough Excavator	80
Converse-Mckenty	98	Hough Fascia Knife	46
Cottle Single Hook	64	Hough Incision	46
Cottle Skin Elevator	99	Hough Pick	54,55
Curette	72,78	Hough Whirly Bird	46,47
Day Ear Hook	61	House Alligator	26,30
Derlacki Mobilizer	75	House Barbara Pick	57
Derlacki Ossecle	31	House Crabtree Dissector	85
Disposable Blades	43	House Crabtree	59
Double-End Curette	40	House Curette	39,40
Double Edge Knife	70	House Detachable Blade	43
Drum Elevators	70	House Elevator	94,101,102
Drum Scraper	75	House Hook	63
Duckbill Ear Forcep	26	House Lancet	50
Duckbill Knife Elevator	66	House Miniature	27,28,29
Dunning Elevator	97	House Obtuse Pick	57
Ear Knife Handle	43	House Oval Window Pick	57,58
Ear Knife Set	42	House Oval	28
Elevator Chisel	76	House Pick	58,59,84
Elevator Type Knife	70	House Sickle	50
Elevator	76	House Stapes	40
Endaural Curette	38	House Strut Pick	58
Excavtor	73,66	House Teflon Coated	102
Farrior	20	House Tympanoplasty	40,49,85
Fascia Knife	75	House Forceps	26
Fenestra Hook	93	House-Billeau Ear Loop	38
Flsch	22,23,24	House-Buck Curette	38,37
Folding Ear Knife	42	House-Dieters	29
Footplate Chisel	81	Incudostapedial Joint	77,89
Footplate Hook	92	Jordan Calal Incision	49
Footplate Pick	87	Jordon Canal	101
Fowler Double-End	40	Joseph Double Hook	65
Fracture Elevator	99	Joseph Single Hook	65
Freer Septum	96	Joseph Periosteal Elevator	98
Friesner Ear Knife	43	Juers Crimper	30
Gillies Zygoma	99	Kleinsasser	24
Goldman Displacer	99	Knife Curette	90
Goldman Single Hook	64	Kos Crimper	30
Goodhill Ear	26	Kos Pick	59



Lancet Knife	77	Storz Otologic	28
Lempert Elevator	95	Struempel Ear	28
Levers Vent Tube	31	Strut Caliper	80
Lillie Ear Hook	62	Strut Pick	83
Louisville Elevator	100	Suction Elevator	96
Louisville Skin Hook	41	Suraci Zygoma Hook-elevator	61
Lucae	33	Tabb Flap	47
Martin	32	Tenzel Double-periosteal	94
Martin-Alloy	32	Toynbee	18
Maumenee Iris Hook	61	Typanlasty Knife	76
Mccollough Elevator	95	Tyrell Iris Hook	61
Mcgee Crimper	30	Whirly Bird Knife	74
Mcgee Footplate Pick	52	Wire Strut Guide	81
Mcgee Malleus	30	Woodson Elevator	96
Mckenty Periosteal	98	Wullstein	23
Micro Bone Curette	69	Wundepreizer, Volle	24
Morgenstein Pick	60	Zollner	21
Myringotomy Knife	44,77	Zygoma Elevator	99
Needle	78,79,91,92		
Noyes Ear Forcep	26		
Oval Window Pick	82,83,88		
Paparella Elevator	101		
Paparella Incudostapedial	45		
Pennington Elevator	97		
Perforating Needle	88		
Perkins Elevator	100		
Pick	66,68,72,86,90,91		
Pierce Elevator	97		
Politzer	18		
Rasp-Caliper	89		
Roger Elevator	97		
Roller Knife	70,92		
Rosen Knife-Curette	46		
Rosen Pick	53		
Royce Ear Knife	44		
Scheer Crimper	30		
Scheer Pick	60		
Schuknecht Crimper	30		
Schuknecht Footplate	62		
Schuknecht Roller	48		
Sempert Elevator	96		
Shambaugh	63,101		
Shapleigh Curette	36		
Shea Fenestra Hook	62		
Shea Pick	56		
Sheath Handle Ear Knife	42		
Sheehy Knife Curette	79,80		
Sheehy Ossicle	31		
Sheehy-House Knife-Curette	39,51		
Sheehy-House	51		
Sickle Kinfe	51,77,89		
Silverstein Sickle	51		
Single Skin Hook	63		
Stapes Curette	78		
Stapes Hook	92		
Stapes Pick	70		
Storz Ear	27,28		
Storz Ear	27		



POLITZER
Ohrtrichter für Kinder

- **BC-21-101** Ø 2.0 mm
- **BC-21-102** Ø 3.0 mm
- **BC-21-103** Ø 4.0 mm
- **BC-21-104** Ø 5.0 mm



HARTMANN
Ohrtrichter

- **BC-21-105** Ø 1.0 mm
- **BC-21-106** Ø 2.0 mm
- **BC-21-107** Ø 3.0 mm
- **BC-21-108** Ø 3.5 mm
- **BC-21-109** Ø 4.0 mm
- **BC-21-110** Ø 5.0 mm
- **BC-21-111** Ø 5.5 mm
- **BC-21-112** Ø 6.0 mm
- **BC-21-113** Ø 6.5 mm
- **BC-21-114** Ø 7.0 mm



TOYNBEE
Ohrtrichter

- **BC-21-115** Ø 2.0 mm
- **BC-21-116** Ø 3.0 mm
- **BC-21-117** Ø 4.0 mm
- **BC-21-118** Ø 5.0 mm








BOUCHERON
Ohrtrichter

- **BC-21-119** Ø 3.0 mm
- **BC-21-120** Ø 4.0 mm
- **BC-21-121** Ø 5.0 mm
- **BC-21-122** Ø 6.0 mm
- **BC-21-123** Ø 7.0 mm
- **BC-21-124** Ø 8.0 mm
- **BC-21-125** Ø 9.0 mm





HARTMANN
Ohrtrichter , oval

-  **BC-21-126** 1.5 x 2.5 mm
-  **BC-21-127** 2.5 x 3.5 mm
-  **BC-21-128** 4.5 x 5.5 mm
-  **BC-21-129** 5.0 x 6.0 mm
-  **BC-21-130** 5.5 x 6.5 mm









HARTMANN
Ohrtrichter zur
Nachbehandlung

-  **BC-21-131** Ø 8.0 mm
-  **BC-21-132** Ø 9.0 mm



HARTMANN
Ohrtrichter, schwarz

-  **BC-21-133** 3.0 mm
-  **BC-21-134** 4.0 mm
-  **BC-21-135** 5.0 mm
-  **BC-21-136** 7.0 mm
-  **BC-21-137** 8.0 mm
-  **BC-21-138** 9.0 mm



FARRIOR
Ohrtrichter

BC-21-139	○	Ø 3.3 mm
BC-21-140	○	Ø 4.3 mm
BC-21-141	○	Ø 4.8 mm
BC-21-142	○	Ø 5.3 mm
BC-21-143	○	Ø 5.8 mm
BC-21-144	○	Ø 6.3 mm
BC-21-145	○	Ø 6.7 mm
BC-21-146	○	Ø 7,2 mm
BC-21-147	○	Ø 7,7 mm



FARRIOR
Ohrtrichter, oval

BC-21-148	○	3.0 x 4.0 mm
BC-21-149	○	4.0 x 5.0 mm
BC-21-150	○	4.5 x 5.5 mm
BC-21-151	○	5.0 x 6.0 mm
BC-21-152	○	5.5 x 6.5 mm
BC-21-153	○	6.0 x 7.0 mm
BC-21-154	○	6.5 x 7.5 mm
BC-21-155	○	7.0 x 8.0 mm
BC-21-156	○	7.5 x 8.5 mm



FARRIOR
Ohrtrichter, oval,
20 schrage

BC-21-157	○	3.0 x 4.0 mm
BC-21-158	○	4.0 x 5.0 mm
BC-21-159	○	4.5 x 5.5 mm
BC-21-160	○	5.0 x 6.0 mm
BC-21-161	○	5.5 x 6.5 mm
BC-21-162	○	6.0 x 70 mm
BC-21-163	○	6.5 x 7.5 mm
BC-21-164	○	7.0 x 80 mm
BC-21-165	○	7.5 x 8.5 mm



FARRIOR
Ohrtrichter, oval,
35 schrage

BC-21-166	○	3.0 x 4.0 mm
BC-21-167	○	4.0 x 5.0 mm
BC-21-168	○	4.5 x 5.5 mm
BC-21-169	○	5.0 x 6.0 mm
BC-21-170	○	5.5 x 6.5 mm
BC-21-171	○	6.0 x 7.0 mm
BC-21-172	○	6.5 x 7.5 mm
BC-21-173	○	7.0 x 8.0 mm
BC-21-174	○	7.5 x 8.5 mm

BC-21-175  Ø 4 mm

BC-21-176  Ø 5 mm

BC-21-177  Ø 6 mm


BC-21-178  Ø 2.0 mm

BC-21-179  Ø 4.0 mm

BC-21-180  Ø 6.0 mm



HOLMGREEN
Endaural-Spekulum,
selbsthaltend



**HOLMGREEN-
RICHNAU**
Ohrspekulum

BC-21-181  4.5 x 5.5 mm

BC-21-182  5.5 x 6.5 mm

BC-21-183  6.5 x 7.5 mm

BC-21-184  7.5 x 8.5 mm

BC-21-185  8.5 x 9.5 mm

BC-21-186  9.5 x 10.5 mm



ZOLLNER
Ohrtrichter,
oval, schwarz eloxiert



FISCH
BC-21-188



FISCH
BC-21-189
3 x 4 Teeth
25 cm



WULLSTEIN
BC-21-190
3 x 3
13 cm



FISCH
BC-21-191
2 x 2 Teeth sharp
13 cm



WULLSTEIN
BC-21-192
3 x 3 Teeth sharp
11 cm



BC-21-193
13 cm 2 x 2 Teeth sharp



FISCH
BC-21-194
3 x 3 Teeth
13 cm



FISCH
BC-21-195
Modell, 3 x 3 Teeth
11 cm



KLEINSASSER
BC-21-196
16 cm



Wundspreizer, volle
BC-21-197
16 cm



Wundspreizer
BC-21-198
16 cm




FISCH
BC-21-199
3 x 3 Teeth
16 cm



FISCH
BC-21-200
3 x 3 Teeth
16 cm



**Beyer
BC-21-201**
2 mm cut



**Beyer
BC-21-202**
1 mm cut



**Citelli
BC-21-203**
3 mm cut



NOYES EAR FORCEPS

BC-21-204
6 mm serrated jaws
Length of shaft 77 mm
overall length 140 mm, 5.5 inches



NOYES EAR FORCEPS

BC-21-205
1 x 2 teeth
7 mm serrated jaws
Length of shaft 77 mm
Overall length 140 mm, 5.5 inches



HOUSE STRUT

BC-21-209
6 mm smooth jaws
Length of shaft 73 mm
Overall length 135 mm, 5.3 inches



DUCKBILL EAR FORCEPS

BC-21-206
Extra flat
8 mm serrated jaw with rounded tip
Length of shaft 77 mm
Overall length 138 mm, 5.4 inches



HOUSE FORCEPS

BC-21-210
6 mm serrated jaws
side opening to right
Length of shaft 73 mm
Overall length 132 mm, 5.2 inches



BELLUCCI EAR FORCEPS

BC-21-207
8 mm serrated jaws
Length of shaft 75 mm
Overall length 135 mm, 5.3 inches



GOODHILL EAR FORCEPS

BC-21-208
5 mm serrated jaws
Length of shaft 76 mm
Overall length 135 mm, 5.3 inches



HOUSE ALLIGATOR AND CRIMPER FORCEPS

BC-21-211
Very flat 8 mm serrated tapered jaws to a dull point
Length of shaft 73 mm
Overall length 130 mm, 5.1 inches



STORZ EAR FORCEPS

BC-21-212
5 mm finely serrated jaws
Length of shaft 74 mm
Overall length 136 mm, 5.4 inches



STORZ EAR FORCEPS

BC-21-215
4 mm finely serrated jaws
Angled up
Length of shaft 72 mm
Overall length 134 mm, 5.3 inches



HOUSE MINIATURE FORCEPS

BC-21-213
Extremely delicate
4 mm finely serrated jaws
Length of shaft 73 mm
Overall length 134 mm, 5.3 inches



GREVEN EAR FORCEPS

BC-21-216
5 mm jaws with 2.5 mm of tips serrated and touching
Opening behind tip
Length of shaft 73 mm
Overall length 135 mm, 5.3 inches



STORZ OTOLOGIC MICRO-ALLIGATOR FORCEPS

BC-21-214
Miniature ear forceps
3.5mm jaws





BC-21-217
2.5 mm x 6 mm flat cups
1mm x 3mm curette opening
Length of shaft 73 mm
Overall length 137 mm, 5.4 inches

**STRUEMPEL EAR
FORCEPS**



BC-21-218
1.5 mm cups
Length of shaft 71 mm
Overall length 130 mm, 5.1 inches

**HARTMAN-HERZFELD
EAR FORCEPS**



BC-21-222
0.9 mm cups
Angled left 15 degrees
Length of shaft 71 mm
Overall length 134 mm, 5.3 inches

**HOUSE OVAL
CUP FORCEPS**



BC-21-219
3 mm cups
Length of shaft 78 mm
Overall length 137 mm, 5.4 inches

**HARTMAN-HERZFELD
EAR FORCEPS**



BC-21-223
Extra fine
0.9 mm x 2 mm cups
Length of shaft 72 mm
Overall length 134 mm, 5.3 inches

STORZ EAR FORCEPS



BC-21-220
2 mm unbiting cups
Length of shaft 78 mm
Overall length 138 mm, 5.4 inches

**HARTMAN-HERZFELD
EAR FORCEPS**



BC-21-224
0.7 mm x 1 mm fenestrated cups
Length of shaft 73 mm
Overall length 134 mm, 5.3 inches

STORZ EAR FORCEPS



BC-21-221
0.9 mm cups
Angled right 15 degrees
Length of shaft 71 mm
Overall length 130 mm, 5.1 inches

**HOUSE OVAL
CUP FORCEPS**



BC-21-225
0.6 mm x 1 mm cups
Length of shaft 73 mm
Overall length 130 mm, 5.1 inches

**HOUSE MINIATURE
EAR FORCEPS**



BC-21-226
0.6 mm x 1 mm cups
Angled left 15 degrees
Length of shaft 72 mm
Overall length 130 mm, 5.1 inches

**HOUSE MINIATURE
EAR FORCEPS**



BC-21-228
0.75 mm cups
Length of shaft 71 mm
Overall length 133 mm, 5.2 inches

STORZ EAR FORCEPS



BC-21-227
0.6 mm x 1 mm cups
Angled right 15 degrees
Length of shaft 72 mm
Overall length 130 mm, 5.1 inches

**HOUSE MINIATURE
EAR FORCEPS**

Micro-Nipper



BC-21-229
3 mm jaws
Angled right 45 degrees
Length of shaft 71mm
Overall length 130 mm, 5.1 inches

**HOUGH CRURA
NIPPER FORCEPS**



BC-21-230
3 mm jaws
Angled right 45 degrees
Length of shaft 71 mm
Overall length 130 mm, 5.1 inches

**HOUGH CRURA
NIPPER FORCEPS**



BC-21-231
2 mm jaws
Length of shaft 75 mm
Overall length 127 mm, 5.0 inches

**HOUSE-DIETERS
MALLEUS NIPPER**



BC-21-232
1.5 mm jaws
Length of shaft 67 mm
Overall length 126 mm, 5.0 inches

**HOUSE-DIETERS
MALLEUS NIPPER**



BC-21-233
Short 3 mm jaws angled slightly downward
Limited closure
Length of shaft 72 mm
Overall length 130mm, 5.1 inches

**McGEE CRIMPER OR
WIRE CLOSURE
FORCEPS**



BC-21-234
3 mm jaws
Designed to crimp incus replacement
prosthesis around the malleus
Length of shaft 74 mm
Overall length 133 mm, 5.2 inches

**McGEE MALLEUS
CRIMPER OR WIRE
CLOSURE FORCEPS**



BC-21-238
7 mm serrated jaws
Length of shaft 72 mm
Overall length 133 mm, 5.2 inches

**GUILFORD-WRIGHT
LOOP SPREADER
FORCEPS**



BC-21-235
6 mm serrated jaws with concave tip to aid
in manipulating strut
Offset to left for improved vision
Length of shaft 75 mm
Overall length 134 mm, 5.2 inches

**JUERS CRIMPER
FORCEPS**



BC-21-239
Very flat 8 mm serrated jaws tapered to
a dull point
Length of shaft 73 mm
Overall length 130 mm, 5.1 inches

**HOUSE ALLIGATOR
AND CRIMPER
FORCEPS**



BC-21-236
6 mm serrated jaws with concave tip to aid
in manipulating strut
Offset to right for improved vision
Length of shaft 75 mm
Overall length 135 mm, 5.3 inches

**JUERS CRIMPER
FORCEPS**



BC-21-240
3 mm jaws
Limited closure
Length of shaft 72 mm
Overall length 132 mm, 5.2 inches

**KOS CRIMPER
FORCEPS**




BC-21-237
6 mm smooth jaws
Length of shaft 73 mm
Overall length 132 mm, 5.2 inches

**SCHUKNECHT WIRE
CRIMPER**



BC-21-241
5 mm smooth jaws have a single length-
wise groove in the lower jaw and a "V"
notch in the tip of both jaws to aid in
positioning wire prosthesis
Length of shaft 72 mm
Overall length 130 mm, 5.1 inches

**SCHEER CRIMPER
FORCEPS**




BC-21-242
 1 mm teeth
 1.5 mm x 4 mm vertical groove
 1.5 mm x 2 mm horizontal groove

DERLACKI OSSICLE HOLDING CLAMP



Thumbscrew



BC-21-243
 5 mm Teeth
 1.5 mm x 4 mm vertical grooves
 8mm from tip to end of serrations

SHEEHY OSSICLE HOLDING CLEMP



Thumbscrew
 Overall length 121mm, 4.8 inches

Micro-Otological Vent Tube Introducers



LEIVERS VENT TUBE INTRODUCER FORCEPS

BC-21-244
 Delicate
 Lower jaw in the form of a 2.5 mm rod to fit into the vent tube
 6 mm movable upper jaw is concave on the underside to securely but gently grasp vent tube
 Length of shaft 73 mm
 Overall length 136 mm, 5.4 inches





Martin
BC-21-245
c 128



Martin
BC-21 -246
C¹ 256



Martin-Alloy
BC-21-247
C² 512



Martin-Alloy
BC-21-248
C³ 1024



Martin-Alloy
WI-21-249
C⁴ 2048



Hartmann
BC-21-250
C 128



Hartmann
BC-21-251
C¹ 256



Hartmann
BC-21-252
C² 512



Hartmann
BC-21-253
C³ 1024



Hartmann
WI-21-254
C⁴ 2048



Hartmann
WI-21-255
C⁵ 4096



Lucae
BC-21-256
C 128



Lucae
BC21-257
C 256



Lucae
BC-21-258
C² 512



Lucae
BC-21-259
C³ 1024



Lucae
BC-21-260
C⁴ 2048



Lucae
BC-21-261
C⁵ 4096




Lucae



BC-21-262
C 128
movable up to h



Hartmann



BC-21-263
C 128
movable up to h



Hartmann
BC21-264
C 256

Blunt



BC-21-265
1 mm diameter
slightly angled tip



Round handle
Overall length 151 mm, 6.0 inches

BUCK EAR CURETTE

Blunt




BC-21-266
1.5 mm diameter
slightly angled tip
Round handle
Overall length 152 mm, 6.0 inches



Round handle
Overall length 151 mm, 6.0 inches

BUCK EAR CURETTE

Blunt



BC-21-267
3 mm diameter
slightly angled tip
Round handle
Overall length 153 mm, 6.0 inches

BUCK EAR CURETTE

Blunt




BC-21-268
2 mm diameter
slightly angled tip
Round handle
Overall length 152 mm, 6.0 inches



Round handle
Overall length 151mm, 6.0 inches


BUCK EAR CURETTE

Blunt



BC-21-269
2.5 mm diameter
slightly angled tip
Round handle
Overall length 153 mm, 6.0 inches

BUCK EAR CURETTE




BC-21-270
1 mm diameter
slightly angled tip



Round handle
Overall length 151 mm, 6.0 inches

BUCK EAR CURETTE




BC-21-271
1.5 mm diameter
slightly angled tip
Round handle
Overall length 152 mm, 6.0 inches




Round handle
Overall length 151 mm, 6.0 inches

BUCK EAR CURETTE



BC-21-272
2 mm diameter
slightly angled tip
Round handle
Overall length 152 mm, 6.0 inches

BUCK EAR CURETTE




BC-21-273
2.5 mm diameter
slightly angled tip
Round handle
Overall length 153 mm, 6.0 inches




Round handle
Overall length 151 mm, 6.0 inches

BUCK EAR CURETTE



BC-21-274
3 mm diameter
slightly angled tip
Round handle
Overall length 153 mm, 6.0 inches

BUCK EAR CURETTE




**SHAPLEIGH EAR
CURETTE**

BC-21-275
Small
Very delicate 1.5 mm serrated loop
Flexible shaft



Flat handle
Overall length 163 mm, 6.4 inches




**SHAPLEIGH EAR
CURETTE**

BC-21-276
Large
Very delicate 2 mm serrated loop
Flexible shaft



Flat handle
Overall length 163 mm, 6.4 inches




**SPRAGUE EAR
CURETTE**

BC-21-277
Small
2 mm curved loop



Round handle
Overall length 174 mm, 6.9 inches




**SPRAGUE EAR
CURETTE**

BC-21-278
Large
3 mm curved loop



Round handle
Overall length 174 mm, 6.9 inches




BC-21-279
Small
2 mm x 3 mm loop
Flexible shaft

BILLEAU EAR LOOP



Round handle
Overall length 161 mm, 6.3 inches




BC-21-280
Medium
3 mm x 5 mm loop
Flexible shaft
Round handle
Overall length 162 mm, 6.4 inches

BILLEAU EAR LOOP




Round handle
Overall length 161 mm, 6.3 inches



BC-21-281
Large
4 mm x 7 mm loop
Flexible shaft
Round handle
Overall length 163 mm, 6.3 inches

BILLEAU LOOP



BC-21-282
No.1
2 mm diameter
Dull curette tip

HOUSE-BUCK CURETTE



Round handle
Overall length 152 mm, 6.0 inches



BC-21-283
No.2
2.5 mm diameter
Dull curette tip

HOUSE- BUCK CURETTE



Round handle
Overall length 153 mm, 6.0 inches




BC-21-284
3.5 mm oval
Teflon coated

**HOUSE-BILLEAU
EAR LOOP**



Overall length 164 mm, 6.5 inches

Otological Curettes, Endaural




BC-21-285
1.5 mm x 3 mm cup

ENDAURAL CURETTE




Overall length 205 mm, 8.2 inches




BC-21-286
2.5 mm x 4 mm cup
Overall length 205 mm, 8.2 inches

ENDAURAL CURETTE



BC-21-287
3 mm x 4.5 mm cup
Overall length 205 mm, 8.2 inches

ENDAURAL CURETTE




BC-21-288
Large
2.25 mm x 3 mm serrated cup

HOUGH CURETTE



Strong angle 35 degrees
Overall length 146 mm, 5.8 inches



BC-21-289
Small
1.5 mm x 2 mm cup
Light angle, 18 degrees

**GUILFORD-WRIGHT
CURETTE**




Flat handle
Overall length 134 mm, 5.3 inches



BC-21-290
Medium
2.25 mm x 3 mm cup
Light angle, 18 degrees

**GUILFORD-WRIGHT
CURETTE**


Flat handle
Overall length 134 mm, 5.3 inches



BC-21-291
Large
2 mm x 2.5 mm cup
strong angle, 35 degrees

**GUILFORD-WRIGHT
CURETTE**

Flat handle
Overall length 134 mm,
5.3 inches



BC-21-292
2 mm diameter blade, semi sharp
straight shaft, malleable
Round handle

**SHEEHY-HOUSE
KNIFE-CURETTE**





Referred to as "Small Weapon"
Overall length 168 mm, 6.6 inches



BC-21-293
2.5 mm diameter blade, semi sharp
straight shaft, malleable
Round handle

**SHEEHY-HOUSE
KNIFE-CURETTE**




BC-21-294
Double-ended
2 mm x 2.5 mm cups and
2.25 mm x 3 mm cups

HOUSE CURETTE



Light angle cups
Flat handle with serrations
Overall length 178 mm, 7.0 inches




BC-21-295
Double-ended
2 mm x 2.5 mm cups and
2.25 mm x 3 mm cups

HOUSE CURETTE



Angled 35 degrees
Flat handle with serrations
Overall length 180 mm, 7.2 inches



BC-21-296
Double-ended
2 mm x 2.5 mm cups and
2.25 mm x 3 mm cups

**HOUSE STAPES
CURETTE**



Angled 35 degrees
Flat handle with serrations
Overall length 152 mm, 6.0 inches




BC-21-297
2 mm cups
One end angled

**FOWLER DOUBLE-END
CURETTE**



Overall length 179 mm, 7.0 inches



BC-21-298
1 mm cup
Angled 35 degrees
Angled shaft, 10 degrees

**HOUSE
TYMPANOPLASTY
CURETTE**




Round hexagonal handle
56 mm from angle to tip
Overall length 168 mm, 6.6 inches

Plastic Surgery Curettes



Overall length 205 mm, 8.0 inches

BC-21-299
DOUBLE-END CURETTE 6 mm x 9 mm and 5 mm x 8 mm cups



BC-21-300
Fine

**LOUISVILLE
SKIN HOOK**



Overall length 162 mm, 6.5 inches




BC-21-301
Medium

**LOUISVILLE
SKIN HOOK**



Overall length 162 mm, 6.3 inches



BC-21-302

**LOUISVILLE
BONE HOOK**



Overall length 152 mm, 6.0 inches




BC-21-303
Blunt

**LOUISVILLE
VEIN HOOK**



Overall length 162 mm, 6.5 inches




BC-21-304
Straight 9 mm downcutting Royce blade
Can be placed in closed handle for protection
Handle length folded 109 mm, 4.3 inches



Blade can be set automatically and rigidly at any point from straight to 40 degree angle without use of screws
Overall length 181 mm, 7.1 inches

**FOLDING
EAR KNIFE**



BC-21-305
Straight 8 mm downcutting Royce blade
Can be placed in closed handle for protection
Rigid Handle when folded and locked




Straight 8 mm downcutting blade
Handle length folded 109 mm, 4.3 inches
Overall length 190 mm, 7.5 inches

**FOLDING
HANDLE EAR KNIFE**



BC-21-306




Straight 9 mm downcutting Royce blade
Can be placed in closed handle for protection
Angular handle
Length of sheath handle 105 mm, 4.2 inches
Overall length 165 mm, 6.5 inches

**SHEATH
HANDLE EAR KNIFE**




BC-21-307



Straight 6 mm Lance blade
Can be placed in closed handle for protection
Angular handle
Length of sheath handle 105 mm, 4.2 inches
Overall length 165 mm, 6.5 inches

**SHEATH
HANDLE EAR KNIFE**



BC-21-308
Consists of 6 disposable bevel blades
and chuck handle



Malleable shafts
Overall assembled length 162 mm, 6.4 inches


**EAR
KNIFE SET**



**EAR
KNIFE HANDLE
BC-21-309**



Overall length 92 mm, 3.6 inches



BC-21-310
3 mm disposable bevel blades
Malleable shaft

DISPOSABLE BLADES



Blades packaged 6 per package
Overall length 82 mm, 3.2 inches




BC-21-311
4 mm disposable bevel blade
Malleable shaft

**DISPOSABLE
EAR KNIFE**



Flat spot on handle designed to guide knife in per
Overall length 165 mm, 6.5 inches




BC-21-312
3 mm disposable Lance blade
Can fit into the NO 521 Storz handle

**HOUSE DETACHABLE
BLADE**



Overall length 82 mm, 3.2 inches




BC-21-313
9 mm blade

**FRIESNER EAR KNIFE
WITH PROTECTOR**



Overall length 152 mm, 6.0 inches




BC-21-314
7 mm
upcutting Royce blade
Bayonet style shaft

ROYCE EAR KNIFE



63 mm from angle to tip
Overall length 185 mm, 7.3, inches



BC-21-315
7 mm downcutting Royce blade
Bayonet style shaft

ROYCE EAR KNIFE




63 mm from angle to tip
Overall length 185 mm, 7.3, inches




BC-21-316

ROYCE EAR KNIFE



8 mm down cutting Royce blade
Angled shaft
60 mm from angle to tip
Overall length 158 mm, 6.2 inches




BC-21-317
8 mm straight blade
Hexagonal handle

MYRINGOTOMY KNIFE



Overall length 166 mm, 6.5 inches



BC-21-318
9 mm slightly curved blade
Hexagonal handle

MYRINGOTOMY KNIFE



Overall length 168 mm, 6.6 inches




BC-21-319
2 mm wide curved tip
Round handle

**AUSTIN DUCKBILL
KNIFE ELEVATOR**



Overall length 160 mm, 6.3 inches




BC-21-320
3 mm curved blade
Round handle

AUSTIN SICKLE KNIFE



Overall length 160 mm, 6.3 inches




BC-21-321
0.5 mm wide x 1.5 mm long
Angled 20 degrees

**PAPARELLA
INCUDOSTAPEDIAL
JOINT KNIFE**



Octagonal handle
Overall length 163 mm, 6.4 inches




BC-21-322
5.5 mm pointed curved blade
Octagonal handle

**PAPARELLA
SICKLE KNIFE**



Overall length 163 mm, 6.4 inches



BC-21-323
Duckbill type sharp
3 mm wide x 8 mm long curved blade

**PAPARELLA
ELEVATOR KNIFE**



Octagonal handle
Overall length 163 mm, 6.4 inches



BC-21-324
Small
2 mm diameter blade
Angled up 45 degrees

ROSEN KNIFE-CURETTE



Round handle
Overall length 168 mm, 6.6 inches




BC-21-325
Large
2.5 mm diameter blade
Angled up 45 degrees

ROSEN KNIFE-CURETTE



Round handle
Overall length 168 mm, 6.6 inches



BC-21-326
2 mm wide x 3 mm long, with rounded
triangular tip
Angled 45 degrees

**HOUGH INCISION
KNIFE**



Hexagonal handle
Overall length 150 mm, 5.9 inches



BC-21-327
5 mm diameter blade
Angled up 45 degrees

HOUGH FASCIA KNIFE



Hexagonal handle
Overall length 150 mm, 5.9 inches



BC-21-328
1 mm wide x 5 mm long, strongly curved tip,
offset to the right

**HOUGH "WHIRLY BIRD"
KNIFE**



Hexagonal handle
Overall length 150 mm, 5.9 inches




BC-21-329
No. 2
1 mm wide x 5 mm long, strongly curved tip,
offset to the left

**HOUGH "WHIRLY BIRD"
KNIFE**



Hexagonal handle
Overall length 150 mm, 5.9 inches



BC-21-330
No.3
0.75 mm wide x 5 mm long, lightly curved tip,
offset to the right

**HOUGH "WHIRLY BIRD"
KNIFE**



Hexagonal handle
Overall length 150 mm, 5.9 inches




BC-21-331
No.4
0.75 mm x 5 mm long, lightly curved tip,
offset to the left

**HOUGH "WHIRLY BIRD"
KNIFE**



Hexagonal handle
Overall length 150 mm, 5.9 inches



BC-21-332
Double-ended
0.9 mm x 1 mm blade
Angled at 25 and 45 degrees

TABB FLAP KNIFE



Round handle
Overall length 160 mm, 6.3 inches




BC-21-333
Double-ended
0.9 mm x 1 mm blade
Angled at 90 and 120 degrees

TABB FLAP KNIFE



Round handle
Overall length 160 mm, 6.3 inches






BC-21-334
3 mm diameter roller blade in line
with handle
Curved shaft

**SCHUKNECHT
ROLLER KNIFE**



Flat handle
15 mm from curve to tip
Overall length 192 mm, 7.6 inches




BC-21-335
3 mm diameter roller blade
angle 90 degrees to handle
Curved shaft

**SCHUKNECHT
ROLLER KNIFE**



Flat handle
15 mm from angle to tip
Overall length 192 mm, 7.6 inches




BC-21-336
3 mm diameter roller blade
15 mm from angle to tip

**GUILFORD-WRIGHT
ROLLER KNIFE**



Round handle
Overall length 140 mm, 5.5 inches




BC-21-337
3 mm wide tip
Flat handle

**GUILFORD-WRIGHT
ELEVATOR KNIFE**



Overall length 140 mm, 5.5 inches



BC-21-338
2 mm wide tip
Flat handle

**GUILFORD-WRIGHT
ELEVATOR KNIFE**



Overall length 140 mm, 5.5 inches




BC-21-339
7 mm strong curved blade
Round handle

**GUILFORD-WRIGHT
DOUBLE EDGE KNIFE**



Overall length 140 mm, 5.5 inches




BC-21-340
1.5 mm long blade
Square handle

**GUILFORD-WRIGHT
INCUDOSTAPEDIAL KNIFE**



Overall length 140 mm, 5.5 inches




BC-21-341
2 mm x 5 mm blade
Round handle

**GUILFORD-WRIGHT
FLAP KNIFE**



Overall length 168 mm, 6.6 inches




BC-21-342
2 mm wide x 3 mm long, tapered handle
Angled 45 degrees

**JORDAN CANAL
INCISION KNIFE**



Round handle
Overall length 166mm, 6.5 inches




BC-21-343
7 mm curved blade
Round handle

**HOUSE
TYMPANOPLASTY KNIFE**



Overall length 168 mm, 6.6 inches




BC-21-344
2 mm x 3 mm tapered blade
Referred to as "Number Two Knife"
Angled 45 degrees

**HOUSE
LANCET KNIFE**



Round handle
Overall length 168 mm, 6.6 inches




BC-21-345
0.5 mm wide x 1.5 mm long
Angled 20 degrees

**HOUSE
INCUDOSTAPEDIAL
JOINT KNIFE**



Round handle
Overall length 168 mm, 6.6 inches



BC-21-346
0.75 mm diameter blade
Angled 60 degrees

**HOUSE
INCUDOSTAPEDIAL
JOINT KNIFE**



Round handle
Overall length 168 mm, 6.6 inches




BC-21-347
7.5 mm curved blade
Referred to as "Number One Knife"

**HOUSE
SICKLE KNIFE**



Round handle
Overall length 168 mm, 6.6 inches



BC-21-348
3.5 mm blade
Round handle

**HOUSE
MYRINGOTOMY KNIFE**



Overall length 168 mm, 6.6 inches




BC-21-349
5 mm curved blade
Round handle

SICKLE KNIFE



Overall length 168 mm, 6.6 inches




BC-21-350
2 mm diameter semi-sharp blade
Referred to as "Small Weapon"
Malleable shaft

**SHEEHY-HOUSE
KNIFE-CURETTE**



Round handle
Overall length 168 mm, 6.6 inches



BC-21-351
2.5 mm diameter semi-sharp blade
Referred to as "Large Weapon"
Malleable shaft

**SHEEHY-HOUSE
KNIFE-CURETTE**



Round handle
Overall length 168 mm, 6.6 inches




BC-21-352
5 mm blade with slightly rounded tip
Round handle

**SILVERSTEIN
SICKLE KNIFE**



Overall length 180 mm, 7.1 inches




BC-21-353
1 mm diameter
Angled 33 degrees

**SILVERSTEIN
ROUND KNIFE**




Round handle
Overall length 180 mm, 7.1 inches




BC-21-354
Anterior
0.2 mm long point
Angled up 90 degrees
shaft angled 20 degrees

MCGEE FOOTPLATE PICK




60 mm from angle to tip
Round handle
Overall length 165 mm, 6.5 inches




BC-21-355
Posterior
0.2 mm long point
Angled down 90 degrees
shaft angled 20 degrees

MCGEE FOOTPLATE PICK



60 mm from angle to tip
Round handle
Overall length 165 mm, 6.5 inches



BC-21-356
Interior
1 mm long point
Angled down 90 degrees
shaft angled 20 degrees

MCGEE FOOTPLATE PICK




60 mm from angle to tip
Round handle
Overall length 165 mm, 6.5 inches




BC-21-357
Posterior
1 mm long point
Angled down 90 degrees
shaft angled 20 degrees

MCGEE FOOTPLATE PICK



60 mm from angle to tip
Round handle
Overall length 165 mm, 6.5 inches



BC-21-358
Slightly curved tip
Round handle

AUSTIN PICK



Overall length 160 mm, 6.3 inches



BC-21-359
1 mm long tip
Angled 45 degrees

AUSTIN PICK



Round handle
Overall length 160mm, 6.3 inches




BC-21-360
3 mm long point
Angled 90 degrees

AUSTIN PICK



Round handle
Overall length 160mm, 6.3 inches




BC-21-361
1 mm long point
Angled 90 degrees

AUSTIN PICK



Round handle
Overall length 160mm, 6.3 inches




BC-21-362
Curved
3 mm curve from tip

ROSEN PICK



Round handle
Overall length 168mm, 6.6 inches




BC-21-363
Curved
1.5 mm curve from tip

ROSEN PICK



Round handle
Overall length 168mm, 6.6 inches




BC-21-364
Straight point
Hexagonal handle

HOUGH PICK



Overall length 150 mm, 5.9 inches




BC-21-365
Slightly curved tip
10 mm curve from tip

HOUGH PICK



Hexagonal handle
Overall length 150 mm, 5.9 inches




BC-21-366
Light
Slightly curved tip
10 mm curve to tip

HOUGH PICK



Hexagonal handle
Overall length 150 mm, 5.9 inches




BC-21-367
0.33 mm long tip
Angled 45 degrees
50 mm from angle to tip

HOUGH PICK



Hexagonal handle, angled 15 degrees
Overall length 150 mm, 5.9 inches



BC-21-368
0.33 mm long tip
Angled 90 degrees
50 mm from angle to tip

HOUGH PICK




Hexagonal handle, angled 15 degrees
Overall length 150 mm, 5.9 inches




BC-21-369
0.66 mm long tip
Angled 90 degrees
50 mm from angle to tips

HOUGH PICK



Hexagonal handle, angled 15 degrees
Overall length 150 mm, 5.9 inches



BC-21-370
Strong curved point
Medium shaft
5 mm curve from tip

BELLUCCI PICK



Round handle
Overall length 168 mm, 6.6 inches




BC-21-371
Short curved point
Fine shaft
3 mm curve from tip

BELLUCCI PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-372
Medium curved point
Fine shaft
7 mm curve from tip

BELLUCCI PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-373
0.5 mm long point
Angled 90 degrees
Fine shaft

BELLUCCI PICK



Round handle
Overall length 168 mm, 6.6 inches




BC-21-374
1.0 mm long point
Angled 90 degrees
Fine shaft

BELLUCCI PICK



Round handle
Overall length 168 mm, 6.6 inches




BC-21-375
Sharp point
3 mm curve from tip

SHEA PICK



Round handle
Overall length 170 mm, 6.7 inches




BC-21-376
Blunt 2 mm long tip
Angled 35 degrees

SHEA PICK



Round handle
Overall length 171 mm, 6.7 inches




BC-21-377
1.25 mm long tip
Angled 40 degrees

**GUILFORD-WRIGHT
STAPES PICK**




Round handle
Overall length 140 mm, 5.5 inches




BC-21-378
0.5 mm long tip
Angled 45 degrees

**GUILFORD-WRIGHT
FOOTPLATE PICK**



Square handle
Overall length 140 mm, 5.5 inches



BC-21-379
Straight point
Malleable shaft

HOUSE BARBARA PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-380
0.25 mm long point
Angled 30 degrees
Malleable shaft

HOUSE OBTUSE PICK



Round handle
Overall length 160 mm, 6.3 inches




BC-21-381
0.33 mm long point
Angled 30 degrees
Malleable shaft

HOUSE OVAL WINDOW PICK



Round handle
Overall length 168 mm, 6.6 inches




BC-21-382
0.66 mm long point
Angled 30 degrees
Malleable shaft

HOUSE OVAL WINDOW PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-383
0.33 mm long point
Angled 90 degrees
Malleable fine shaft

HOUSE OVAL WINDOW PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-384
0.33 mm long point
Angled 90 degrees
Malleable fine shaft

**HOUSE OVAL
WINDOW PICK**



Round handle
Overall length 168 mm, 6.6 inches



BC-21-385
0.5 mm long point
Angled 90 degrees
Malleable medium shaft

HOUSE STRUT PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-386
0.66 mm long point
Angled 90 degrees
medium shaft

HOUSE STRUT PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-387
1 mm long point
Angled 90 degrees
Malleable heavy shaft

HOUSE PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-388
1.5 mm long point
Angled 90 degrees
Malleable heavy shaft

HOUSE PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-389
3mm long tip
Angled 90 degrees
Malleable heavy shaft

HOUSE PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-390
Original design
0.5 mm wide x 2.5 mm long flattened tip
Angled 90 degrees

**HOUSE CRABTREE
DISSECTOR PICK**



Referred to as the "Jimmy"
Round handle
Overall length 168 mm, 6.6 inches



BC-21-391
Small 0.5 mm wide x 1.5 long, flattened
tip
Angled 90 degrees

**HOUSE CRABTREE
DISSECTOR**



Round handle
Overall length 168 mm, 6.6 inches




BC-21-392
Straight point
85 mm from angle to tip

KOS PICK



Bayonet flat handle
Overall length 172 mm, 6.8 inches




BC-21-393
0.33 mm long point
Angled up 90 degrees
85 mm from angle to tip

KOS PICK



Bayonet flat handle
Overall length 172 mm, 6.8 inches



BC-21-394
0.3 mm long point
Angled 90 degrees
Malleable shaft

SCHEER PICK



Hexagonal handle
Overall length 162 mm, 6.4 inches




BC-21-395
1 mm long point
Angled 90 degrees
Malleable shaft

SCHEER PICK



Hexagonal handle
Overall length 162 mm, 6.4 inches



BC-21-396
2 mm blunt tip
Angled 90 degrees

MORGENSTEIN PICK




Round handle
Overall length 164 mm, 6.5 inches



BC-21-397
Flat serrated handle

**TYRELL
IRIS HOOK**

Overall length 122 mm, 4.8 inches



BC-21-398
Flat serrated handle

**MAUMENEE
IRIS HOOK**

Overall length 125 mm, 4.8 inches




ORAL MAXILLOFACIAL

**SURACI ZYGOMA
HOOK-ELEVATOR**



BC-21-399
Overall length 228 mm, 9.0 inches


OTOLOGIC, BASIC



BC-21-400
Extra small
1 mm blunt tip
Angled 90 degrees


DAY EAR HOOK

Round handle
Overall length 159 mm, 6.3 inches




BC-21-401
Small
2 mm blunt tip
Angled 90 degrees
Round handle
Overall length 159 mm, 6.3 inches

DAY EAR HOOK



BC-21-402
Large
2.5 mm blunt tip
Angled 90 degrees
Round handle
Overall length 159 mm, 6.3 inches

DAY EAR HOOK



BC-21-403
1 mm diameter ball tip
Angled 90 degrees
2 mm from angle to tip

LILLIE EAR HOOK



Flat handle
Overall length 173 mm, 6.8 inches

OTOLOGIC, MICRO



BC-21-404
Anterior
0.3 mm long tip
Slight upward curve
shaft angled 18 degrees

**SCHUKNECHT
FOOTPLATE HOOK**




55 mm from angle to tip
Round handle
Overall length 165 mm, 6.5 inches




BC-21-405
Posterior
0.6 mm long tip
Full downward curve
shaft angled 18 degrees

**SCHUKNECHT
FOOTPLATE HOOK**




55 mm from angle to tip
Round handle
Overall length 165 mm, 6.5 inches



BC-21-406
1 mm long point
shaft angled 18 degrees
55 mm from angle to tip

**SCHUKNECHT
STAPES HOOK**



Round handle
Overall length 165 mm, 6.5 inches




BC-21-407
1 mm long point
Angled 45 degrees

**SHEA FENESTRA
HOOK**



Round handle
Overall length 171 mm, 6.7 inches



BC-21-408
1 mm tip
Angled 90 degrees

**SHAMBAUGH
AND DERLACKI
MICRO HOOK**



Flat handle
Overall length 156 mm, 6.1 inches



BC-21-409
1 mm tip
Angled left 90 degrees
Shaft angled 8 degrees

HOUSE HOOK



Round handle, malleable shaft
53 mm from angle to tip
Overall length 168 mm, 6.6 inches




BC-21-410
1 mm long tip
Angled right 90 degrees
Shaft angled 8 degrees

HOUSE HOOK



Round handle, malleable shaft
53 mm from angle to tip
Overall length 168 mm, 6.6 inches




BC-21-411
3 mm tip angled 90 degrees

**ALLPORT
INCUS HOOK**



Overall length 172 mm, 6.8 inches

SOFT TISSUE



BC-21-412
Fine
Short Working end

**SINGLE
SKIN HOOK**



Overall length 125 mm, 4.9 inches




BC-21-413
Shallow curve
Small

COTTLE SINGLE HOOK



Overall length 133 mm, 5.2 inches




BC-21-414
Deep curve
Small

COTTLE SINGLE HOOK



Overall length 133 mm, 5.2 inches




BC-21-415
Deep curve
Large

COTTLE SINGLE HOOK



Overall length 134 mm, 5.3 inches




BC-21-416

**GOLDMAN
SINGLE HOOK**



Overall length 136 mm, 5.4 inches



BC-21-417
Small

**CONVERSE
SKIN HOOK**



Overall length 188 mm, 7.4 inches




BC-21-418
Large

**CONVERSE
SKIN HOOK**



Overall length 186 mm, 7.3 inches



BC-21-419

**JOSEPH
SINGLE HOOK**



Overall length 156 mm, 6.1 inches

**JOSEPH
DOUBLE HOOK
BC-21-420**
2 mm wide



Overall length 156 mm, 6.1 inches

**JOSEPH
DOUBLE HOOK
BC-21-421**
5 mm wide



Overall length 158 mm, 6.2 inches

**JOSEPH
DOUBLE HOOK
BC-21-422**
7 mm wide



Overall length 158 mm, 6.2 inches

**JOSEPH
DOUBLE HOOK
BC-21-423**
10 mm wide



Overall length 158 mm, 6.2 inches



BC-21-424
2 mm wide curved tip
Round handle
Overall length 160 mm, 6.3 inches

DUCKBILL KNIFE ELEVATOR




BC-21-425
0.75 mm wide tip
Angled to left
Round handle
Overall length 160 mm, 6.3 inches

EXCAVATOR



BC-21-426
0.75 mm wide tip
Angled to right
Round handle
Overall length 160 mm, 6.3 inches

EXCAVATOR




BC-21-427
Slightly curved tip
Round handle
Overall length 160 mm, 6.3 inches

PICK




BC-21-428
1 mm long tip
Angled 45 degrees
Round handle
Overall length 160mm, 6.3 inches

PICK




BC-21-429
 3 mm long point
 Angled 90 degrees
 Round handle
 Overall length 160 mm, 6.3 inches

PICK




BC-21-430
 1 mm long point
 Angled 90 degrees
 Round handle
 Overall length 160 mm, 6.3 inches

PICK




BC-21-431
 5 mm flattened tip
 Angled 90 degrees
 16 mm from angled to tip
 Round handle
 Overall length 160 mm,
 6.3 inches

ATTIC DISSECTOR




BC-21-432
 3.5 mm, 4 mm, and 4.5 mm from tip
 6 mm from angle to tip
 Round handle
 Overall length 160 mm, 6.3 inches

STRUT CALIPER




BC-21-433
 3 mm curved blade
 Round handle
 Overall length 160mm, 6.3 inches

SICKLE KNIFE




BC-21-434
Strong curved point
Medium shaft
5 mm curve from tip
Round handle
Overall length 168 mm, 6.6 inches

PICK




BC-21-435
Short curved point
Fine shaft
4 mm curve from tip
Round handle
Overall length 168 mm, 6.6 inches

PICK




BC-21-436
Medium curved point
Fine shaft
7 mm curve from tip
Round handle
Overall length 168 mm, 6.6 inches

PICK




BC-21-437
0.5 mm long point
Angled 90 degrees
Medium shaft
Hexagonal handle
Overall length 162 mm, 6.4 inches

PICK



BC-21-438
0.5 mm long point
Angled 90 degrees
Fine shaft
Round handle
Overall length 168 mm, 6.6 inches

PICK




BC-21-439
1.0 mm long point
Angled 90 degrees
Fine shaft
Round handle
Overall length 168 mm, 6.6 inches

PICK



Small




BC-21-441
1.5 mm x 2 mm cup
Angled 18 degrees
Flat handle
Overall length 134 mm, 5.3 inches

MICRO BONE CURETTE



Large




BC-21-442
2.25 mm x 3 mm cup
Angled 18 degrees
Flat handle
Overall length 134 mm, 5.3 inches

MICRO BONE CURETTE




Medium




BC-21-443
2. mm x 2.5 mm cup
Angled 35 degrees
Flat handle
Overall length 134 mm, 5.3 inches

MICRO BONE CURETTE




BC-21-444
3 mm diameter blade
15 mm from angle to tip
Round handle
Overall length 140 mm, 5.5 inches

ROLLER KNIFE




BC-21-445
3 mm wide tip
Round handle
Overall length 140 mm, 5.5 inches

ELEVATOR TYPE KNIFE




BC-21-446
2 mm wide tip
Round handle
Overall length 140 mm, 5.5 inches

ELEVATOR TYPE KNIFE




BC-21-447
7 mm long curved, pointed blade
Round handle
Overall length 140 mm, 5.5 inches

DOUBLE EDGE KNIFE




BC-21-448
1.25 mm long tip
Angled 40 degrees
Round handle
Overall length 140 mm, 5.5 inches

STAPES PICK



BC-21-449
2 mm wide
15 mm from angle to tip
Round handle
Overall length 140 mm, 5.5 inches

DRUM ELEVATOR




BC-21-450
1.5 mm blade
Square handle

**GUILFORD- WRIGHT
INCUDOSTAPEDIAL
KNIF**



Overall length 140 mm, 5.5 inches




BC-21-451
2 mm x 5mm blade
Round handle

**GUILFORD- WRIGHT
ELEVATOR FLAP KNIFE**



Overall length 168 mm, 6.7 inches




BC-21-452
0.5 mm long tip
Angled 45 degrees

**GUILFORD- WRIGHT
FOOTPLATE PICK**



Square handle
Overall length 140 mm, 5.5 inches



BC-21-453
1 mm wide x 4 mm long measuring
plate with 1 mm markings on both sides
Angled 65 degrees


**GUILFORD- WRIGHT
FENESTROMETER**



Malleable shaft
Round handle
Overall length 140 mm, 5.5 inches




Large




BC-21-454
2.25 mm x 3 mm serrated cup
Angled 35 degrees
Hexagonal handle
Overall length 143 mm, 5.6 inches

CURETTE



BC-21-455
Straight point
Hexagonal handle
Overall length 150 mm, 5.9 inches

PICK




BC-21-456
Slightly curved tip
10 mm curved from tip
Hexagonal handle
Overall length 150 mm, 5.9 inches

PICK




Light



BC-21-457
Slightly curved tip
10 mm curved from tip
Hexagonal handle
Overall length 150 mm, 5.9 inches

PICK




BC-21-458
 0.33 mm long tip
 Angled 45 degrees
 50 mm from angle to tip
 Hexagonal handle, angled 15 degrees
 Overall length 150 mm, 5.9 inches

PICK




BC-21-459
 0.5 mm long flattened tip
 Angled 45 degrees
 Hexagonal handle, angled 15 degrees
 Overall length 150 mm, 5.9 inches

EXCAVATOR




BC-21-460
 Hoe
 0.5 mm long tip
 Angled 90 degrees
 50 mm from angle to tip
 Hexagonal handle, angled 15 degrees
 Overall length 150 mm, 5.9 inches

EXCAVATOR




BC-21-461
 0.33 mm long tip
 Angled 90 degrees
 50 mm from angle to tip
 Hexagonal handle, angled 15 degrees
 Overall length 150 mm, 5.9 inches

PICK




BC-21-462
 0.66 mm long tip
 Angled 90 degrees
 50 mm from angle to tip
 Hexagonal handle, angled 15 degrees
 Overall length 150 mm, 5.9 inches

PICK




BC-21-463
No.1
1 mm wide x 5 mm long strong
curved tip, offset to the right
Hexagonal handle
Overall length 150 mm, 5.9 inches

“WHIRLY BIRD” KNIFE




BC-21-464
No.2
1 mm wide x 5 mm long strong
curved tip, offset to the left
Hexagonal handle
Overall length 150 mm, 5.9 inches

“WHIRLY BIRD” KNIFE




BC-21-465
No.3
0.75 mm wide x 5 mm long lightly
curved tip, offset to the right
Hexagonal handle
Overall length 150 mm, 5.9 inches

“WHIRLY BIRD” KNIFE




BC-21-466
No.4
0.75 mm wide x 5 mm long lightly
curved tip, offset to the left
Hexagonal handle
Overall length 150 mm, 5.9 inches

“WHIRLY BIRD” KNIFE




BC-21-467
 4 prong, 2 mm diameter, curved tip
 Hexagonal handle
 Overall length 150 mm, 5.9 inches

DERLACKI MOBILIZER




BC-21-468
 2 mm wide x 3 mm long, rounded
 triangular tip
 Angled 45 degrees
 Overall length 150 mm, 5.9 inches

INCISION KNIFE




BC-21-469
 5 mm diameter blade
 Angled 45 degrees
 Hexagonal handle
 Overall length 150 mm, 5.9 inches

FASCIA KNIFE




BC-21-470
 0.5 mm long point on 45 degree angle,
 2.5 mm long tip, angled 60 degrees
 50 mm from angle to tip
 Hexagonal handle, angled 15 degrees
 Overall length 150 mm, 5.9 inches

DRUM SCRAPER




BC-21-471
1mm wide x 12 mm long tip
Angled 15 degrees
Referred to as the "Gimmick"
Round handle
Overall length 168 mm, 6.6 inches

ELEVATOR




BC-21-472
1mm wide x 12 mm long tip
Angled 15 degrees
Teflon(R) Coated
With diathermy attachment
Overall length 168 mm, 6.6 inches

COATED ELEVATOR




BC-21-473
1.5 mm wide, curved scraping tip
88 mm working length
Flat bayonet handle
Overall length 170 mm, 6.7 inches

ELEVATOR & CHISEL




BC-21-474
1 mm wide, curved scraping tip
88 mm working length
Flat bayonet handle
Overall length 170 mm, 6.7 inches

ELEVATOR & CHISEL




BC-21-475
7 mm curved blade
Round handle
Overall length 168 mm, 6.6 inches

TYPANOPLASTY KNIFE




BC-21-476
 2 mm x 3 mm tapered tip
 Angled 45 degrees
 Referred to as " Number Two Knife"
 Round handle
 Overall length 168 mm, 6.6 inches

LANCET KNIFE




BC-21-477
 0.75 mm diameter blade
 Angled 60 degrees
 Round handle
 Overall length 168 mm, 6.6 inches

INCUDOSTAPEDIAL JOINT KNIFE




BC-21-478
 0.75 mm curved blade
 Referred to as " Number One Knife"
 Round handle
 Overall length 168 mm, 6.6 inches

SICKLE KNIFE




BC-21-479
 0.5 mm wide x 1.5 mm long blade
 Angled 20 degrees
 Round handle
 Overall length 168 mm, 6.6 inches

INCUDOSTAPEDIAL JOINT KNIFE




BC-21-480
 3.5 mm blade
 Round handle
 Overall length 168 mm, 6.6 inches

MYRINGOTOMY KNIFE




BC-21-481
Light angle
Double-ended
2 mm and 2.5 mm oval cups
Flat handle
Overall length 178 mm, 7.0 inches

CURETTE




BC-21-982
Strong angle
Double-ended
2.25 mm x 3 mm and 2 mm 2.5 mm cups
Flat handle
Overall length 178 mm, 7.0 inches

CURETTE




BC-21-483
Double-ended
2.25 mm x 3 mm and 2 mm x 2.5 mm cups
Flat handle
Overall length 152 mm, 6.0 inches

STAPES CURETTE




BC-21-484
Curved tip
Semi-sharp point
Malleable shaft
Round handle
Overall length 168 mm, 6.6 inches

NEEDLE




BC-21-485
 Rosen type curved and slightly
 dulled tip
 Malleable shaft
 Overall length 168 mm, 6.6 inches

NEEDLE




BC-21-486
 1 mm long tip
 Angled left 90 degrees
 Malleable shaft, angled 8 degrees
 53 mm from angle to tip
 Overall length 168 mm, 6.6 inches

HOOK




BC-21-487
 1 mm long point
 Angled right 90 degrees
 Malleable shaft, angled 8 degrees
 53 mm from angle to tip
 Overall length 168 mm, 6.6 inches

HOOK




BC-21-488
 2 mm diameter semi-sharp blade
 Referred to as "Small Weapon"
 Malleable shaft
 Overall length 168 mm, 6.6 inches

SHEEHY KNIFE CURETTE



BC-21-489
2.5 mm diameter semi-sharp blade
Referred to as "Large Weapon"
Malleable shaft
Overall length 168 mm, 6.6 inches

SHEEHY KNIFE CURETTE




BC-21-490
Rosen point (semi-sharp)
3mm curve from tip
Malleable shaft
Overall length 160 mm, 6.3 inches

NEEDLE




BC-21-491
0.33 mm long tip
Angled 90 degrees
Malleable shaft
Overall length 168 mm, 6.6 inches

HOUGH EXCAVATOR




BC-21-492
2 mm diameter
Angled 45 degrees
Overall length 168 mm, 6.6 inches

MIDDLE EAR MIRROR




BC-21-493
Marker 4.5 mm from tip
Malleable shaft
Overall length 168 mm, 6.6 inches

STRUT CALIPER



BC-21-494
Marker 3.5 mm from tip
Malleable shaft
Overall length 168 mm, 6.6 inches

STRUT CALIPER




BC-21-495
Marker 4 mm from tip
Malleables length 168 mm, 6.6 inches
Overall length 168mm, 6.6 inches

STRUT CALIPER




Straight




BC-21-496
0.6 mm tip
Malleable shaft, angled 10 degrees
62 mm from angle to tip
Overall length 156 mm, 6.2 inches

FOOTPLATE CHISEL




BC-21-497
0.6 mm long tip with notched end
Malleable shaft, angled 10 degrees
62 mm from angle to tip
Overall length 156 mm, 6.2 inches

WIRE STRUT GUIDE




BC-21-498
Marker 5 mm from tip
Malleable shaft
Overall length 168 mm, 6.6 inches

STRUT CALIPER




BC-21-499
Straight point
Malleable shaft
Overall length 168mm, 6.6 inches

BARBARA PICK




BC-21-500
0.25 mm long point
Angled 30 degrees
Malleable shaft
Overall length 160mm, 6.3 inches

OBTUSE PICK




BC-21-501
0.33 mm long point
Angled 30 degrees
Malleable shaft
Overall length 168 mm, 6.6 inches

OVAL WINDOW PICK




BC-21-502
0.66 mm long point
Angled 30 degrees
Malleable shaft
Overall length 168 mm, 6.6 inches

OVAL WINDOW PICK




BC-21-503
0.33 mm long point
Angled 90 degrees
Malleable shaft
Overall length 168 mm, 6.6 inches

OVAL WINDOW PICK




BC-21-504
0.33 mm long point
Angled 90 degrees
Malleable medium shaft
Overall length 168 mm, 6.6 inches

OVAL WINDOW PICK




BC-21-505
0.5 mm long point
Angled 90 degrees
Malleable medium shaft
Overall length 168 mm, 6.6 inches

STRUT PICK



BC-21-506
0.66 mm long tip
Angled 90 degrees
Malleable medium shaft
Overall length 168 mm, 6.6 inches

STRUT PICK




BC-21-507
1 mm long tip
Angled 90 degrees
Malleable heavy shaft

HOUSE PICK



Round handle
Overall length 168 mm, 6.6 inches




BC-21-508
1.5 mm long point
Angled 90 degrees
Malleable heavy shaft
Overall length 168 mm, 6.6 inches

HOUSE PICK



Round handle
Overall length 168 mm, 6.6 inches



BC-21-509
3 mm long tip
Angled 90 degrees
Malleable heavy shaft
Overall length 168 mm, 6.6 inches


HOUSE PICK



Round handle
Overall length 168 mm, 6.6 inches



Referred to as the "Jimmy"
Round handle
Overall length 168 mm, 6.6 inches




BC-21-510
Original design
0.5 mm wide x 2.5 mm long, flattened tip
Angled 90 degrees

**HOUSE CRABTREE
DISSECTOR PICK**



Referred to as the "Jimmy"
Round handle
Overall length 168 mm, 6.6 inches



BC-21-511
Small 0.5 mm wide x 1.5 mm long, flattened tip
Angled 90 degrees
Overall length 168 mm, 6.6 inches

**HOUSE CRABTREE
DISSECTOR PICK**




56 mm from angle to tip
Round handle
Overall length 168 mm, 6.6 inches




BC-21-512
Small
1 mm cup, angled 35 degrees
Shaft angled 10 degrees

**HOUSE
TYMPANOPLASTY
CURETTE**




BC-21-513
2 mm wide x 3 mm long
tapered blade
Angled 45 degrees
Overall length 166 mm, 6.5 inches

CANAL INCISION KNIFE




BC-21-514
2 mm diameter blade
Angled 60 degrees
Overall length 166 mm, 6.5 inches

CANAL ELEVATOR



BC-21-515
Straight point
85 mm, from angle to tip
Bayonet flat handle
Overall length 172 mm, 6.8 inches

PICK




BC-21-516
0.33 mm long point
Angled up 90 degrees
85 mm, from angle to tip
Bayonet flat handle
Overall length 172 mm, 6.8 inches

PICK



Anterior




BC-21-520
 0.2 mm long point
 Angled up 90 degrees
 shaft angled 20 degrees
 60 mm from angle to tip
 Overall length 165 mm, 6.5 inches

FOOTPLATE PICK



Posterior




BC-21-521
 0.2 mm long point
 Angled down 90 degrees
 Shaft angled 20 degrees
 60 mm from angle to tip
 Overall length 165 mm, 6.5 inches

FOOTPLATE PICK



Anterior




BC-21-522
 1 mm long point
 Angled down 90 degrees
 Shaft angled 20 degrees
 60 mm from angle to tip
 Overall length 165 mm, 6.5 inches

FOOTPLATE PICK



Posterior




BC-21-523
 1 mm long point
 Angled down 90 degrees
 Shaft angled 20 degrees
 60 mm from angle to tip
 Overall length 165 mm, 6.5 inches

FOOTPLATE PICK



Small




BC-21-524
0.4 mm diameter rasp
Shaft angled 20 degrees
60 mm from angle to tip
Overall length 165 mm, 6.5 inches

OVAL WINDOW RASP




Large




BC-21-525
0.7 mm diameter rasp
Shaft angled 20 degrees
60 mm from angle to tip
Overall length 165 mm, 6.5 inches

OVAL WINDOW RASP




BC-21-526
Straight point
Shaft angled 20 degrees
60 mm from angle to tip
Overall length 165 mm, 6.5 inches

PERFORATING NEEDLE




BC-21-527
 0.5 wide x 1.5 mm long blade
 Angled 20 degrees
 Octagonal handle
 Overall length 163 mm, 6.4 inches

INCUDOSTAPEDIAL JOINT KNIFE




BC-21-528
 Mushroom type rasp
 Markers at 3.5 mm and 4.5mm
 Malleable handle
 Overall length 163 mm, 6.4 inches

RASP-CALIPER




BC-21-529
 1 mm wide x 4 mm long blade,
 with 1 mm, markings on one side and
 0.5 mm markings on the other side
 Malleable shaft
 Octagonal handle
 Overall length 163 mm, 6.4 inches

FENESTROMETER




BC-21-530
 5.5 mm curved, pointed blade
 Octagonal handle
 Overall length 163 mm, 6.4 inches

SICKLE KNIFE




BC-21-531
 3mm wide x 8 mm long, curved
 duckbill blade
 Octagonal handle
 Overall length 163 mm, 6.4 inches

ELEVATOR KNIFE



BC-21-532
Curved
3mm curve from tip
Overall length 168 mm, 6.6 inches

PICK




BC-21-533
Curved
1.5mm curve from tip
Overall length 168 mm, 6.6 inches

PICK



Small




BC-21-534
2mm diameter blade, angled up
45 degrees
Overall length 168 mm, 6.6 inches

PICK




Large




BC-21-535
2.5mm diameter blade, angled up
45 degrees
Overall length 168 mm, 6.6 inches

KNIFE CURETTE




BC-21-536
Malleable shaft
Hexagonal handle
Overall length 162 mm, 6.4 inches

NEEDLE



BC-21-537
Slightly curved
6 mm curve from tip
Malleable shaft
Hexagonal handle
Overall length 162 mm, 6.4 inches

NEEDLE




BC-21-538
0.3 mm long point
Angled 90 degrees
Malleable shaft
Hexagonal handle
Overall length 162 mm, 6.4 inches

PICK




BC-21-539
1 mm long point
Angled 90 degrees
Malleable light shaft
Hexagonal handle
Overall length 162 mm, 6.4 inches

PICK



BC-21-540
3 mm diameter blade in line
with handle
Curved shaft
15 mm curve from tip
Overall length 192 mm, 7.6 inches

ROLLER KNIFE




BC-21-541
3 mm diameter roller blade,
angled 90 degrees to handle
Curved shaft
15 mm curve from tip
Overall length 192 mm, 7.6 inches

ROLLER KNIFE



Anterior




BC-21-542
0.3 mm long tip
Slight upward curve
shaft angled 18 degrees
55 mm from angle to tip
Overall length 165 mm, 6.5 inches

FOOTPLATE HOOK




Posterior




BC-21-543
0.6 mm long tip
Full downward curve
shaft angled 18 degrees
55 mm from angle to tip
Overall length 165 mm, 6.5 inches

FOOTPLATE HOOK




BC-21-544
1 mm long point
shaft angled 18 degrees
55 mm from angle to tip
Overall length 165 mm, 6.5 inches

STAPES HOOK




BC-21-545
shaft angled 18 degrees
55 mm from angle to tip
Overall length 165 mm, 6.5 inches

NEEDLE




BC-21-546
Sharp point
3 mm curve from tip
Overall length 171 mm, 6.7 inches

PICK




BC-21-547
1 mm long point
Angled 45 degrees
Overall length 171 mm, 6.7 inches

FENESTRA HOOK




BC-21-548
Blunt 2 mm long tip
Angled 35 degrees
Overall length 171 mm, 6.7 inches

PICK



BC-21-549
4 mm x 6mm dissector end with shaft
marked at 10 mm, 20 mm, 30,mm from tip
6 mm, wide curved elevator



TENZEL DOUBLE-ENDED
PERIOSTEAL ELEVATOR



Round knurled handle
Overall length 195 mm, 7.7 inches

BC-21-550
ELEVATOR
Heavy
7.5 mm curved tip



Overall length 202 mm, 7.9 inches

BC-21-551
LEMPERT ELEVATOR
7 mm curved tip



Overall length 183 mm, 7.2 inches

BC-21-552
HOUSE ELEVATOR
3 mm wide
Octagonal handle



Overall length 180 mm, 7.0 inches

BC-21-553
LEMPERT ELEVATOR
Slightly angled, 2 mm wide



Overall length 175 mm, 6.8 inches



BC-21-554
LEMPERT ELEVATOR
Slight angled, 3mm wide



Overall length 175 mm, 6.8 inches

BC-21-555
HENNER ELEVATOR
Double-ended
2 mm wide



Overall length 189 mm, 7.4 inches



BC-21-556
Double-ended
3 mm round semi sharp end and
2.75 mm elevator end

MCCOLLOUGH
ELEVATOR



Angled 45 degrees
Overall length 185 mm, 7.3 inches

BC-21-557
HALLE ELEVATOR
2 mm tip



Overall length 180 mm, 7.0 inches



BC-21-558
3 mm tip

GORNEY SEPTUM
SUCTION ELEVATOR



Overall length 216 mm, 8.5 inches

BC-21-559
HABERMAN
SUCTION ELEVATOR
2.5 mm x 4 mm opening
Hexagonal handle



Overall length 177 mm, 7.0 inches

BC-21-560
FREER
SEPTUM ELEVATOR
Double-ended
Sharp and blunt 4.5 mm tips



Overall length 190 mm, 7.5 inches

BC-21-561
FREER
SEPTUM ELEVATOR
Double-ended
Sharp and blunt 6.5 mm tips



Overall length 195 mm, 7.6 inches




BC-21-562
Double-ended
4 mm x 10 mm flat edge like scraper and
2 mm x 11 mm semi sharp top

WOODSON ELEVATOR



Overall length 176 mm, 6.9 inches



BC-21-563
Double-ended

BALL END ELEVATOR



Overall length 190 mm, 7.4 inches



BC-21-564

Double-ended
6 mm x 9 mm and 6 mm x 11 mm tip
PENNINGTON
ELEVATOR



Overall length 215 mm, 8.4 inches

BC-21-565

Double-ended
Semi sharp
Curved 3 mm tips
ROGER ELEVATOR




Angled right and left
Overall length 210 mm, 8.2 inches

BC-21-566

Double-ended
Blunt and semi sharp 5 mm tips
HAJEK-BALLENGER
ELEVATOR



Overall length 197 mm, 7.7 inches



BC-21-567


Double-ended
3 mm x 5 mm curette end and 3 mm
spatula end



**DUNNING ELEVATOR
& CURETTE**




Overall length 200 mm, 7.8 inches



BC-21-568

Double-ended
5 mm x 7 mm tips



PIERCE ELEVATOR



Overall length 188 mm, 7.4 inches

BC-21-569
MCKENTY
PERIOSTEAL ELEVATOR
4 mm wide tip



Overall length 124 mm, 4.8 inches

BC-21-570
MCKENTY
PERIOSTEAL ELEVATOR
5 mm wide tip



Overall length 125 mm, 4.9 inches

BC-21-571
CONVERSE-MCKENTY
PERIOSTEAL ELEVATOR
5 mm wide tip



Overall length 145 mm, 5.7 inches

BC-21-572
JOSEPH
PERIOSTEAL ELEVATOR
4 mm wide tip



Overall length 170 mm, 6.6 inches

BC-21-573
CONVERSE
PERIOSTEAL ELEVATOR
4 mm wide tip



Overall length 180 mm, 7.0 inches



BC-21-574
 GILLIES
 ZYGOMA ELEVATOR
 9 mm wide tip
 35 mm from angle to tip



Overall length 225 mm, 8.8 inches

BC-21-575
 COTTLE SKIN ELEVATOR
 10 mm wide tip



Overall length 190 mm, 7.4 inches

BC-21-576
 BOIES NASAL
 FRACTURE ELEVATOR
 9 mm wide tip




Overall length 195 mm, 7.6 inches


BC-21-577
 GOLDMAN DISPLACER
 SEPTUM ELEVATOR
 6 mm wide tip
 Gentle tapered blade



Overall length 190 mm, 7.4 inches



BC-21-578
 Double-ended
 4 mm x 4 mm round end and 3 mm
 spatula end



**COTTLE
 SEPTUM ELEVATOR**



Overall length 225 mm, 8.8 inches

BC-21-579
CLEFT PALATE
ELEVATOR
4 mm x 9 mm tip



Overall length 160 mm, 6.2 inches



BC-21-580
Small
Double-ended

LOUISVILLE ELEVATOR




3 mm x 5 mm and 4 mm x 5 mm tips
Overall length 170 mm, 6.7 inches




BC-21-581
Large
Double-ended

LOUISVILLE ELEVATOR



3 mm pointed end and 5 mm angled end
Overall length 170 mm, 6.7 inches

OTOLOGIC, MICRO



BC-21-582
2 mm curved tip
Round handle

**AUSTIN DUCKBILL
KNIFE-ELEVATOR**



Overall length 159 mm, 6.3 inches



BC-21-583
1.25 mm x 2.5 mm blade
135 degrees back scratcher

PERKINS ELEVATOR



Round handle
Overall length 165 mm, 6.5 inches




BC-21-584
Duckbill type, sharp
3 mm x 8 mm curved blade

**PAPARELLA
ELEVATOR-KNIFE**



Octagonal handle
Overall length 163 mm, 6.5 inches




BC-21-585
2 mm
Angled 13 degrees, 15 mm from angled
to tip

**GUILFORD-WRIGHT
DRUM ELEVATOR**



Round handle
Overall length 140 mm, 5.6 inches




BC-21-586
Duckbill
3.5 mm tip

**SHAMBAUGH &
DERLACKI ELEVATOR**



Light curve
Overall length 172 mm, 6.7 inches




BC-21-587
2 mm diameter blade
Angled 60 degrees
Sharp

**JORDON
CANAL ELEVATOR**



Round handle
Overall length 166 mm, 6.6 inches



BC-21-588
1 mm x 12 mm
Angled 15 degrees
Referred to as the "Gimmick"

HOUSE ELEVATOR



Straight shaft
Round handle
Overall length 168 mm, 6.7 inches




BC-21-589
1.5 mm x 12 mm tip
Angled 15 degrees
coated

**HOUSE TEFLON
COATED ELEVATOR**



With diathermy attachment
Overall length 168 mm, 6.6 inches




BC-21-590
1.5 curved tip
Flat handle, bayonet type

HOUSE T ELEVATOR



88 mm from band to tip
Overall length 170 mm, 6.7 inches



BC-21-591
1 curved tip
Flat handle, bayonet type

HOUSE ELEVATOR



88 mm from band to tip
Overall length 170 mm, 6.7 inches



**WELMED Instrument Company
Defence Road, Near Wapda Grid Station
Sialkot 51310, Pakistan.**

Phone : +92 52 3554782, 3558704

Fax : +92 52 3550296

📞 Mobile : +92 333 8680640

Email : welmedpk@gmail.com

info@welmedinst.com

wmm@brain.net.pk

welmed@cyber.net.pk

Website : www.welmedinst.com