



Serving Humanity Since 1908

WELMED
INSTRUMENTS CO.

**ENT Catalogue
Group 20, 21, 22, 23, 24 & 25**



Manufacturers & Exporters of Quality Surgical Instruments

**Anaesthesia, Otolaryngology, Rhinology, Tonsillectomy,
Laryngo-Bronchoscopy, Oral Maxillo-Facial Surgery
and Cranio-Maxillo-Facial Surgery Instruments**



Tonsil 23

WELMED Instruments Company
Defence Road
Near Wapda Grid Colony
Sialkot 51310, Pakistan.

Phone : +92 52 3554782, 3558704
Fax : +92 52 3550296
Cell : +92 333 8680640

Email : info@welmedinst.com
welmedpk@gmail.com
wmm@brain.net.pk
welmed@cyber.net.pk

Website : www.welmedinst.com

About Us

Quality begins at **WELMED** Instruments Company where master craftsmen forge and manufacture each instrument to its exact specifications and measurements. Instruments must undergo a series of quality control inspections in order to be worthy of **WELMED** name and warranty. We offer a wide selection of medical instruments including Surgical and Dental instruments.

We have acquired a unique and unmatched understanding of manufacturing through a large experience, working since 1908. Years of experience & state of the art manufacturing facilities helped us to keep our prices low. Our Surgical & Dental instruments are created by expert craftsmen, providing unfaltering quality you can rely on.

We have launched a fully interactive website where you can get online and track your orders 24 hours a day, 7 days in a week, and 365 days a year. Ask for our free e-catalogue on CD-Rom to save time and enjoy hassle free ordering.

We provide our customers with most reliable customer service as your satisfaction is our top priority. Our customers always enjoy our wonderful customer service and quality tools at the best price.

Lastly we invite you to have a work experience with **WELMED** Instruments Company and we will prove ourselves as one of the most valuable source you have.

You can find us at

WELMED Instruments Company
Defence Road, Near Wapda Grid Colony
Sialkot 51310, Pakistan.

Phone : +92 52 3554782, 3558704
Fax : +92 52 3550296
Cell : +92 333 8680640

Emails :
Info@welmedinst.com
welmedpk@gmail.com
wmm@brain.net.pk
welmed@cyber.net.pk

Website : www.welmedinst.com



History

WELMED Instruments Company was founded in 1995. Our Great Grandfather **Mian Muhammad Fazal Elahi** was the founder of M/s M.F. Elahi & Company (E.P.) Pvt. Ltd, established in 1908. WELMED is a sister concern of M.F. Elahi & Co. It is our pride to be a member of M.F. Elahi Family.

WELMED Instruments Company is located in the heart of beautiful city, Sialkot, Pakistan. We are manufacturers and exporters of quality Surgical, Dental instruments and Scissors of all sorts in various categories. We are closely attached with highly skilled surgeons and other technicians of the surgical industry to make the new instruments development possible.

WELMED is in a continues process of development and testing of new instruments for their precision and perfect functioning. We are also working at imported forging supply by foreign customers.

Certificates

Welmed Instruments Company has achieved the certificates for ISO-9001, ISO-13485, cGMP (FDA) and CE Marking.

Quality Control

All material use in the manufacturing of our instruments are tested from relevant laboratories. Our continuous quality control is our best warranty so that you built confidence in our products.

Deliveries

Deliveries are quick and in-time.

Mode Of Payment

Terms of Payment is ADVANCE or LETTER OF CREDIT (Irrevocable)

Quality Policy

WELMED Instruments Company aims to produce quality Surgical, Dental instruments and Scissors of all sorts at reasonable prices to become an exemplary manufacturer. Our total emphasize is to achieve customer's satisfaction by clearly understanding their requirements and following a quality management system which is a primary concern of all our work force. We do our best to produce instruments of finest quality.



Warranty

All WELMED Instruments are produced from finest stainless steel available to be free from any defect in material or workmanship. Any product that proves to be defective in material or workmanship will either be repaired or replaced without any charge. This warranty is void if the instrument is abused, altered, modified or is not used in the manner for which it is intended.

Tonsil 23



WELMED
INSTRUMENTS CO.

Adson	159
Andrew-pynchon	167
Blohmke	152
Bruenings (Tube)	155
Bruenings (Tube, Curved)	155
Bruenings	153,156
Byrd	167
Carpenter	149
Chevalier-Jackson	162
Cogswell	168
Cooly	166
Cordes	156
Cottle-Joseph Double Hook	157
Coupland	164
Eves	153
Fraenkel	156
Frazier (Fergusson)	160
Frazier	159,160
Hartmann	154
Henke	150
Huber	155
Hurd	149,150
King Adenold Punch Forcep	154
Krause	156
Krishaber	162
Lange	156
Laryngeal Mirrors, Handles	163
Luer	162
Luniatschek	151
Millard Double Hook Forcep	157
Millard Double Hook	157
Poole	164
Nuboer	168
Rosenberg	156
Scheinmann	156
Schmidt	152
Schumacher	156
Seiffert	156
Tonsilectomes	158
Tyding	153
White	152
Yankauer	165,167



BC-23-101
20.5 cm / 8 "
for right hand

Hurd




BC-23-102
20.5 cm / 8 "
for left hand


Hurd




BC-23-103
25 cm / 9 3/4 "
4.4 & 8.8 mm



Carpenter
Tonsil Knife / Enucleator




BC-23-104
22.5 cm / 8 3/4"
11.0 & 14.0 mm




Hurd



Henke



BC-23-105
22.5 cm / 8 3/4"
10.5 & 5.0 mm




BC-23-106
22.5 cm / 8 3/4"
16.0 & 5.0 mm



Henke





BC-23-107
17.5 cm / 7"
2 mm

BC-23-108
17.5 cm / 7"
3 mm

BC-23-109
17.5 cm / 7"
2 mm

round handle

Luniatschek

Luniatschek



White



1/2

BC-23-110
18 cm / 7"



Blohmke



1/2

BC-23-111
20 cm / 8"



1:2




Schnidt



1/2

BC-23-112
19 cm / 7 1/2"



BC-23-113
18.5 cm / 7 1/4"



Schnidt



1/2

BC-23-114
19 cm / 7 1/2"



BC-23-115
18.5 cm / 7 1/4"



BC-23-116
28.5 cm / 11 1/4"
shaft 6 mm ø

Without ratchet

Eves



BC-23-117
28.5 cm / 11 1/4"
shaft 6 mm ø

With ratchet


Eves



BC-23-118
29 cm / 11 1/2"
shaft end 9.5 mm wide

With double ratchet


Brunings



Tyding
BC-23-119
with one straight tip




120 mm



BC-23-120
12 cm / 4 3/4"
cut

13 mm \varnothing
through cutting



BC-23-121
12 cm / 4 3/4"
cut

15 mm \varnothing
through cutting

Hartmann

Tonsil and Adenoid



BC-23-122
No. 1
Triangular
7 mm lower basket, 6 mm upper jaw
strong curve
45 mm from curve to tip
length of shaft 110 mm
Overall length 179 mm, 7.0 inches




King Adenoid
Punch Forceps



Huber
BC-23-123
Universal handle
for interchangeable tubes ref.



Bruening
Interchangeable tube for handle ref.
BC-23-124 20 cm / 8"
BC-23-125 30 cm / 12"
BC-23-126 40 cm / 16"



Bruening
Laryngeal Tube, Curved
BC-23-127 20 cm / 8 "



B-23-128



Fraenkel



BC-23-129



Seiffert



BC-23-130



Fraenkel



BC-23-131



Scheinmann



BC-23-132



Scheinmann



BC-23-133



BC-23-134



Bruenings



BC-23-135



BC-23-136



Bruenings



BC-23-137



Lange



BC-23-138



BC-23-139



Krause



BC-23-141



BC-23-142



BC-23-143



BC-23-144



BC-23-140



Krause



Cordes



BC-23-145



Schumacher



BC-23-146



BC-23-147



Cordes



BC-23-148



Cordes



BC-23-149




Rosenberg



Overall length 137mm, 5.4 inches

**Cottle-Joseph
Double Hook
BC-23-150**
12 mm wide



Millard Double
Hook Forceps

BC-23-151
with swivel thimble
1 x 2 teeth




10 mm wide double hook overall length
101 mm, 4.0 inches


**Millard
Double Hook
BC-23-153**
10 mm wide double hook thimble




Short flat handles
Overall length 84mm, 3.3 inches




BC-23-154
Fig. 1



BC-23-155
Fig. 2



BC-23-156
Fig. 3



Komplett mit a aufsatzen



Adson
B-23-157
ø 4 mm, Charr.12



Adson
B-23-158
ø 4 mm, Charr.12



Adson
B-23-159
ø 4 mm, Charr.12



Frazier
BC-23-160 Charr. 8
BC-23-161 Charr. 10
BC-23-162 Charr. 12
12.5 cm



Frazier
BC-23-163 Charr. 8
BC-23-164 Charr. 10
BC-23-165 Charr. 12
12.5 cm



Frazier
BC-23-166 Charr. 6
BC-23-167 Charr. 8
BC-23-168 Charr. 10
BC-23-169 Charr. 12
 12.5 cm



Frazier (Fergusson)
BC-23-170 Charr. 6
BC-23-171 Charr. 8
BC-23-172 Charr. 10
BC-23-173 Charr. 12
 12.5 cm



Frazier
BC-23-174 Charr. 6
BC-23-175 Charr. 8
BC-23-176 Charr. 10
BC-23-177 Charr. 12
 18 cm



Frazier (Fergusson)
BC-23-178 Charr. 6
BC-23-179 Charr. 8
BC-23-180 Charr. 10
BC-23-181 Charr. 12
 18 cm



International scale of Trachea Tube Sizes

Size	0	1	2	3	4	5	6	7	8	9
ø mm	4	5	6	7	8	9	10	11	12	13
French	12	15	18	21	24	27	30	33	36	39
Length mm	45	50	55	60	65	70	75	80	85	90



Luer
Trachea Tube,
standard pattern

- | | | |
|----------------------|------------------|------------------------|
| Nickel silver | | Sterling silver |
| BC-23-182 | Size 0-- ø 4 mm | BC-23-192 |
| BC-23-183 | Size 1-- ø 5 mm | BC-23-193 |
| BC-23-184 | Size 2-- ø 6 mm | BC-23-194 |
| BC-23-185 | Size 3-- ø 7 mm | BC-23-195 |
| BC-23-186 | Size 4-- ø 8 mm | BC-23-196 |
| BC-23-187 | Size 5-- ø 9 mm | BC-23-197 |
| BC-23-188 | Size 6-- ø 10 mm | BC-23-198 |
| BC-23-189 | Size 7-- ø 11 mm | BC-23-199 |
| BC-23-190 | Size 8-- ø 12 mm | BC-23-200 |
| BC-23-191 | Size 9-- ø 13 mm | BC-23-201 |



Luer
Trachea Tube,
with opening in inner
and outer tube

- | | |
|----------------------|------------------|
| Nickel silver | |
| BC-23-202 | Size 0-- ø 4 mm |
| BC-23-203 | Size 1-- ø 5 mm |
| BC-23-204 | Size 2-- ø 6 mm |
| BC-23-205 | Size 3-- ø 7 mm |
| BC-23-206 | Size 4-- ø 8 mm |
| BC-23-207 | Size 5-- ø 9 mm |
| BC-23-208 | Size 6-- ø 10 mm |
| BC-23-209 | Size 7-- ø 11 mm |
| BC-23-210 | Size 8-- ø 12 mm |
| BC-23-211 | Size 9-- ø 13 mm |



Krishaber
Trachea Tube,
obturator has 4 fenestrations

Nickel --silver	
BC-23-212	Size 0-- ø 4 mm
BC-23-213	Size 1-- ø 5 mm
BC-23-214	Size 2-- ø 6 mm
BC-23-215	Size 3-- ø 7 mm
BC-23-216	Size 4-- ø 8 mm
BC-23-217	Size 5-- ø 9 mm
BC-23-218	Size 6-- ø 10 mm
BC-23-219	Size 7-- ø 11 mm
BC-23-220	Size 8-- ø 12 mm



Chevalier-Jackson
Trachea Tube,
Nickel silver

BC-23-221	Size 0-- ø 4 mm
BC-23-222	Size 1-- ø 5 mm
BC-23-223	Size 2-- ø 6 mm
BC-23-224	Size 3-- ø 7 mm
BC-23-225	Size 4-- ø 8 mm
BC-23-226	Size 5-- ø 9 mm
BC-23-227	Size 6-- ø 10 mm
BC-23-228	Size 7-- ø 11 mm
BC-23-229	Size 8-- ø 12 mm
BC-23-230	Size 9-- ø 13 mm

Sterling silver	
BC-23-231	
BC-23-232	
BC-23-333	
BC-23-234	
BC-23-235	
BC-23-236	
BC-23-237	
BC-23-238	
BC-23-239	
BC-23-340	



Luer
Speaking Tube,
inner tube with sliding valve
Nickel silver

BC-23-241	Size 4-- ø 8 mm
BC-23-242	Size 5-- ø 9 mm
BC-23-243	Size 6-- ø 10 mm
BC-23-244	Size 7-- ø 11 mm



BC-23-245
18.5 cm / 7 1/4"
Plane, with fixes handle

BC-23-246
18.5 cm / 7 1/4"
rhodium plated



BC-23-247
13.5 cm / 5 1/4"
Plane, without handle

BC-23-248
13.5 cm / 5 1/4"
rhodium plated



BC-23-249
Fig. 1
12 mmø



BC-23-250
Fig. 2
14 mmø



BC-23-251
Fig. 3
16 mmø



BC-23-252
Fig. 4
18 mmø



BC-23-253
10 cm / 4"



BC-23-254
10 cm / 4"



BC-23-255
12.5 cm / 5"

Michel
25 cm / 10"
moveable mirror



Nuboer
BC-23-256
21 cm / 8 1/4"
6 mm ø
round opening



Nuboer
BC-23-257
21 cm / 8 1/4"
6 mm ø
oval opening



BC-23-259
tip laterally
perforated



Coupland
BC-23-258
17 cm / 6 3/4"
Suction Handle with tube connector

BC-23-260
9 cm / 3 1/2"
1.5 mm ø

BC-23-261
9 cm / 3 1/2"
2.0 mm ø

BC-23-262
9 cm / 3 1/2"
2.5 mm ø

BC-23-263
9 cm / 3 1/2"
3.5 mm ø



**Yankauer
BC-23-264**
31 cm / 12 1/4"
complete with suction tip
and connector

**Yankauer
BC-23-265**
31 cm / 12 1/4"
complete with handle, 4 tubes, suction
tip and tube connector

BC-23-272
Suction tip



BC-23-268
Fig. 0
10 cm / 4"



BC-23-269
Fig. 1
10 cm / 4"



BC-23-270
Fig. 2
10 cm / 4"

BC-23-271
Fig. 3
13.5 cm / 5 1/4"



BC-23-266
15 cm / 6"
handle



BC-23-267
2.5 cm / 1"
Olive



1/2

**Cooley
BC-23-273**
36 cm / 14 1/4"
8 mm ø

1/3

**Cooley
BC-23-274**
31 cm / 12 1/4"
10 mm ø
with perforated
screw tip

1/2

**Cooley
BC-23-275**
31 cm / 12 1/4"
7 mm ø
with perforated
screw tip

1/2

**Cooley
BC-23-276**
36 cm / 14 1/4"
8 mm ø
with perforated
screw tip

130 mm

**Cooley
BC-23-277**
25 cm



Yankauer
BC-23-278
29 cm

Ø 6 mm



Andrew-Pynchon
BC-23-279
23 cm / 9"



Byrd
21.5 cm / 8 1/2"

○
BC-23-280
1.5 mm ø

○
BC-23-281
2.5 mm ø

○
BC-23-282
3.5 mm ø



Cogswell
BC-23-283
24.5 cm / 9 3/4"
2.0 mm \varnothing
with finger cut-off



Poole
BC-23-284
22 cm / 8 3/4"
10 mm \varnothing
(Charr. 30)



Poole
BC-23-285
22.5 cm / 8 3/4"
8 mm \varnothing
(Charr. 23)



**WELMED Instrument Company
Defence Road, Near Wapda Grid Station
Sialkot 51310, Pakistan.**

Phone : +92 52 3554782, 3558704

Fax : +92 52 3550296

☎ Mobile : +92 333 8680640

Email : welmedpk@gmail.com

info@welmedinst.com

wmm@brain.net.pk

welmed@cyber.net.pk

Website : www.welmedinst.com