



Serving Humanity Since 1908

WELMED
INSTRUMENTS CO.

**Gynecology Catalogue
Group 9 & 10**



Manufacturers & Exporters of Quality Surgical Instruments

Gynecology & Obstetrics Instruments



Surgical Instruments

Gynecology

9

Obstractics

10

WELMED Instruments Company
Defence Road
Near Wapda Grid Colony
Sialkot 51310, Pakistan.

Phone : +92 52 3554782, 3558704
Fax : +92 52 3550296
Cell : +92 333 8680640

Email : info@welmedinst.com
welmedpk@gmail.com
wmm@brain.net.pk
welmed@cyber.net.pk

Website : www.welmedinst.com

About Us

Quality begins at **WELMED** Instruments Company where master craftsmen forge and manufacture each instrument to its exact specifications and measurements. Instruments must undergo a series of quality control inspections in order to be worthy of **WELMED** name and warranty. We offer a wide selection of medical instruments including Surgical and Dental instruments.

We have acquired a unique and unmatched understanding of manufacturing through a large experience, working since 1908. Years of experience & state of the art manufacturing facilities helped us to keep our prices low. Our Surgical & Dental instruments are created by expert craftsmen, providing unfaltering quality you can rely on.

We have launched a fully interactive website where you can get online and track your orders 24 hours a day, 7 days in a week, and 365 days a year. Ask for our free e-catalogue on CD-Rom to save time and enjoy hassle free ordering.

We provide our customers with most reliable customer service as your satisfaction is our top priority. Our customers always enjoy our wonderful customer service and quality tools at the best price.

Lastly we invite you to have a work experience with **WELMED** Instruments Company and we will prove ourselves as one of the most valuable source you have.

You can find us at

WELMED Instruments Company
Defence Road, Near Wapda Grid Colony
Sialkot 51310, Pakistan.

Phone : +92 52 3554782, 3558704
Fax : +92 52 3550296
Cell : +92 333 8680640

Emails :
Info@welmedinst.com
welmedpk@gmail.com
wmm@brain.net.pk
welmed@cyber.net.pk

Website : www.welmedinst.com



History

WELMED Instruments Company was founded in 1995. Our Great Grandfather **Mian Muhammad Fazal Elahi** was the founder of M/s M.F. Elahi & Company (E.P.) Pvt. Ltd, established in 1908. WELMED is a sister concern of M.F. Elahi & Co. It is our pride to be a member of M.F. Elahi Family.

WELMED Instruments Company is located in the heart of beautiful city, Sialkot, Pakistan. We are manufacturers and exporters of quality Surgical, Dental instruments and Scissors of all sorts in various categories. We are closely attached with highly skilled surgeons and other technicians of the surgical industry to make the new instruments development possible.

WELMED is in a continuous process of development and testing of new instruments for their precision and perfect functioning. We are also working at imported forging supply by foreign customers.

Certificates

WELMED Instruments Company has achieved the certificates for ISO-9001, ISO-13485, cGMP (FDA) and CE Marking.

Quality Control

All material used in the manufacturing of our instruments are tested from relevant laboratories. Our continuous quality control is our best warranty so that you built confidence in our products.

Deliveries

Deliveries are quick and in-time.

Mode Of Payment

Terms of Payment is ADVANCE or LETTER OF CREDIT (Irrevocable)

Quality Policy

WELMED Instruments Company aims to produce quality Surgical, Dental instruments and Scissors of all sorts at reasonable prices to become an exemplary manufacturer. Our total emphasize is to achieve customer's satisfaction by clearly understanding their requirements and following a quality management system which is a primary concern of all our work force. We do our best to produce instruments of finest quality.



Warranty

All WELMED Instruments are produced from finest stainless steel available to be free from any defect in material or workmanship. Any product that proves to be defective in material or workmanship will either be repaired or replaced without any charge. This warranty is void if the instrument is abused, altered, modified or is not used in the manner for which it is intended.

Gynecology 09



WELMED
INSTRUMENTS CO.

Adair	29	Simon	1
Alexander	41	Simpson	27, 1
Aspirating	6	Sims	2, 11, 27
Auvar	19	Somer	37
Auvar-Berlind	19	Staude-Moore	29
Baby-Tischler	42	Teale	34
Barrett	29	Thoms-Gaylor	38
Berger	40	Tischler	42
Blake	1	Tischler-Morgan	40
Braun	30, 31	Townsend	42
Breisky	18	Trelat	8
Bruenings	43	Utrine	4
Bumm	5	Valleix	27
Burke	42	Van Doren	39
Collin	8, 37	Wittner	40
Collin-Landau	8	Schubert	39, 38
Cusco	7, 8		
Czerny	33		
Doyen	16, 17, 35		
Duplay	28		
Emmet	26		
Gourdet	1		
Grave Modif	9		
Grave	9, 10		
Hank	23		
Hartmann	40		
Hegar	21, 22		
Heywood-Smith	35		
Huber Handle	43		
Hulka	28		
Hulka-Kenwick	28		
Jacobs	34		
Kallmorgen	15		
Kelly	35		
Kevourian-Pacific	42		
Kevourkian-Young	36, 41, 44		
Kogan	24, 25		
Krause	41		
Kristeller	11-14		
Lahey	31		
Landau	18		
Laser Surgery Instruments	26		
Mathieu	17		
Matzenauer	18		
Mod. Munchen	36		
Museux	34		
Noto	35		
Novak	36		
Pederson	9, 10		
Pozzi	28		
Pratt	33		
Randall	36		
Recamier	03		
Scheinmann	41		
Scherback	20		
Schroder	1, 30, 31, 32		
Schumacher	41, 44		
Semm	9		



1
BC-09-101
 6.8 mm



2
BC-09-102
 8.5 mm



3
BC-09-103
 10.0 mm



4
BC-09-104
 11.5 mm



5
BC-09-105
 13.0 mm



6
BC-09-106
 14.5 mm



Simon
 25 cm
 Uterine scoops, sharp / rigid
 Cucharillas uterinas, cortant / rigidas
 Curettes uterines, tranchantes / rigides



1
BC-09-107
 6.8 mm



2
BC-09-108
 8.5 mm



3
BC-09-109
 10.0 mm



4
BC-09-110
 11.5 mm



5
BC-09-111
 13.0 mm



6
BC-09-112
 14.5 mm



Schroder
 30 cm
 Uterine scoops, sharp / rigid
 Cucharillas uterinas, cortant / rigidas
 Curettes uterines, tranchantes, rigides



BC-09-113
 6.0 mm



BC-09-114
 8.0 mm



BC-09-115
 10.0 mm



BC-09-116
 12.0 mm



BC-09-117
 14.0 mm



Gourdet
 28 cm
 Uterine scoops, blunt / rigid
 Cucharillas uterinas, romas / rigidas
 Curettes uterines, mousses, rigides

13 mm



Blunt
 romas
 mousse




Sharp
 cortantes
 tranchantes



8 mm

Blake
BC-09-118
 27.5 cm
 Uterine curette
 Cureta uterina
 Curette uterine



Sims

BC-09-119 Fig 1 6.5 mm
BC-09-120 Fig 2 7.5 mm
BC-09-121 Fig 3 9.0 mm
BC-09-122 Fig 4 11.0 mm
BC-09-123 Fig 5 12.0 mm
BC-09-124 Fig 6 14.0 mm
24.5 - 26.0 cm
sharp, rigid
cortantes, rígidas
tranchantes, rigides

Sims

BC-09-125 Fig 1 6.5 mm
BC-09-126 Fig 2 7.5 mm
BC-09-127 Fig 3 9.0 mm
BC-09-128 Fig 4 11.0 mm
BC-09-129 Fig 5 12.0 mm
BC-09-130 Fig 6 14.0 mm
24.5 - 26.0 cm
sharp, malleable
cortantes, maleables
tranchantes, malleables

Sims

BC-09-131 Fig 1 6.5 mm
BC-09-132 Fig 2 7.5 mm
BC-09-133 Fig 3 9.0 mm
BC-09-134 Fig 4 11.0 mm
BC-09-135 Fig 5 12.0 mm
BC-09-136 Fig 6 14.0 mm
24.5 - 26.0 cm
blunt, rigid
romas, rígidas
mousses, rigides

Sims

BC-09-137 Fig 1 6.5 mm
BC-09-138 Fig 2 7.5 mm
BC-09-139 Fig 3 9.0 mm
BC-09-140 Fig 4 11.0 mm
BC-09-141 Fig 5 12.0 mm
BC-09-142 Fig 6 14.0 mm
24.5 - 26.0 cm
blunt, malleable
romas, maleables
mousses, malleables



Recamier
BC-09-143 Fig 00 6.0 mm
BC-09-144 Fig 0 7.5 mm
BC-09-145 Fig 1 8.5 mm
BC-09-146 Fig 2 9.5 mm
BC-09-147 Fig 3 10.5 mm
BC-09-148 Fig 4 11.5 mm
BC-09-149 Fig 5 12.5 mm
BC-09-150 Fig 6 14.0 mm
BC-09-151 Fig 7 15.0 mm
BC-09-152 Fig 8 16.0 mm
BC-09-153 Fig 9 19.0 mm
BC-09-154 Fig 10 20.0 mm
 29.5 - 32.5 cm
 sharp, rigid
 cortantes, rígidas
 tranchantes, rigides



Recamier
BC-09-155 Fig 00 6.0 mm
BC-09-156 Fig 0 7.5 mm
BC-09-157 Fig 1 8.5 mm
BC-09-158 Fig 2 9.5 mm
BC-09-159 Fig 3 10.5 mm
BC-09-160 Fig 4 11.5 mm
BC-09-161 Fig 5 12.5 mm
BC-09-162 Fig 6 14.0 mm
BC-09-163 Fig 7 15.0 mm
BC-09-164 Fig 8 16.0 mm
BC-09-165 Fig 9 19.0 mm
BC-09-166 Fig 10 20.0 mm
 29.5 - 32.5 cm
 sharp, malleable
 cortantes, maleables
 tranchantes, malleables



Recamier
BC-09-167 Fig 00 6.0 mm
BC-09-168 Fig 0 7.5 mm
BC-09-169 Fig 1 8.5 mm
BC-09-170 Fig 2 9.5 mm
BC-09-171 Fig 3 10.5 mm
BC-09-172 Fig 4 11.5 mm
BC-09-173 Fig 5 12.5 mm
BC-09-174 Fig 6 14.0 mm
BC-09-175 Fig 7 15.0 mm
BC-09-176 Fig 8 16.0 mm
BC-09-177 Fig 9 19.0 mm
BC-09-178 Fig 10 20.0 mm
 29.5 - 32.5 cm
 blunt, rigid
 romas, rígidas
 mousses, rigides



Recamier
BC-09-179 Fig 00 6.0 mm
BC-09-180 Fig 0 7.5 mm
BC-09-181 Fig 1 8.5 mm
BC-09-182 Fig 2 9.5 mm
BC-09-183 Fig 3 10.5 mm
BC-09-184 Fig 4 11.5 mm
BC-09-185 Fig 5 12.5 mm
BC-09-186 Fig 6 14.0 mm
BC-09-187 Fig 7 15.0 mm
BC-09-188 Fig 8 16.0 mm
BC-09-189 Fig 9 19.0 mm
BC-09-190 Fig 10 20.0 mm
 29.5 - 32.5 cm
 blunt, malleable
 romas, maleables
 mousses, malleables

6.0
mm



00
BC-09-143
BC-09-155
BC-09-167
BC-09-179

7.5
mm



0
BC-09-144
BC-09-156
BC-09-168
BC-09-180

8.5
mm



1
BC-09-145
BC-09-157
BC-09-169
BC-09-181

9.5
mm



2
BC-09-146
BC-09-158
BC-09-170
BC-09-182

10.5
mm



3
BC-09-147
BC-09-159
BC-09-171
BC-09-183

11.5
mm



4
BC-09-148
BC-09-160
BC-09-172
BC-09-184

12.5
mm



5
BC-09-149
BC-09-161
BC-09-173
BC-09-185

14
mm



6
BC-09-150
BC-09-162
BC-09-174
BC-09-186

15
mm



7
BC-09-151
BC-09-163
BC-09-175
BC-09-187

16
mm



8
BC-09-152
BC-09-164
BC-09-176
BC-09-188

19
mm



9
BC-09-153
BC-09-165
BC-09-177
BC-09-189

20
mm



10
BC-09-154
BC-09-166
BC-09-178
BC-09-190



24 mm

22 mm

20 mm

18 mm

16 mm

14 mm

13 mm



6
 BC-09-197
 BC-09-204

5
 BC-09-196
 BC-09-203

4
 BC-09-195
 BC-09-202

3
 BC-09-194
 BC-09-201

2
 BC-09-193
 BC-09-200

1
 BC-09-192
 BC-09-199



0
 BC-09-191
 BC-09-198



Bumm
 BC-09-191 Fig 0 13 mm
 BC-09-192 Fig 1 14 mm
 BC-09-193 Fig 2 16 mm
 BC-09-194 Fig 3 18 mm
 BC-09-195 Fig 4 20 mm
 BC-09-196 Fig 5 22 mm
 BC-09-197 Fig 6 24 mm
 28 - 29 cm
 sharp, malleable
 cortantes, maleables
 tranchantes, malleables



Bumm
 BC-09-198 Fig 0 13 mm
 BC-09-199 Fig 1 14 mm
 BC-09-200 Fig 2 16 mm
 BC-09-201 Fig 3 18 mm
 BC-09-202 Fig 4 20 mm
 BC-09-203 Fig 5 22 mm
 BC-09-204 Fig 6 24 mm
 28 - 29 cm
 blunt, malleable
 romas, maleables
 mouses, malleables



BC-09-212
Coupling piece for the connection of aspirating curettes BC-09-207 - BC-09-211 to flexible tube having an interior diameter of 11 mm

Adaptador de tubo para curetas de aspiracion BC-09-207 - Bc-09-211 para tubos de aspiracion, \emptyset interior de 11 mm
Embout adaptable aux curettes BC-09-207 - BC-09-211 Pour tube d' aspiration, \emptyset intérieur á 11 mm

BC-09-205 28.5 cm **BC-09-206** 28.5 cm **BC-09-207** 28.5 cm **BC-09-208** 28.5 cm **BC-09-209** 28.5 cm **BC-09-210** 28.5 cm **BC-09-211** 28.5 cm

Specimen and aspirating curettes
Curetats de biopsia y de aspiracion
Curettes á biopsie et d'aspiration

Aspirating curettes
Curetats de aspiracion
Curettes d' aspiration



Cusco
BC-09-213
75 x 17 mm (AxB)
Vaginal speculum for virgins with lateral fixation screw,
tilt-down handle
Especulo vaginal para virgenes, con tornillo de
sujeto lateral, mango plegable
Spéculum vaginal pour vierges, avec fixation latérale á vis



Cusco
BC-09-214 75 x 32 mm (AXB)
BC-09-215 85 x 32 mm (AXB)
BC-09-216 100 x 35 mm (AXB)
Vaginal speculum for virgins with lateral fixation screw,
tilt-down handle
Especulo vaginal para virgenes, con tornillo de
sujeto lateral, mango plegable
Spéculum vaginal pour vierges, avec fixation latérale á vis



Cusco
BC-09-217
75 x 17 mm (AxB)
Vaginal speculum for virgins with central fixation screw,
tilt-down handle
Especulo vaginal para virgenes, con tornillo de
sujeto central, mango plegable
Spéculum vaginal pour vierges, avec fixation centrale á vis



Cusco
BC-09-218 75 x 32 mm (AXB)
BC-09-219 85 x 35 mm (AXB)
BC-09-220 100 x 35 mm (AXB)
Vaginal speculum for virgins with central fixation screw,
tilt-down handle
Especulo vaginal para virgenes, con tornillo de
sujeto central, mango plegable
Spéculum vaginal pour vierges, avec fixation centrale á vis



Cusco

BC-09-221 90 x 22-25 mm (AxB)

BC-09-222 100 x 25-27 mm (AxB)

BC-09-223 105 x 27-30 mm (AxB)

Vaginal specula, Swiss pattern, cone shaped, with lateral fixation screw, tilt-down handle

Especulos vaginales, modelo suizo (cónico), con tornillo de sujeto lateral, mango plegable

Spéculums vaginaux, modèle suisse (conique), avec fixation latéral a vis



Trelat

BC-09-224 85 x 33 mm (AXB)

BC-09-225 95 x 35 mm (AXB)

BC-09-226 115 x 43 mm (AXB)

Vaginal specula, parallel opening, tilt-down handle

Especulos vaginales, apertura paralela, mango plegable

Spéculum vaginaux, à valves parallèles



Collin

BC-09-227 90 x 28 mm (AxB)

BC-09-228 100 x 35 mm (AxB)

BC-09-229 110 x 40 mm (AxB)

Vaginal specula

Especulos vaginales

Spéculums vaginaux



Collin-Landau

BC-09-230 90 x 28 mm (AXB)

BC-09-231 100 x 35 mm (AXB)

BC-09-232 110 x 40 mm (AXB)

Vaginal specula, upper blade of specula shortened

Especulos vaginales, valvas superior más corta

Spéculum vaginaux, valve supérieure plus courte



Semm

BC-09-233 90 x 14-18 mm (AXB)

BC-09-234 100 x 17-20 mm (AXB)

BC-09-235 100 x 25-30 mm (AXB)

Vaginal specula with double lock,
one-hand operated
Especulos vaginales con cierre doble,
para el uso con una mano
Spéculums vaginalux, double ouverture
utilisé à la main



Pederson

BC-09-236 75 x 13 mm (AXB)

BC-09-237 100 x 16 mm (AXB)

BC-09-238 100 x 22 mm (AXB)

BC-09-239 120 x 25 mm (AXB)

Vaginal specula
Especulos vaginales
Spéculums vaginalux



Grave

BC-09-240 75 x 20 mm (AXB)

BC-09-241 95 x 35 mm (AXB)

BC-09-242 115 x 35 mm (AXB)

Vaginal specula
Especulos vaginales
Spéculums vaginalux



Grave modif.

BC-09-243

100 x 35 mm (AXB)

Vaginal speculum, right side open
Especulo vaginal, parte derecha abierta
Spéculum vaginalux, ouverts à droite

Surface black-finish
Superficie negra
Surface noire



Pederson

BC-09-244 75 x 13 mm (AXB)

BC-09-245 100 x 22 mm (AXB)

BC-09-246 120 x 25 mm (AXB)

Vaginal specula for laser surgery,
suctioning from upper blade

Especulos vaginales para la cirugía con laser,
tubo de succión en la valva superior

Spéculums vaginaux pour chirurgie laser, avec
embout aspiration sur la valve supérieure

Pederson

BC-09-247 75 x 13 mm (AXB)

BC-09-248 100 x 22 mm (AXB)

BC-09-249 120 x 25 mm (AXB)

Vaginal specula for laser surgery,
suctioning for lower blade

Especulos vaginales para la cirugía con laser,
tubo de succión en la valva inferior

Spéculums vaginaux pour chirurgie laser, avec
embout aspiration sur la valve inférieure

Surface black-finish
Superficie negra
Surface noire



Grave

BC-09-250 75 x 20 mm (AXB)

BC-09-251 95 x 35 mm (AXB)

BC-09-252 115 x 35 mm (AXB)

Vaginal specula for laser surgery,
suctioning for lower blade

Especulos vaginales para la cirugía con laser,
tubo de succión en la valva superior

Spéculums vaginaux pour chirurgie laser, avec
embout aspiration sur la valve supérieure

Grave

BC-09-253 75 x 20 mm (AXB)

BC-09-254 95 x 35 mm (AXB)

BC-09-255 115 x 35 mm (AXB)

Vaginal specula for laser surgery,
suctioning for lower blade

Especulos vaginales para la cirugía con laser,
tubo de succión en la valva inferior

Spéculums vaginaux pour chirurgie laser, avec
embout aspiration sur la valve inférieure



1/2

Sims
BC-09-256 65 x 25 mm / 70 x 30 mm
BC-09-257 70 x 30 mm / 70 x 35 mm
BC-09-258 85 x 35 mm / 90 x 40 mm



Kristeller
BC-09-259 Fig. 1 70 x 27 mm
BC-09-260 Fig. 2 80 x 30 mm
BC-09-261 Fig. 3 90 x 36 mm



Kristeller
BC-09-262 Fig. 1 75 x 23 mm
BC-09-263 Fig. 2 85 x 26 mm
BC-09-264 Fig. 3 100 x 32 mm



Kristeller
BC-09-265 = Consisting of
BC-09-266 = Consiste de
BC-09-267 = Se composant de
 Sets **BC-09-259 + BC-09-262**
BC-09-260 + BC-09-263
BC-09-261 + BC-09-264

Slender pattern
modelo delgado
modèle gracile



Kristeller
BC-09-268 Fig. 1 105 x 18 mm
BC-09-269 Fig. 2 105 x 21 mm

Kristeller
BC-09-270 Fig. 1 110 x 14 mm
BC-09-271 Fig. 2 110 x 17 mm

Kristeller
BC-09-272 = Consisting of
BC-09-273 = Consiste de
Sets = Se composant de
BC-09-270 + BC-09-268
BC-09-261 + BC-09-269



Flat speculum
 espéculo plano
 spéculum plat



Kristeller
BC-09-274 Fig. 1 105 x 26 mm
BC-09-275 Fig. 2 105 x 28 mm
BC-09-276 Fig. 3 105 x 33 mm



Kristeller
BC-09-277 Fig. 1 115 x 22 mm
BC-09-278 Fig. 2 115 x 24 mm
BC-09-279 Fig. 3 115 x 29 mm



Kristeller
BC-09-280 = Consisting of
BC-09-281 = Consiste de
BC-09-282 = Se composant de
 Sets **BC-09-274 + BC-09-277**
BC-09-275 + BC-09-278
BC-09-276 + BC-09-279



Deep speculum
espéculeo profundo
spéculum profond



Kristeller
BC-09-283 Fig. 1 110 x 27 mm
BC-09-284 Fig. 2 110 x 30 mm
BC-09-285 Fig. 3 110 x 36 mm



Kristeller
BC-09-286 Fig. 1 115 x 22 mm
BC-09-287 Fig. 2 115 x 24 mm
BC-09-288 Fig. 3 115 x 29 mm



Kristeller
BC-09-289 = Consisting of
BC-09-290 = Consiste de
BC-09-291 = Se composant de
Sets = BC-09-283 + BC-09-286
BC-09-284 + BC-09-287
BC-09-285 + BC-09-288




Kallmorgen
BC-09-292 Fig. 1 70 x 40 mm
BC-09-293 Fig. 2 90 x 40 mm



Kallmorgen
BC-09-294 Fig. 1 70 x 33 mm
BC-09-295 Fig. 2 90 x 33 mm



Kallmorgen
BC-09-296 = Consisting of
BC-09-297 = Consiste de
Sets = Se composant de
 = **BC-09-292 + BC-09-294**
 = **BC-09-293 + BC-09-295**



Slightly concave
ligeramente concave
légèrement concave

BC-09-298 Fig. 1 63 x 35 mm
BC-09-299 Fig. 2 94 x 35 mm
BC-09-300 Fig. 3 124 x 35 mm

BC-09-301 Fig. 1 66 x 47 mm
BC-09-302 Fig. 2 95 x 47 mm
BC-09-303 Fig. 3 126 x 47 mm
BC-09-304 Fig. 4 160 x 47 mm

BC-09-305 Fig. 1 55 x 60 mm
BC-09-306 Fig. 2 95 x 60 mm
BC-09-307 Fig. 3 130 x 60 mm
BC-09-308 Fig. 4 175 x 60 mm

Doyen



Mathieu
BC-09-309 95 x 13-23 mm (AXB)
BC-09-310 95 x 16-30 mm (AXB)
Specula only
Especulos sueltos
Spéculums seuls



Mathieu
BC-09-311 100 x 11-21 mm (AXB)
BC-09-312 100 x 14-27 mm (AXB)
Retractors only
Valvas sueltas
Valves seules



Doyen
BC-09-313 60 x 30 mm (AXB)
BC-09-314 70 x 30 mm (AXB)
BC-09-315 80 x 30 mm (AXB)



Matzenauer

- BC-09-316** 125 x 15 mm (AXB)
- BC-09-317** 125 x 20 mm (AXB)
- BC-09-318** 125 x 25 mm (AXB)
- BC-09-319** 125 x 30 mm (AXB)
- BC-09-320** 125 x 35 mm (AXB)
- BC-09-321** 125 x 40 mm (AXB)



Landau

- BC-09-322** 100 x 28 mm (AXB)
- BC-09-323** 110 x 38 mm (AXB)



Breisky

- BC-09-324** 100 x 15 mm (AXB)
- BC-09-325** 100 x 20 mm (AXB)
- BC-09-326** 100 x 25 mm (AXB)
- BC-09-327** 100 x 30 mm (AXB)
- BC-09-328** 100 x 35 mm (AXB)
- BC-09-329** 100 x 40 mm (AXB)
- BC-09-330** 130 x 15 mm (AXB)
- BC-09-331** 130 x 20 mm (AXB)
- BC-09-332** 130 x 25 mm (AXB)
- BC-09-333** 130 x 30 mm (AXB)
- BC-09-334** 130 x 35 mm (AXB)
- BC-09-335** 130 x 40 mm (AXB)
- BC-09-336** 160 x 40 mm (AXB)
- BC-09-337** 180 x 40 mm (AXB)



**Auvard-Berlind
BC-09-338**
102 x 45 mm (AXB)
1.14 kg
Vaginal speculum with fixed
weight and holes in blade
Especulo vaginal con peso
fijo y valva perforada
Spéculum vaginal avec poids
fixe et valve à trous



**Auvard
BC-09-339**
83 x 38 mm (AXB)
0.91 kg
BC-09-340
83 x 38 mm (AXB)
1.14 kg
BC-09-341
83 x 38 mm (AXB)
1.36 kg
Vaginal speculum with fixed weight
Especulo vaginal con peso fijo
Spéculum vaginal avec poids fixe



**Auvard
BC-09-342**
102 x 45 mm (AXB)
1.14 kg
Vaginal speculum with detachable weight
Especulo vaginal con peso desmontable
Spéculum vaginal avec poids démontable



**Scherback
BC-09-343**

Set, consisting of:

- BC-09-344** 1 Handle
- BC-09-345** 1 Weight 0.4 kg
- BC-09-347** 1 Blade 80 x 30 mm
- BC-09-348** 1 Blade 85 x 35 mm
- BC-09-349** 1 Blade 85 x 40 mm
- BC-09-350** 1 Blade 55 x 40-60 mm
- BC-09-351** 1 Blade 70 x 45-70 mm

Jeu, composé de:

- BC-09-344** 1 Manche
- BC-09-345** 1 Poids 0.4 kg
- BC-09-347** 1 Valve 80 x 30 mm
- BC-09-348** 1 Valve 85 x 35 mm
- BC-09-349** 1 Valve 85 x 40 mm
- BC-09-350** 1 Valve 55 x 40-60 mm
- BC-09-351** 1 Valve 70 x 45-70 mm

Juego, compuesto de:

- BC-09-344** 1 Mango
- BC-09-345** 1 Peso 0.4 kg
- BC-09-347** 1 Valva 80 x 30 mm
- BC-09-348** 1 Valva 85 x 35 mm
- BC-09-349** 1 Valva 85 x 40 mm
- BC-09-350** 1 Valva 55 x 40-60 mm
- BC-09-351** 1 Valva 70 x 45-70 mm



Hegar

- BC-09-352 ø 1.0 mm
- BC-09-353 ø 2.0 mm
- BC-09-354 ø 3.0 mm
- BC-09-355 ø 4.0 mm
- BC-09-356 ø 5.0 mm
- BC-09-357 ø 6.0 mm
- BC-09-358 ø 7.0 mm
- BC-09-359 ø 8.0 mm
- BC-09-360 ø 9.0 mm
- BC-09-361 ø 10.0 mm
- BC-09-362 ø 11.0 mm
- BC-09-363 ø 12.0 mm
- BC-09-364 ø 13.0 mm
- BC-09-365 ø 14.0 mm
- BC-09-366 ø 15.0 mm
- BC-09-367 ø 16.0 mm
- BC-09-368 ø 17.0 mm
- BC-09-369 ø 18.0 mm
- BC-09-370 ø 19.0 mm
- BC-09-371 ø 20.0 mm
- BC-09-372 ø 21.0 mm
- BC-09-373 ø 22.0 mm
- BC-09-374 ø 23.0 mm
- BC-09-375 ø 24.0 mm
- BC-09-376 ø 25.0 mm
- BC-09-377 ø 26.0 mm

18.5 cm
 with sloped handle
 con mango aplanado
 á manche aplati



Hegar

- BC-09-378 ø 1.5 mm
- BC-09-379 ø 2.5 mm
- BC-09-380 ø 3.5 mm
- BC-09-381 ø 4.5 mm
- BC-09-382 ø 5.5 mm
- BC-09-383 ø 6.5 mm
- BC-09-384 ø 7.5 mm
- BC-09-385 ø 8.5 mm
- BC-09-386 ø 9.5 mm
- BC-09-387 ø 10.5 mm
- BC-09-388 ø 11.5 mm
- BC-09-389 ø 12.5 mm
- BC-09-390 ø 13.5 mm
- BC-09-391 ø 14.5 mm
- BC-09-392 ø 15.5 mm
- BC-09-393 ø 16.5 mm
- BC-09-394 ø 17.5 mm

18.5 cm
 with sloped handle
 con mango aplanado
 á manche aplati



Hegar

- BC-09-395 ø 2.0 / 3.0 mm
- BC-09-396 ø 3.0 / 4.0 mm
- BC-09-397 ø 5.0 / 6.0 mm
- BC-09-398 ø 7.0 / 8.0 mm
- BC-09-399 ø 9.0 / 10.0 mm
- BC-09-400 ø 11.0 / 12.0 mm
- BC-09-401 ø 13.0 / 14.0 mm
- BC-09-402 ø 15.0 / 16.0 mm
- BC-09-403 ø 17.0 / 18.0 mm
- BC-09-404 ø 19.0 / 20.0 mm
- BC-09-405 ø 21.0 / 22.0 mm
- BC-09-406 ø 23.0 / 24.0 mm
- BC-09-407 ø 25.0 / 26.0 mm

20 cm



With sloped handle
con mango aplanado
á manche aplati

Hegar

- BC-09-408**
- BC-09-409** 1 Rack for dilators \varnothing 4-12 mm
- BC-09-410** 1 Rack for dilators \varnothing 13-17 mm
- BC-09-411** 1 Metal case with hinged lid 205 x 105 x 45 mm (wxdxh)

- BC-09-409** 1 Soporte para dilatadores \varnothing 4-12 mm
- BC-09-410** 1 Soporte para dilatadores de \varnothing 13-17 mm
- BC-09-411** 1 Estuche de metal con tapa de bisagra 205 x 105 x 45 mm

- BC-09-409** 1 Chevalet pour dilateurs \varnothing 4-12 mm
- BC-09-410** 1 Chevalet pour dilateurs \varnothing 13-17 mm
- BC-09-411** 1 Boîte métal avec couvercle á charnière 205 x 105 x 45 mm



Hegar

- BC-09-412**
- BC-09-413** 1 Rack for dilators \varnothing 1-17 mm
- BC-09-414** 1 Rack for dilators \varnothing 18-26 mm
- BC-09-415** 1 Metal case with hinged lid 310 x 210 x 55 mm (wxdxh)


- BC-09-413** 1 Soporte para dilatadores de \varnothing 1-17 mm
- BC-09-414** 1 Soporte para dilatadores de \varnothing 18-26 mm
- BC-09-415** 1 Estuche de metal con tapa de bisagra 310 x 210 x 55 mm

- BC-09-413** 1 chevalet pour dilateurs de \varnothing 1-17 mm
- BC-09-414** 1 chevalet pour dilateurs de \varnothing 18-26 mm
- BC-09-415** 1 Boîte métal avec couvercle á charnière 310 x 210 x 55 mm



Hegar

- BC-09-416**
- BC-09-417** 1 Boilable canvas roll
- BC-09-417** 1 Bolsa de tela esterilizable por ebullición
- BC-09-417** 1 Trousse stérilisable par ébullition




Hank
BC-09-418
set of 6 piezas
in boilable canvas roll
juego de 6 piezas
en bolsa de tela esterilizable por ebullición
jeu de 6 pièces en pochette en treillis stériliable par ébullition

Hank
BC-09-419 4.5 / 5.0 mm \emptyset
BC-09-420 5.5 / 6.0 mm \emptyset
BC-09-421 6.5 / 7.0 mm \emptyset
BC-09-422 7.5 / 8.0 mm \emptyset
BC-09-423 8.5 / 9.0 mm \emptyset
BC-09-424 9.5 / 10.0 mm \emptyset

American scale
escala americana
échelle américaine



BC-09-427
3 mm



BC-09-428
5 mm

5 mm



Kogan
BC-09-425
24 cm
Endospeculum for cervix and urethra, without catch
Endospéculo para uretra y cerviz, sin cierre
Endospeculum pour urètre et col utérin, sans verrou

5 mm



Kogan
BC-09-426
24 cm
Endospeculum for cervix and urethra, with
scale-ratchet and fixing screws
Endospéculo para uretra y cerviz, con cierre
graduado tornillo fijador
Endospeculum pour uretre et col utérin
crémaille à vis



Kogan
BC-09-427
BC-09-428
27 cm
Endospecula for cervix and
urethra, with catch
Endoespeculos para uretra
y cerviz, sin cierre
Endospeculums pour uretre
et col utérin avec crémaille,
avec verrou



BC-09-430
3 mm



BC-09-431
5 mm



BC-09-429
25 cm
Lateral speculum
Espécule latéral
Speculum latéral



KOGAN
BC-09-430
BC-09-431
27 cm
Endospecula for cervix and urethra
Endoespécules para cerviz y uretra
Speculums pour urètre et col utérin



BC-09-432
24,5 cm
Hook
Gancho
Crochet

EMMET
BC-09-433
24,5 cm
Hook
Gancho
Crochet

BC-09-434
24,5 cm
Hook
Gancho
Crochet

BC-09-435
23,5 cm
Hook
Gancho
Crochet

BC-09-436
26 cm
Hook
Gancho
Crochet

Ø 4.0 mm



**Sims
BC-09-437**
34 cm
graduated in inch,
malleable, silve-plated

graduado en pulgadas,
maleable, plateado

gradu , en pouces,
mall able, argent

Ø 4.0 mm



**Sims
BC-09-438**
34 cm
graduated in centimeters,
malleable, silver-plated

graduado en centimetros,
maleable, plateado

gradu , en cm,
mall able, argent

Ø 4.0 mm



**Sims
BC-09-439**
34 cm
graduated in centimeters,
malleable, 18/8 steel

graduado en centimetros,
maleable, acero 18/8

gradu , en cm,
mall able, acier inox

Ø 4.0 mm



**Simpson
BC-09-440**
34 cm
graduated in centimeters,
malleable, silver-plated

graduado en centimetros,
maleable, plateado

gradu , en cm,
mall able, argent


Ø 4.0 mm



**Valleix
BC-09-441**
29 cm
graduated in centimeters, with silder,
malleable, chromium plated

graduado en centimetros, con un cursor,
maleable, cromado

gradu , en cm avec curseur,
mall able, chrom 



Pozzi
BC-09-442
25.5 cm
heavy pattern
modelo fuerte
modèle fort

Duplay
BC-09-443
28.5 cm

Hulka
BC-09-444
29 cm

Hulka-Kenwick
BC-09-445
29 cm



Staude-Moore
BC-09-446
17 cm



Barrett
BC-09-447
18 cm



Adair
Bc-09-448
18.5 cm



Schröder
BC-09-449
25.5 cm



Schröder
BC-09-450
25 cm



Braun
BC-09-451
25 cm



Lahey
BC-09-452
15.5 cm



Braun
BC-09-453
15.5 cm



Schröder
BC-09-454
25 cm
Straight
recta
droits



Schröder
BC-09-455
25 cm
curved on flat
curva sobre el plano
courbe sur le plat



Schröder
BC-09-456
25 cm
curved sideways
curva lateral
courbe sur le côté



Czerny
BC-09-457
21 cm



Pratt
BC-09-458
26 cm



Pratt
BC-09-459
27.5 cm



Jacobs
BC-09-460
21 cm
Straight
recta
droit



Teale
BC-09-461
BC-09-462
Curved side ways
curva lateral
courbe sur le côté



BC-09-461
23 cm



BC-09-462
23 cm



Museux
BC-09-463
BC-09-464
Straight
rectas
droits



BC-09-463
23.5 cm



BC-09-464
23.5 cm




Museux
BC-09-465
BC-09-466
curved sideways
curva lateral
courbe sur le côté



BC-09-465
23.5 cm



BC-09-466
23.5 cm




Doyen
BC-09-467
18.5 cm

Heywood-Smith
BC-09-468
25 cm

Noto
BC-09-469
27 cm

Kelly
BC-09-470
31 cm



**Randall
BC-09-471**
25 cm
Endometrical suction
biopsy curette

Cureta aspiradora para
biopsia endo-uterina

Curette aspiratrice pour
biopsie de l'endomètre

**Novak
BC-09-472**
25 cm
Endometrical suction
biopsy curette

Cureta aspiradora para
biopsia endo-uterina

Curette aspiratrice pour
biopsie de l'endomètre

**Mod.München
BC-09-473**
27 cm
Endometrical biopsy curette

Cureta para biopsia endo-uterina

Curette biopsie de l'endomètre

**Kevourkian-Young
BC-09-474**
28cm
Endocervical biopsy
curette

Cureta para biopsia
endo-cerviz

Curette á biopsie
endo-cervicale

**Kevourkian-Young
BC-09-475**
28cm
Endocervical biopsy
curette

Cureta para biopsia
endo-cerviz

Curette á biopsie
endo-cervicale



With movable jaws
con mandibula movable
avec mors mobiles



Somer
BC-09-476
23 cm
Uterine elevating forceps
Pinza para el útero
Pince hystérolabe



Collin
BC-09-477
26.5 cm
Uterine elevating forceps
Pinza para el útero
Pince hystérolabe



Thoms-Gaylor
BC-09-478
24 cm
Cervical biopsy and specimen forceps
Pinza saca-bocados para biopsia
Pince emporte-pieces á biopsie






Schubert
BC-09-479
28.5 cm
Uterine Biopsy Forceps






Schubert
BC-09-480
26 cm
Uterine Biopsy Forceps
Handle Angled To Side




Van Doren
BC-09-481
25 cm
Uterine Biopsy Forceps
S-Shaped







Burke
BC-09-492 Shaft size 18 cm
BC-09-492A Shaft size 20 cm
BC-09-492B Shaft size 25 cm



Kevourkian-Pacific
BC-09-493 Shaft size 18 cm
BC-09-493A Shaft size 20 cm
BC-09-493B Shaft size 25 cm

Baby-Tischler
BC-09-494 Shaft size 18 cm
BC-09-494A Shaft size 20 cm
BC-09-494B Shaft size 25 cm

Tischler-Morgan
BC-09-495 Shaft size 18 cm
BC-09-495A Shaft size 20 cm
BC-09-495B Shaft size 25 cm




BC-09-496
 140 mm
 Universal Handle
 With complete



Huber BC-09-497
 140 mm
 Universal Handle
 With complete




Bruenings
BC-09-498 200 mm
BC-09-499 250 mm




Townsendmini
BC-09-500

2.3 x 4.2 mm
 straight



Tischler kevorkian
BC-09-501

3 x 9.5 mm
 straight



Baggish
BC-09-502

4 x 6.5 mm
 straight



4.6 mm / A
5.4 mm / B

BC-09-503
25 cm
complete
completo
complet

Tip only
Punta suelta
Embout, seul

Handle
mango
manche

Schumacher

A



3.7 x 5 mm

**Townsend
BC-09-504**
Straight
recta
droit



3.7 x 5 mm

**Townsend
BC-09-505**
angled upward
angulada hacia arriba
angulé au dessus



3.7 x 5 mm

**Townsend
BC-09-505**
angled downward
angulada hacia abajo
angulé au dessous




3 x 7 mm

**Tischler
BC-09-507**
Straight
recta
droit



3 x 7 mm

**Tischler
BC-09-508**
angled upward
angulada hacia arriba
angulé au dessus



Shaft size 20 cm
shaft size 23 cm
shaft size 25 cm

Handle
mango
manche

B



3 x 7 mm

**Tischler
BC-09-509**
angled downward
angulada hacia abajo
angulé en dessous



3 x 8 mm / A
3 x 9 mm / B

**Kevourkian-Younger
BC-09-510**
straight
recta
droit



Simpson-Braun
BC-10-101
31 cm



Simpson-Braun
BC-10-102
36 cm



Naegele
BC-10-103
35.5 cm



Naegele
BC-10-104
40 cm



Kieland
BC-10-105
41 cm



Kieland-Luikart
BC-10-106
41 cm



Bill
BC-10-107
Traction handle to be used with all forceps with standard type handle
Tractor adecuado para todas las fórceps con mango standard
Tracteurs pour tous les forceps avec manches ordinaires



Wrigley
BC-10-108
28 cm
Forceps for caesarean section
Fórceps para operación cesárea
Forceps pour césarienne



Piper
BC-10-109
44.5 cm
For delivery of the after coming head
in breech presentation
especialmente para partos en
presentación de nalgas
Spatules pour dégager la tête



Tarnier
BC-10-110
40 cm
Forceps with traction rod
Fórceps con dispositivo de tracción
Forceps avec dispositif de traction



Sellheim
BC-10-111
32 cm
Obstetric lever
Separador para la salida de cabeza
Levier obstétrique



Green-Armytage
BC-10-112
20.5 cm
Caesarean section haemostatic forceps
Pinza-clamp uterina para hemostasis
Pince-clamp utérine pour hémostase



Willett
BC-10-113
19 cm
Scalp flap forceps
Pinza para pericráneo
Pince pour le cuir chevelu



Murless
BC-10-114
30 cm
Head extractor for use in lower segment in caesarean operation
Elevador para cabeza usándose en le segmento inferior del útero en casos de operaciones cesáreas
Levier á tête utilisé en cas de césarienne



Winter

Without ratchet
sin cremallera
sans crémaillère



BC-10-116
28.5 cm



BC-10-117
28.5 cm



BC-10-118
28.5 cm



BC-10-119
28.5 cm



BC-10-120
28.5 cm



BC-10-121
28.5 cm



Winter

Without ratchet
sin cremallera
sans crémaillère



Without ratchet
sin cremallera
sans crémaillère



Saenger
BC-10-122
29 cm



With ratchet
con cremallera
avec crémaillère



Saenger
BC-10-123
27.5 cm



With ratchet
con cremallera
avec crémaillère



Saenger
BC-10-124
27.5 cm



14 mm



BC-10-125
14 mm

18 mm



BC-10-126
18 mm

21 mm



BC-10-127
21 mm

Cuzzi
BC-10-125 14 mm
BC-10-126 18 mm
BC-10-127 21 mm
30 cm

Placenta & ovum scoops
Cucharillas para aborto y placenta
Curettes á avortement et placenta



Bumm
BC-10-128 24 mm
BC-10-129 26 mm
BC-10-130 28 mm
 31 cm
 sharp, malleable
 cortantes, maleable
 trachentes, malléables

Bumm
BC-10-131 24 mm
BC-10-132 26 mm
BC-10-133 28 mm
 31 cm
 blunt, malleable
 romas, maleable
 mouses, malléables

Recamir
BC-10-134 35 mm
BC-10-135 40 mm
 33.5 cm
 Placenta and ovum curettes, blunts / rigid
 Curetas para aborto y placenta, romas / rigidas
 Curettes á avortement et placenta, mousse / rigide



**WELMED Instrument Company
Defence Road, Near Wapda Grid Station
Sialkot 51310, Pakistan.**

Phone : +92 52 3554782, 3558704

Fax : +92 52 3550296

📞 Mobile : +92 333 8680640

Email : welmedpk@gmail.com

info@welmedinst.com

wmm@brain.net.pk

welmed@cyber.net.pk

Website : www.welmedinst.com