



WELMED

INSTRUMENTS CO.

MANUFACTURERS OF ORTHODONTIC INSTRUMENTS

ORTHODONTICS



www.welmedinst.com



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodónticos
Ortodonti Aletleri



DI-17-101 TC
Wire Cutter upto
.022" x .028" Ligature Wires
(0.22 mm x 0.51 mm)



DI-17-102 TC
Hard Wire Cutter
.022" x .028"
(0.56 mm x 0.71 mm)



DI-17-103 TC
Side Nipper for Lock Pins and
Ligatures for Soft Thin Wire
Ø .016" (0.4 mm)



DI-17-104 TC
Wire Cutter Angled 45°
upto .015" (0.38 mm)



DI-17-105 TC
Micro Light Wire Cutter upto
Max 0.15 mm , 14 cm



DI-17-106 TC
Micro Distal & Cutter
Safety Hold Hard Wire
.022" x .025"
(0.55 mm x 0.64 mm)



DI-17-107 TC
Wire Cutter upto
.014" (0.36 mm)

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-108 TC
Distal End Cutter Safety Hold
Re : .022" x .028" (0.55 mm x 0.71 mm)
Ro : .020" x .016" (0.51 mm x 0.41 mm)



DI-17-111 TC
Distal End Cutter
Safety Holding (Round Shape)
Re : 0.56 mm x 0.75 mm
Ro : 0.51 mm



DI-17-109 TC
Micro Distal End Cutter - Safety Hold
Re : upto .022" x .028" (0.55 mm x 0.71 mm)
Ro : .020" (0.50 mm)



DI-17-112
Posterior Band Remover



DI-17-113
Direct Bonding Bracket
Remover with Pad



DI-17-110 TC
Micro Distal End Cutter
Safety Hold
upto 0.21 mm x 0.70 mm
14 cm



DI-17-114 TC
Detailing Step Plier Special
Double Side Offset Head
Forms Bayonet Bands,
Forms Right Hand or Left Hand
Step-up or Step-Down
Bends In 1 mm Increments.
Remover with Pad



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-115
Direct Bonding Attachment
Bracket Remover (Angled) with Stop Screw



DI-17-116
Direct Bonding Bracket
Bracket Remover Narrow, with Stop Screw



DI-17-117
TC and without TC
Tweed Angle Arch Banding
Blade 14 mm



DI-17-118
TC and without TC
Tweed Angle Arch Banding
Blade 10 mm



DI-17-119
Bend Seating Plier Will Firmly
Grip The Bracket's Tie Wing Aiding
In The Seating of The Band. Designed
with Serrated Tips For Maximum Grip
That Prevent Slippage During use.



DI-17-120
TC and without TC
Tweed Angle Arch Banding
Blade 6.5 mm



DI-17-121
TC and without TC
Arch Forming Plier



DI-17-122
Legature Arch Banding Plier

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-123
TC and without TC
Utility Plier Slim Tip



DI-17-124
TC and without TC
Bend Seating Plier (Dir)



DI-17-125
TC and without TC
How Plier with a 45° Angle
of Tips (Cur)



DI-17-126
TC and without TC
Weingart Well Grip Plier
with Serrated




DI-17-127
Tweed Plier Loop Forming
with Stop Screw
Re : 0.56 mm x 0.71 mm



DI-17-128
Tweed Plier Loop Forming
with Stop Screw
Re : 0.56 mm x 0.71 mm



DI-17-129
Angled Bell Hook Crimping Plier



DI-17-130
Bell Hook Crimping Plier




Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodónticos
Ortodonti Aletleri



DI-17-131
ADERER
Three Jaw Wire Bending Plier
(Thin Large Points)




DI-17-132
ADERER
Three Jaw Wire Bending Plier
(Medium Points)




DI-17-133
YOUNG
Loop Bending Plier




DI-17-134
JOHNSON
Contouring Plier



DI-17-135
Nance Closing Loop Plier



DI-17-136
ADERER
Three Jaw Wire Bending Plier




DI-17-137
Light Wire Pliers with Cutter Bezz



DI-17-138
Clasp Bending Plier



DI-17-139
Arch Forming Plier



DI-17-140
Arch Forming Plier (Plain)

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-141
Palatal Arch Plier
Allow Easy Placement or Removal



DI-17-142
Flat Nose Plier with Smooth
Breaks for Bending




DI-17-143
Adams Plier for Wire Bending




DI-17-144
Round Nose Plier
for Wire Bending



DI-17-145
Double Rounded Jaw
Wire Bending Plier



DI-17-146
Mckellops
16 cm



DI-17-147
Mckellops
14 cm



DI-17-148
Reynolds



DI-17-149
Bird Beak Plier Bending Round,
Square or Rectangular



DI-17-150
Bird Beak Plier Bending Round,
Square or Rectangular



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodónticos
Ortodonti Aletleri



DI-17-151
Universal




DI-17-152
New Universal



DI-17-153
Nance 15 cm




DI-17-154
Schwarz



DI-17-155
Waldschs 15 cm




DI-17-156
Heavy Wire Cutter
for Hard Wires upto \varnothing 1.5 mm



DI-17-157
Niti Tip Back Plier




DI-17-158
Aderer 14 cm




DI-17-159
Versatile Light Wire Bending Plier.
Grooves Assist with Precise Bending
and Closing of Loops. Serrated Surface
Grip Wire upto .020" (51 mm)

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri




DI-17-160 **DI-17-161** **DI-17-162** **DI-17-163** **DI-17-164**



DI-17-165
 Distal Bender Simple to use
 Bends tails of springs Works on all Stainless
 Steel and most NiTi and Beta III Wires



DI-17-166
 NITI Ice Instruments Turns NITI wire soft and
 flexible Immediately Designed for Intraoral use
 Durable Stainless Steel design



DI-17-167
 Positive Positioner Precise Placement
 of any direct bond attachment
 Light weight design Hard Anodized aluminum



DI-17-169
 Bracket Positioning Gauge



DI-17-168
 Taylor Tweezer 15 cm



DI-17-170 Twirl-On Director Applicator

DI-17-171 Legatur Director

DI-17-172 Legatur Director

DI-17-173 Band Pusher and Schure Scaler

DI-17-174 Band Pusher and Schure Scaler



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodónticos
Ortodonti Aletleri




DI-17-176
Beebee
(Cur)

DI-17-175
Beebee
(Dir)



DI-17-177
Universal



DI-17-179
Mosquito Forceps
(Cur)

DI-17-178
Mosquito Forceps
(Dir)

DI-17-180 A .018"
DI-17-180 B .022"

DI-17-180 A
Aluminum Height Gauge



DI-17-181
Mershon Band
Pusher



DI-17-182
Bend Setter with
Tin Insert 15 cm



DI-17-183
Bend Pusher
Plastic Handle



DI-17-184
Mathieu Without TC
14 cm



DI-17-185
Mathieu
14 cm



Mathieu With TC
14 cm



DI-17-184T
THIN



DI-17-184K
KOCHER



DI-17-184B
HEAVY

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-186
Castroviejo



DI-17-187
Bone Caliper Double Scale
15.5 cm



DI-17-190
Boley



DI-17-188
Separating Plier for Separators



DI-17-189
Ligature Tying Plier



DI-17-191
Beerendonk



WELMED

INSTRUMENTS CO.

Defence Road, Near Wapda Grid Colony

Sialkot, 51310 - Pakistan.

Tel : +92-52-3554782 - 3558704

Fax: +92-52-3550296

Cell: +92-333-8680640

E-mail: wmm@brain.net.pk

: welmed@cyber.net.pk

: info@welmedinst.com

Web : www.welmedinst.com

ISO 13485 / 9001

CE Mark

FDA USA (certified)