



Serving Humanity Since 1908

WELMED
INSTRUMENTS CO.

**TC-Gold Catalogue
Group 2, 3 & 4**



Manufacturers & Exporters of Quality Surgical Instruments

Tungsten Carbide Instruments



TC Forceps 03

WELMED Instruments Company
Defence Road
Near Wapda Grid Colony
Sialkot 51310, Pakistan.

Phone : +92 52 3554782, 3558704
Fax : +92 52 3550296
Cell : +92 333 8680640

Email : info@welmedinst.com
welmedpk@gmail.com
wmm@brain.net.pk
welmed@cyber.net.pk

Website : www.welmedinst.com

About Us

Quality begins at **WELMED** Instruments Company where master craftsmen forge and manufacture each instrument to its exact specifications and measurements. Instruments must undergo a series of quality control inspections in order to be worthy of **WELMED** name and warranty. We offer a wide selection of medical instruments including Surgical and Dental instruments.

We have acquired a unique and unmatched understanding of manufacturing through a large experience, working since 1908. Years of experience & state of the art manufacturing facilities helped us to keep our prices low. Our Surgical & Dental instruments are created by expert craftsmen, providing unfaltering quality you can rely on.

We have launched a fully interactive website where you can get online and track your orders 24 hours a day, 7 days in a week, and 365 days a year. Ask for our free e-catalogue on CD-Rom to save time and enjoy hassle free ordering.

We provide our customers with most reliable customer service as your satisfaction is our top priority. Our customers always enjoy our wonderful customer service and quality tools at the best price.

Lastly we invite you to have a work experience with **WELMED** Instruments Company and we will prove ourselves as one of the most valuable source you have.

You can find us at

WELMED Instruments Company
Defence Road, Near Wapda Grid Colony
Sialkot 51310, Pakistan.

Phone : +92 52 3554782, 3558704
Fax : +92 52 3550296
Cell : +92 333 8680640

Emails :
Info@welmedinst.com
welmedpk@gmail.com
wmm@brain.net.pk
welmed@cyber.net.pk

Website : www.welmedinst.com



History

WELMED Instruments Company was founded in 1995. Our Great Grandfather **Mian Muhammad Fazal Elahi** was the founder of M/s M.F. Elahi & Company (E.P.) Pvt. Ltd, established in 1908. WELMED is a sister concern of M.F. Elahi & Co. It is our pride to be a member of M.F. Elahi Family.

WELMED Instruments Company is located in the heart of beautiful city, Sialkot, Pakistan. We are manufacturers and exporters of quality Surgical, Dental instruments and Scissors of all sorts in various categories. We are closely attached with highly skilled surgeons and other technicians of the surgical industry to make the new instruments development possible.

WELMED is in a continuous process of development and testing of new instruments for their precision and perfect functioning. We are also working at imported forging supply by foreign customers.

Certificates

Welmed Instruments Company has achieved the certificates for ISO-9001, ISO-13485, cGMP (FDA) and CE Marking.

Quality Control

All material used in the manufacturing of our instruments are tested from relevant laboratories. Our continuous quality control is our best warranty so that you built confidence in our products.

Deliveries

Deliveries are quick and in-time.

Mode Of Payment

Terms of Payment is ADVANCE or LETTER OF CREDIT (Irrevocable)

Quality Policy

WELMED Instruments Company aims to produce quality Surgical, Dental instruments and Scissors of all sorts at reasonable prices to become an exemplary manufacturer. Our total emphasize is to achieve customer's satisfaction by clearly understanding their requirements and following a quality management system which is a primary concern of all our work force. We do our best to produce instruments of finest quality.



Warranty

All WELMED Instruments are produced from finest stainless steel available to be free from any defect in material or workmanship. Any product that proves to be defective in material or workmanship will either be repaired or replaced without any charge. This warranty is void if the instrument is abused, altered, modified or is not used in the manner for which it is intended.

TC Forceps

03




WELMED
INSTRUMENTS CO.


Adson Delicate	20
Adson Micro	20, 25
Adson Par Style	20
Adson	20,25,30
Adson-Brown	25
Baker Mod. (bonney)	28
Bonney	24
Cushing Micro	21
Cushing	21,26
Cushing-Potts	21
Delicate	27
Disectting Forceps	20,24
Evans	23
Ferris-Smith	28
Gerald Milcro	21
Gerald	21,26
Gillies	27
Leibinger	30
Leibinger-Debakey	29, 30
Leibinger-Williams	27
Mc Indoe	23
Obwegeser	23
Potts-Smith Micro	23
Potts-Smith	23
Rees Bayonet	24
Semkin Micro	26
Semkin	26
Skinner Kelicate	27
Skinner Micro	27
Skinner	27
Standard Delicate	22
Standard	22,28
Thorlakson	30
Tickner	24
Tissue Forceps	25,30
Uda Mod.	23
Wangensteen	24
Wheeler	30



- Wide selection of patterns for use in all types of surgical procedures
- safe and secure grasping of finest tissue
- precisely machined and hand-crafted for longlasting performance
- gold-plated ends for easy identification




Serrations: 0.5 mm



REGULAR

This tooth size is ideal for use in standard pattern needle holders like Mayo-Hegar, Mathieu, Bozeman a.o. And in needle-pulling tissue forceps. For holding large needles with sutures sized 3-0 and larger.


Serrations: 0.4 mm



FINE

Best for use in needle holders such as Crile-Wood, DeBakey, Ryder, a.o. And in delicate needle pulling tissue forceps. For holding needles and sutures between 4-0 and 6-0.


Serrations: 0.65 mm



COARSE

This coarsely serrated insert is placed in the most powerful needle holder patterns and wire twisters and assures maximum holding ability when using heavy needles or surgical wire.

Serrations: 0.2 mm



MICRO

This micro size is extremely effective for use in microsurgical needle holders like Castroviejo, Mini-Ryder, Microvascular, Lemmon and other microoriented needle holders. For use with the most delicate needles and sutures from size 7-0 to 10-0. Excellent delicate needle holding ability without minimal damage to needle or suture.



SMOOTH

The smooth insert is designed to hold the smaller suture sizes without damage to the suture or needle. Its needle holding ability is somewhat limited, however. Used primarily in needle holders for plastic surgery and microsurgery.



Adson



Adson Par Style




Adson Delicate




Adson



Adson Micro






1.8 mm

Regular
BC-03-107
 18 cm


Cushing-Potts



1.8 mm


Regular
BC-03-108
 18 cm

Cushing



Fine
BC-03-109
 18 cm


Cushing Micro



1.8 mm

Fine
BC-03-110
 18 cm
 With lock


Cushing



1.8 mm

Regular
BC-03-111
 18 cm

Cushing




Fine
BC-03-112
 17.5 cm

Gerald



Fine
BC-03-113
 17.5 cm

Gerald Micro





Regular
BC-03-122 15 cm

USA Mod.



Regular
BC-03-123 18 cm
BC-03-124 21 cm
BC-03-125 24.5 cm

Potts-smith



Fine
BC-03-126 21 cm

Potts-smith Micro



Regular
BC-03-127 25.5 cm

Evans



Fine
BC-03-128 15 cm

Mc Indoe



BC-03-129 16 cm
BC-03-130 18 cm

Obwegeser



Tickner



Regular
BC-03-131 14.5 cm
BC-03-132 19.5 cm



Bonney



Coarse
BC-03-133 17 cm



Wangenstein



Regular
BC-03-134 15.5 cm
BC-03-135 23 cm
BC-03-136 25.5 cm



Regular



BC-03-137 13.5 cm



Regular



BC-03-138 15.5 cm

Rees Bayonet



Regular



BC-03-139 19.5 cm
BC-03-140 23.5 cm

Rees Bayonet



Regular
BC-03-141 12 cm

Adson



3 x 4 Teeth



Regular
BC-03-142 12 cm

Adson-Brown



BC-03-143 12 cm

Adson Micro



Regular
BC-03-144 15 cm

Adson



3 x 4 Teeth



Regular
BC-03-145 15 cm

Adson-Brown



Fine
BC-03-146 15 cm



Regular

BC-03-147 14.5 cm
BC-03-148 16 cm

Semkin



Fine

BC-03-149 14.5 cm
BC-03-150 16 cm

Semkin Micro



Fine

BC-03-151 18 cm

Semkin Micro



Regular

BC-03-152 17.5 cm

Gerald



Regular

BC-03-153 17.5 cm

Gerald



Regular

BC-03-154 18 cm

Cushing





Regular
BC-03-155 15 cm

Gillies



Regular
BC-03-156 18 cm

Delicate



3 x 4 Teeth



Regular
BC-03-157 20 cm

Leibinaer-Williams



Regular
BC-03-158 23 cm

Skinner




Regular
BC-03-159 23 cm

Skinner Delicate




Regular
BC-03-160 23 cm

Skinner Micro




Coarse
BC-03-161 14.5 cm

Standard




Coarse
BC-03-162 16 cm

Standard




Coarse
BC-03-163 18 cm

Standard



Coarse
BC-03-164 17 cm

Baker Mod. (Bonney)











Coarse
BC-03-165 17 cm

Ferries-Smith



DEBAKEY-LEIBINGER Tissue Forceps with atraumatic serrations at the distal tip and finely cut Tungsten Carbide inserts placed directly behind.
 The serration pattern cut in to the tips will grasp delicate tissues firmly yet with minimal trauma. The additional Tungsten Carbide insert allow the use of these forceps for controlled needle "pulling" as the needle is drawn through the tissue. Particularly use full for delicate vascular and cardiac surgery.

Atraumatic Serrations

1.5 mm	2.0 mm	2.7 mm	3.5 mm
			
Fine	Fine	Fine	Fine
BC-03-166 15 cm BC-03-167 19.5 cm BC-03-168 24 cm	BC-03-169 15 cm BC-03-170 19.5 cm BC-03-171 24 cm	BC-03-172 15 cm BC-03-173 19.5 cm BC-03-174 24 cm	BC-03-175 19.5 cm BC-03-176 24 cm
			
Leibinger-DeBakey	Leibinger-DeBakey	Leibinger-DeBakey	Leibinger-DeBakey

Atraumatic Serrations
Handles are shapes
to allow comfortable
finger placement.



2.0 mm
Fine

Leibinger-DeBakey
BC-03-177 15 cm



2.0 mm
Fine

Leibinger-DeBakey
BC-03-178 19.5 cm



2.0 mm
Fine

Leibinger-DeBakey
BC-03-179 24 cm




Regular
BC-03-180 12.5 cm
BC-03-181 15 cm

Adson



Regular
BC-03-186 29.5 cm



Regular
BC-03-182 13.5 cm
BC-03-183 19.5 cm
BC-03-184 22.5 cm

Wheeler



Fine
BC-03-185 19.5 cm

Leibinger



Thorlakson



**WELMED Instrument Company
Defence Road, Near Wapda Grid Station
Sialkot 51310, Pakistan.**

Phone : +92 52 3554782, 3558704

Fax : +92 52 3550296

☎ Mobile : +92 333 8680640

Email : welmedpk@gmail.com

info@welmedinst.com

wmm@brain.net.pk

welmed@cyber.net.pk

Website : www.welmedinst.com