

VELMED

INSTRUMENTS CO.

Laparoscopy Instruments

Minimally Invasive Surgery



DEFENCE ROAD NEAR WAPDA GRID
STATION SIALKOT, PAKISTAN



(+92) 340 4040061, (+92) 3338680640



W-001-00
Manhes Style Ratchet




W-002-00
Hemostat Style Ratchet



W-003-00
Without Ratchet



W-004-00
Disengageable with Ratchet




W-005-00
Without Ratchet



W-006-00
Without Ratchet



W-007-00
Angulated ring handle With Ratchet



W-008-00
Disengageable With Ratchet

All are in Insulated Handles



W-011-00
Hemostat Style Ratchet.



W-012-00
Handle



W-013-00
Disengageable with
Ratchet.



W-014-00
without Ratchet.



W-015-00
Manhes Style Ratchet.

Dissecting & Grasping Forceps



Operating & Instrumental Information

Rotating, Dismantling,
With Connector Pin For Unipolar Coagulation, Size 2 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 20 cm,
to be use with trocars size 2.5 mm.

All are Double Action Jaws

Working Length	Dia.			
20 cm / 8"	2mm			
 Reddick-Olsen Dissecting and Grasping Forceps				
		W-001-20	W-002-20	W-003-20



Operating & Instrumental Information

Rotating, Dismantling,
Without Connector Pin For Unipolar Coagulation, Size 2 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 20 cm,
to be use with trocars size 3.5 mm.

All are Double Action Jaws

Working Length	Dia.				
20 cm / 8"	2mm				
 Atraumatic Dissecting and Grasping Forceps					
		W-004-20	W-005-20	W-006-20	W-007-20
 Kelly Dissecting and Grasping Forceps, Long.					
		W-008-20	W-009-20	W-010-20	W-011-20
 Kelly Dissecting and Grasping Forceps.					
		W-012-20	W-013-20	W-014-20	W-015-20
 Dissecting and Grasping Forceps, right angled.					
		W-016-20	W-017-20	W-018-20	W-019-20



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 20 cm, to be use with trocars size 3.5 mm.

Dissecting & Grasping Forceps

→ Operating & Instrumental Information

Rotating, Dismantling,
With Connector Pin For Unipolar Coagulation, Size 3 mm.

All are Double Action Jaws

Working Length	Dia.				
20 cm / 8"	2mm				



Reddick-Olsen
Dissecting and Grasping Forceps, Strong.

W-020-20

W-021-20

W-022-20

W-023-20



Grasping Forceps, Atraumatic, Fenestrated.

W-024-20

W-025-20

W-026-20

W-027-20




Alain-Grousseau
Pyloric Spreader.

W-028-20

W-029-20

W-030-20

W-031-20



Tan Pyloric Spreader and Grasping Forceps, in and outside serrated jaws.

W-032-20

W-033-20

W-034-20

W-035-20

All are Single / Double Action Jaws



Babcock Grasping Forceps, Atraumatic, Fenestrated.
Double action jaws.

W-036-20

W-037-20

W-038-20

W-039-20



Nuwayhid Grasping Forceps,
Double action jaws.

W-040-20

W-041-20

W-042-20

W-043-20



Grasping Forceps, with especially fine Atraumatic Serration Fenestrated, Single action jaws

W-044-20

W-045-20

W-046-20

W-047-20

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 20 cm, to be use with trocars size 3.5 mm.

Operating & Instrumental Information

Rotating, Dismantling,
With & Without Connector Pin For Unipolar Coagulation,
Size 3 mm.

All are Single / Double Action Jaws

Working Length	Dia.				
20 cm / 8"	3mm				



Grasping Forceps, Atraumatic, Wavy, Double Action Jaws

W-048-20

W-049-20

W-050-20

W-051-20



Manhes Grasping Forceps, jaws with multi teeth, width of jaws 4.8 mm, for Atraumatic and accurate grasping, single action jaws.

W-052-20

W-053-20

W-054-20

W-055-20



Manhes Grasping Forceps, Cobra- jaws 1 x 2 teeth, double action jaws.

W-056-20

W-057-20

W-058-20

W-059-20



Manhes Grasping Forceps, Tiger jaws 2 x 4 teeth, for grasping and liberation of solid organs, in particular during adhesiolyses, double action jaws.

W-060-20

W-061-20

W-062-20

W-063-20



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, to be use with trocars size 3.5 mm.

Operating & Instrumental Information

Rotating, Dismantling,
Without Connector Pin For Unipolar Coagulation,
Size 3 mm.

All are Single Action Jaws

Working Length	Dia.		
30 cm / 12"	3mm		



KOH Ultramicro I Grasping Forceps, straight. Fine Graspers with sub-millimeter jaws for handling 6-0 to 8-0 suture without cutting the suture.

W-064-30

W-065-30



KOH Ultramicro I Grasping Forceps, curved left. For ipsilateral suturing technique.

W-066-30

W-067-30

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, to be use with trocars size 3.5 mm.

→ **Operating & Instrumental Information**
Rotating, Dismantling,
Without Connector Pin For Unipolar Coagulation,
Size 3 mm.

All are Single Action Jaws

Working Length	Dia.		
30 cm / 12"	3mm		



KOH Ultramicro I Grasping Forceps, curved right. For ipsilateral suturing technique.

W-068-30

W-069-30



KOH Ultramicro II Grasping Forceps.

W-070-30

W-071-30



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, to be use with trocars size 3.9 mm.

→ **Operating & Instrumental Information**
Rotating, Dismantling,
Without Connector Pin For Unipolar Coagulation,
Size 3 mm.

All are Single / Double Action Jaws

Working Length	Dia.				
30 cm / 12"	3.5mm				




Kelly Dissecting and Grasping Forceps.

W-072-30

W-073-30

W-074-30

W-075-30



Reddick-Olsen Dissecting and Grasping Forceps, Strong.

W-076-30

W-077-30

W-078-30

W-079-30



Grasping Forceps, Atraumatic, Wavy.

W-080-30

W-081-30

W-082-30

W-083-30

Dissecting & Grasping Forceps



Operating & Instrumental Information
Rotating, Dismantling,
Without Connector Pin For Unipolar Coagulation,
Size 3.5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
to be use with trocars size 3.9 mm.

All are Single / Double Action Jaws

Working Length	Dia.				
30 cm / 12"	3.5mm				
 Grasping Forceps , with especially fine Atraumatic Serration, Fenestrated, Single Action Jaws.					
W-084-30 W-085-30 W-086-30 W-087-30					
 Babcock Grasping Forceps , Atraumatic, Fenestrated, Double Action Jaws.					
W-088-30 W-089-30 W-090-30 W-091-30					
 Grasping Forceps , Fenestrated, Double Action Jaws.					
W-092-30 W-093-30 W-094-30 W-095-30					

Dissecting & Grasping Forceps



Operating & Instrumental Information

Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 3.5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length					
30 cm / 12" 36 cm / 14" 43 cm / 17"					



Manhes Dissecting and Grasping Forceps, atraumatic, specially designed for grasping and liberation of hollow organs.

W-096-30	W-097-30	W-098-30	W-099-30	W100-30
W-096-36	W-097-36	W-098-36	W-099-36	W100-36
W-096-43	W-097-43	W-098-43	W-099-43	W100-43



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, blunt, specially designed for extended clinching of tissue thus causing mechanical hemostasis length of jaws 15mm.

W-101-30	W-102-30	W-103-30	W-104-30	W105-30
W-101-36	W-102-36	W-103-36	W-104-36	W105-36
W-101-43	W-102-43	W-103-43	W-104-43	W105-43



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, blunt, specially designed for extended clinching of tissue thus causing mechanical hemostasis length of jaws 10mm.

W-106-30	W-107-30	W-108-30	W-109-30	W110-30
W-106-36	W-107-36	W-108-36	W-109-36	W110-36
W-106-43	W-107-43	W-108-43	W-109-43	W110-43



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, blunt, specially designed for partial clinching of tissue thus causing selective hemostasis length of jaws 15mm.

W-111-30	W-112-30	W-113-30	W-114-30	W115-30
W-111-36	W-112-36	W-113-36	W-114-36	W115-36
W-111-43	W-112-43	W-113-43	W-114-43	W115-43

Dissecting & Grasping Forceps



→ **Operating & Instrumental Information** ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 3.5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Manhes Dissecting and Grasping Forceps, atraumatic, specially designed for grasping and liberation of hollow organs.

W-116-30 W-116-36 W-116-43	W-117-30 W-117-36 W-117-43	W-118-30 W-118-36 W-118-43	W-119-30 W-119-36 W-119-43	W-120-30 W-120-36 W-120-43	W-121-30 W-121-36 W-121-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, blunt, specially designed for extended clinching of tissue thus causing mechanical hemostasis length of jaws 15mm.

W-122-30 W-122-36 W-122-43	W-123-30 W-123-36 W-123-43	W-124-30 W-124-36 W-124-43	W-125-30 W-125-36 W-125-43	W-126-30 W-126-36 W-126-43	W-127-30 W-127-36 W-127-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, blunt, specially designed for extended clinching of tissue thus causing mechanical hemostasis length of jaws 10mm.

W-128-30 W-128-36 W-128-43	W-129-30 W-129-36 W-129-43	W-130-30 W-130-36 W-130-43	W-131-30 W-131-36 W-131-43	W-132-30 W-132-36 W-132-43	W-133-30 W-133-36 W-133-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, blunt, specially designed for partial clinching of tissue thus causing selective hemostasis length of jaws 15mm.

W-134-30 W-134-36 W-134-43	W-135-30 W-135-36 W-135-43	W-136-30 W-136-36 W-136-43	W-137-30 W-137-36 W-137-43	W-138-30 W-138-36 W-138-43	W-139-30 W-139-36 W-139-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



→ **Operating & Instrumental Information** ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Manhes Dissecting and Grasping Forceps, atraumatic, specially designed for grasping and liberation of hollow organs.

W-140-30 W-140-36 W-140-43	W-141-30 W-141-36 W-141-43	W-142-30 W-142-36 W-142-43	W-143-30 W-143-36 W-143-43	W-144-30 W-144-36 W-144-43	W-145-30 W-145-36 W-145-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, bulnt, specially designed for extebded clinching of tissue thus causing mechanical hemostasis lenght of jaws 15mm.

W-146-30 W-146-36 W-146-43	W-147-30 W-147-36 W-147-43	W-148-30 W-148-36 W-148-43	W-149-30 W-149-36 W-149-43	W-150-30 W-150-36 W-150-43	W-151-30 W-151-36 W-151-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, bulnt, specially designed for extebded clinching of tissue thus causing mechanical hemostasis lenght of jaws 10mm.

W-152-30 W-152-36 W-152-43	W-153-30 W-153-36 W-153-43	W-154-30 W-154-36 W-154-43	W-155-30 W-155-36 W-155-43	W-156-30 W-156-36 W-156-43	W-157-30 W-157-36 W-157-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Dissecting and Grasping Forceps, Duck-bill jaws, bulnt, specially designed for partial clinching of tissue thus causing selective hemostasis lenght of jaws 15mm.

W-158-30 W-158-36 W-158-43	W-159-30 W-159-36 W-159-43	W-160-30 W-160-36 W-160-43	W-161-30 W-161-36 W-161-43	W-162-30 W-162-36 W-162-43	W-163-30 W-163-36 W-163-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



Operating & Instrumental Information

Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length					
30 cm / 12"					
36 cm / 14"					
43 cm / 17"					



Manhes Dissecting and Grasping Forceps, concical, blunt, specially designed for extended
clishing of tissue thus causing mechanical hemostasis, lenght of jaws 15mm.

W-164-30	W-165-30	W-166-30	W-167-30	W-168-30
W-164-36	W-165-36	W-166-36	W-167-36	W-168-36
W-164-43	W-165-43	W-166-43	W-167-43	W-168-43



Manhes Dissecting and Grasping Forceps, concical, blunt, specially designed for extended
clishing of tissue thus causing mechanical hemostasis, lenght of jaws 10mm.

W-169-30	W-170-30	W-171-30	W-172-30	W-173-30
W-169-36	W-170-36	W-171-36	W-172-36	W-173-36
W-169-43	W-170-43	W-171-43	W-172-43	W-173-43



Manhes Dissecting and Grasping Forceps, concical, blunt, specially designed for partial
clishing of tissue thus causing selective hemostasis, lenght of jaws 15mm.

W-174-30	W-175-30	W-176-30	W-177-30	W-178-30
W-174-36	W-175-36	W-176-36	W-177-36	W-178-36
W-174-43	W-175-43	W-176-43	W-177-43	W-178-43



Dissecting and Grasping Forceps, Dolphin Nose, small.

W-179-30	W-180-30	W-181-30	W-182-30	W-183-30
W-179-36	W-180-36	W-181-36	W-182-36	W-183-36
W-179-43	W-180-43	W-181-43	W-182-43	W-183-43

Dissecting & Grasping Forceps



► Operating & Instrumental Information ◀
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Manhes Dissecting and Grasping Forceps, concical, blunt, specially designed for extended
clipping of tissue thus causing mechanical hemostasis, lenght of jaws 15mm.

W-184-30	W-185-30	W-186-30	W-187-30	W-188-30	W-189-30
W-184-36	W-185-36	W-186-36	W-187-36	W-188-36	W-189-36
W-184-43	W-185-43	W-186-43	W-187-43	W-188-43	W-189-43



Manhes Dissecting and Grasping Forceps, concical, blunt, specially designed for extended
clipping of tissue thus causing mechanical hemostasis, lenght of jaws 10mm.

W-190-30	W-191-30	W-192-30	W-193-30	W-194-30	W-195-30
W-190-36	W-191-36	W-192-36	W-193-36	W-194-36	W-195-36
W-190-43	W-191-43	W-192-43	W-193-43	W-194-43	W-195-43



Manhes Dissecting and Grasping Forceps, concical, blunt, specially designed for partial
clipping of tissue thus causing selective hemostasis, lenght of jaws 15mm.

W-196-30	W-197-30	W-198-30	W-199-30	W-200-30	W-201-30
W-196-36	W-197-36	W-198-36	W-199-36	W-200-36	W-201-36
W-196-43	W-197-43	W-198-43	W-199-43	W-200-43	W-201-43



Dissecting and Grasping Forceps, Dolphin Nose, small.

W-202-30	W-203-30	W-204-30	W-205-30	W-206-30	W-207-30
W-202-36	W-203-36	W-204-36	W-205-36	W-206-36	W-207-36
W-202-43	W-203-43	W-204-43	W-205-43	W-206-43	W-207-43

Dissecting & Grasping Forceps



► Operating & Instrumental Information ◀
Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length					
30 cm / 12"					
36 cm / 14"					
43 cm / 17"					



Dissecting and Grasping Forceps, atraumatic.

W-208-30	W-209-30	W-210-30	W-211-30	W-212-30
W-208-36	W-209-36	W-210-36	W-211-36	W-212-36
W-208-43	W-209-43	W-210-43	W-211-43	W-212-43




Dissecting and Grasping Forceps, atraumatic.

W-213-30	W-214-30	W-215-30	W-216-30	W-217-30
W-213-36	W-214-36	W-215-36	W-216-36	W-217-36
W-213-43	W-214-43	W-215-43	W-216-43	W-217-43



Dissecting and Grasping Forceps, "Dolphin Nose".

W-218-30	W-219-30	W-220-30	W-221-30	W-222-30
W-218-36	W-219-36	W-220-36	W-221-36	W-222-36
W-218-43	W-219-43	W-220-43	W-221-43	W-222-43



Reddick-Olsen Dissecting and Grasping Forceps.

W-223-30	W-224-30	W-225-30	W-226-30	W-227-30
W-223-36	W-224-36	W-225-36	W-226-36	W-227-36
W-223-43	W-224-43	W-225-43	W-226-43	W-227-43

Dissecting & Grasping Forceps




► Operating & Instrumental Information ◀
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.


All are Double Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Reddick-Olsen Dissecting and Grasping Forceps, heavy.

W-228-30	W-229-30	W-230-30	W-231-30	W-232-30	W-233-30
W-228-36	W-229-36	W-230-36	W-231-36	W-232-36	W-233-36
W-228-43	W-229-43	W-230-43	W-231-43	W-232-43	W-233-43



Dissecting and Grasping Forceps, alligator jaws.

W-234-30	W-235-30	W-236-30	W-237-30	W-238-30	W-239-30
W-234-36	W-235-36	W-236-36	W-237-36	W-238-36	W-239-36
W-234-43	W-235-43	W-236-43	W-237-43	W-238-43	W-239-43



Kelly Dissecting and Grasping Forceps.

W-240-30	W-241-30	W-242-30	W-243-30	W-244-30	W-245-30
W-240-36	W-241-36	W-242-36	W-243-36	W-244-36	W-245-36
W-240-43	W-241-43	W-242-43	W-243-43	W-244-43	W-245-43



Kelly Dissecting and Grasping Forceps, long.

W-246-30	W-247-30	W-248-30	W-249-30	W-250-30	W-251-30
W-246-36	W-247-36	W-248-36	W-249-36	W-250-36	W-251-36
W-246-43	W-247-43	W-248-43	W-249-43	W-250-43	W-251-43

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, 36 cm, to be use with trocars size 6 mm.

All are Double Action Jaws

Operating & Instrumental Information

Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Single Puncture Approach.
Operating instruments for single Puncture approach length, 43 cm, to be use with operating laparoscopes with built-in instrument channel.

Working Length					
30 cm / 12" 36 cm / 14" 43 cm / 17"					



Dissecting and Grasping Forceps, right angled.

W-252-30 W-252-36 W-252-43	W-253-30 W-253-36 W-253-43	W-254-30 W-254-36 W-254-43	W-255-30 W-255-36 W-255-43	W-256-30 W-256-36 W-256-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Dissecting and Grasping Forceps, atraumatic, with hollow jaws.

W-257-30 W-257-36 W-257-43	W-258-30 W-258-36 W-258-43	W-259-30 W-259-36 W-259-43	W-260-30 W-260-36 W-260-43	W-261-30 W-261-36 W-261-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, 36 cm, to be use with trocars size 6 mm.

All are Double Action Jaws

Operating & Instrumental Information

Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Single Puncture Approach.
Operating instruments for single Puncture approach length, 43 cm, to be use with operating laparoscopes with built-in instrument channel.

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Dissecting and Grasping Forceps, right angled.

W-262-30 W-262-36 W-262-43	W-263-30 W-263-36 W-263-43	W-364-30 W-264-36 W-264-43	W-265-30 W-265-36 W-265-43	W-266-30 W-266-36 W-266-43	W-267-30 W-267-36 W-267-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Dissecting and Grasping Forceps, atraumatic, with hollow jaws.

W-268-30 W-268-36 W-268-43	W-269-30 W-269-36 W-269-43	W-270-30 W-270-36 W-270-43	W-271-30 W-271-36 W-271-43	W-272-30 W-272-36 W-272-43	W-273-30 W-273-36 W-273-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



Operating & Instrumental Information

Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length					
30 cm / 12"					
36 cm / 14"					
43 cm / 17"					



Dissecting and Grasping Forceps, with especially fine atraumatic serration, fenestrated.

W-274-30	W-275-30	W-276-30	W-277-30	W-278-30
W-274-36	W-275-36	W-276-36	W-277-36	W-278-36
W-274-43	W-275-43	W-276-43	W-277-43	W-278-43



Dissecting and Grasping Forceps, atraumatic, with hollow jaws.

W-279-30	W-280-30	W-281-30	W-282-30	W-283-30
W-279-36	W-280-36	W-281-36	W-282-36	W-283-36
W-279-43	W-280-43	W-281-43	W-282-43	W-283-43



Croce-Olmi Dissecting and Grasping Forceps, atraumatic, fenestrated, curved.

W-284-30	W-285-30	W-286-30	W-287-30	W-288-30
W-284-36	W-285-36	W-286-36	W-287-36	W-288-36
W-284-43	W-285-43	W-286-43	W-287-43	W-288-43

Dissecting & Grasping Forceps



► Operating & Instrumental Information ◀
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Dissecting and Grasping Forceps, with especially fine atraumatic serration, fenestrated.

W-289-30	W-290-30	W-391-30	W-292-30	W-293-30	W-294-30
W-289-36	W-290-36	W-291-36	W-292-36	W-293-36	W-294-36
W-289-43	W-290-43	W-291-43	W-292-43	W-293-43	W-294-43



Dissecting and Grasping Forceps, atraumatic, with hollow jaws.

W-295-30	W-296-30	W-297-30	W-298-30	W-299-30	W-300-30
W-295-36	W-296-36	W-297-36	W-298-36	W-299-36	W-300-36
W-295-43	W-296-43	W-297-43	W-298-43	W-299-43	W-300-43



Croce-Olmi Dissecting and Grasping Forceps, atraumatic, fenestrated, curved.

W-301-30	W-302-30	W-303-30	W-304-30	W-305-30	W-306-30
W-301-36	W-302-36	W-303-36	W-304-36	W-305-36	W-306-36
W-301-43	W-302-43	W-303-43	W-304-43	W-305-43	W-306-43

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length					
30 cm / 12" 36 cm / 14" 43 cm / 17"					



Schneider Lymph Node Dissecting and Grasping Forceps, atraumatic.

W-307-30 W-307-36 W-307-43	W-308-30 W-308-36 W-308-43	W-308-30 W-308-36 W-308-43	W-308-30 W-308-36 W-308-43	W-309-30 W-309-36 W-309-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



DeBaey Dissecting and Grasping Forceps, atraumatic.

W-310-30 W-310-36 W-310-43	W-311-30 W-311-36 W-311-43	W-312-30 W-312-36 W-312-43	W-313-30 W-313-36 W-313-43	W-314-30 W-314-36 W-314-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Grasping Forceps, atraumatic.

W-315-30 W-315-36 W-315-43	W-316-30 W-316-36 W-316-43	W-317-30 W-317-36 W-317-43	W-318-30 W-318-36 W-318-43	W-319-30 W-319-36 W-319-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Babcock Grasping Forceps, atraumatic, jaws with multiple teeth, fenestrated.

W-320-30 W-320-36 W-320-43	W-321-30 W-321-36 W-321-43	W-322-30 W-322-36 W-322-43	W-323-30 W-323-36 W-323-43	W-324-30 W-324-36 W-324-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



► Operating & Instrumental Information ◀
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Schneider Lymph Node Dissecting and Grasping Forceps, atraumatic.

W-325-30 W-325-36 W-325-43	W-326-30 W-326-36 W-326-43	W-327-30 W-327-36 W-327-43	W-328-30 W-328-36 W-328-43	W-329-30 W-329-36 W-329-43	W-330-30 W-330-36 W-330-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



DeBaKey Dissecting and Grasping Forceps, atraumatic.

W-331-30 W-331-36 W-331-43	W-332-30 W-332-36 W-332-43	W-333-30 W-333-36 W-333-43	W-334-30 W-334-36 W-334-43	W-335-30 W-335-36 W-335-43	W-336-30 W-336-36 W-336-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Grasping Forceps, atraumatic.

W-337-30 W-337-36 W-337-43	W-338-30 W-338-36 W-338-43	W-339-30 W-339-36 W-339-43	W-340-30 W-340-36 W-340-43	W-341-30 W-341-36 W-341-43	W-342-30 W-342-36 W-342-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Babcock Grasping Forceps, atraumatic, jaws with multiple teeth, fenestrated.

W-343-30 W-343-36 W-343-43	W-344-30 W-344-36 W-344-43	W-345-30 W-345-36 W-345-43	W-346-30 W-346-36 W-346-43	W-347-30 W-347-36 W-347-43	W-348-30 W-348-36 W-348-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length					
30 cm / 12"					
36 cm / 14"					
43 cm / 17"					



Babcock Grasping Forceps, atraumatic, jaws with multiple teeth, fenestrated, long.

W-349-30	W-350-30	W-351-30	W-352-30	W-353-30
W-349-36	W-350-36	W-351-36	W-352-36	W-353-36
W-349-43	W-350-43	W-351-43	W-352-43	W-353-43



Manhes Grasping Forceps, "Cobra-Jaws", 1 x 2 teeth.

W-354-30	W-355-30	W-356-30	W-357-30	W-358-30
W-354-36	W-355-36	W-356-36	W-357-36	W-358-36
W-354-43	W-355-43	W-356-43	W-357-43	W-358-43



Manhes Grasping Forceps, jaws with multiple teeth, width of jaws 4.8 mm, for atraumatic and accurate grasping.

W-359-30	W-360-30	W-361-30	W-362-30	W-363-30
W-359-36	W-360-36	W-361-36	W-362-36	W-363-36
W-359-43	W-360-43	W-361-43	W-362-43	W-363-43



Manhes Grasping Forceps, "Tiger-Jaws", 2 x 4 teeth, for grasping and liberation of solid organs, in particular during adhesiolyses.

W-364-30	W-365-30	W-366-30	W-367-30	W-368-30
W-364-36	W-365-36	W-366-36	W-367-36	W-368-36
W-364-43	W-365-43	W-366-43	W-367-43	W-368-43

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Babcock Grasping Forceps, atraumatic, jaws with multiple teeth, fenestrated, long.

W-369-30	W-370-30	W-371-30	W-372-30	W-373-30	W-374-30
W-369-36	W-370-36	W-371-36	W-372-36	W-373-36	W-374-36
W-369-43	W-370-43	W-371-43	W-372-43	W-373-43	W-374-43



Manhes Grasping Forceps, "Cobra-Jaws", 1 x 2 teeth.

W-375-30	W-376-30	W-377-30	W-378-30	W-379-30	W-380-30
W-375-36	W-376-36	W-377-36	W-378-36	W-379-36	W-380-36
W-375-43	W-376-43	W-377-43	W-378-43	W-379-43	W-380-43



Manhes Grasping Forceps, jaws with multiple teeth, width of jaws 4.8 mm, for atraumatic and accurate grasping.

W-381-30	W-382-30	W-383-30	W-384-30	W-385-30	W-386-30
W-381-36	W-382-36	W-383-36	W-384-36	W-385-36	W-386-36
W-381-43	W-382-43	W-383-43	W-384-43	W-385-43	W-386-43



Manhes Grasping Forceps, "Tiger-Jaws", 2 x 4 teeth, for grasping and liberation of solid organs, in particular during adhesiolyses.

W-387-30	W-388-30	W-389-30	W-390-30	W-391-30	W-392-30
W-387-36	W-388-36	W-389-36	W-390-36	W-391-36	W-392-36
W-387-43	W-388-43	W-389-43	W-390-43	W-391-43	W-392-43

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, 36 cm, to be use with trocars size 6 mm.

All are Double Action Jaws

Operating & Instrumental Information

Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Single Puncture Approach.
Operating instruments for single Puncture approach length, 43 cm, to be use with operating laparoscopes with built-in instrument channel.

Working Length					
30 cm / 12" 36 cm / 14" 43 cm / 17"					



Manhes Grasping Forceps, jaws with multiple teeth, width of jaws 3 mm, for atraumatic and accurate grasping.

W-393-30 W-393-36 W-393-43	W-394-30 W-394-36 W-394-43	W-395-30 W-395-36 W-395-43	W-396-30 W-396-36 W-396-43	W-397-30 W-397-36 W-397-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Claw Forceps, 2 x 3 teeth.

W-398-30 W-398-36 W-398-43	W-399-30 W-399-36 W-399-43	W-400-30 W-400-36 W-400-43	W-401-30 W-401-36 W-401-43	W-402-30 W-402-36 W-402-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, 36 cm, to be use with trocars size 6 mm.

All are Double Action Jaws

Operating & Instrumental Information

Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Single Puncture Approach.
Operating instruments for single Puncture approach length, 43 cm, to be use with operating laparoscopes with built-in instrument channel.

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Manhes Grasping Forceps, jaws with multiple teeth, width of jaws 3 mm, for atraumatic and accurate grasping.

W-403-30 W-403-36 W-403-43	W-404-30 W-404-36 W-404-43	W-405-30 W-405-36 W-405-43	W-406-30 W-406-36 W-406-43	W-407-30 W-407-36 W-407-43	W-408-30 W-408-36 W-408-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Claw Forceps, 2 x 3 teeth.

W-409-30 W-409-36 W-409-43	W-410-30 W-410-36 W-410-43	W-411-30 W-411-36 W-411-43	W-412-30 W-412-36 W-412-43	W-413-30 W-413-36 W-413-43	W-414-30 W-414-36 W-414-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length					
30 cm / 12"					
36 cm / 14"					
43 cm / 17"					



Grasping Forceps, atraumatic, wavy.

W-415-30	W-416-30	W-417-30	W-418-30	W-419-30
W-415-36	W-416-36	W-417-36	W-418-36	W-419-36
W-415-43	W-416-43	W-417-43	W-418-43	W-419-43



Grasping Forceps, atraumatic.

W-420-30	W-421-30	W-422-30	W-423-30	W-424-30
W-420-36	W-421-36	W-422-36	W-423-36	W-424-36
W-420-43	W-421-43	W-422-43	W-423-43	W-424-43



Grasping Forceps, atraumatic, double spoon, with multiple teeth.

W-425-30	W-426-30	W-427-30	W-428-30	W-429-30
W-425-36	W-426-36	W-427-36	W-428-36	W-429-36
W-425-43	W-426-43	W-427-43	W-428-43	W-429-43



Grasping Forceps, fenestrated.

W-430-30	W-431-30	W-432-30	W-433-30	W-434-30
W-430-36	W-431-36	W-432-36	W-433-36	W-434-36
W-430-43	W-431-43	W-432-43	W-433-43	W-434-43

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Grasping Forceps, atraumatic, wavy.

W-435-30	W-436-30	W-437-30	W-438-30	W-439-30	W-440-30
W-435-36	W-436-36	W-437-36	W-438-36	W-439-36	W-440-36
W-435-43	W-436-43	W-437-43	W-438-43	W-439-43	W-440-43



Grasping Forceps, atraumatic.

W-441-30	W-442-30	W-443-30	W-444-30	W-445-30	W-446-30
W-441-36	W-442-36	W-443-36	W-444-36	W-445-36	W-446-36
W-441-43	W-442-43	W-443-43	W-444-43	W-445-43	W-446-43



Grasping Forceps, atraumatic, double spoon, with multiple teeth.

W-447-30	W-448-30	W-449-30	W-450-30	W-541-30	W-452-30
W-447-36	W-448-36	W-449-36	W-450-36	W-451-36	W-452-36
W-447-43	W-448-43	W-449-43	W-450-43	W-451-43	W-452-43



Grasping Forceps, fenestrated.

W-453-30	W-454-30	W-455-30	W-456-30	W-457-30	W-458-30
W-453-36	W-454-36	W-455-36	W-456-36	W-457-36	W-458-36
W-453-43	W-454-43	W-455-43	W-456-43	W-457-43	W-458-43

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Babcock Grasping Forceps.

W-459-30	W-460-30	W-461-30	W-462-30	W-463-30	W-464-30
W-459-36	W-460-36	W-461-36	W-462-36	W-463-36	W-464-36
W-459-43	W-460-43	W-461-43	W-462-43	W-463-43	W-464-43



Bowel Grasper, fenestrated.

W-465-30	W-466-30	W-467-30	W-468-30	W-469-30	W-470-30
W-465-36	W-466-36	W-467-36	W-468-36	W-469-36	W-470-36
W-465-43	W-466-43	W-467-43	W-468-43	W-469-43	W-470-43



Mangeshikar Grasping Forceps, serrated.

W-471-30	W-472-30	W-473-30	W-474-30	W-575-30	W-476-30
W-471-36	W-472-36	W-473-36	W-474-36	W-475-36	W-476-36
W-471-43	W-472-43	W-473-43	W-474-43	W-475-43	W-476-43



Mangeshikar Grasping Forceps, without serrated.

W-477-30	W-478-30	W-479-30	W-480-30	W-481-30	W-482-30
W-477-36	W-478-36	W-479-36	W-480-36	W-481-36	W-482-36
W-477-43	W-478-43	W-479-43	W-480-43	W-481-43	W-482-43

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length					
30 cm / 12"					
36 cm / 14"					
43 cm / 17"					



Babcock Grasping Forceps.

W-483-30	W-484-30	W-485-30	W-486-30	W-487-30
W-483-36	W-484-36	W-485-36	W-486-36	W-487-36
W-483-43	W-484-43	W-485-43	W-486-43	W-487-43



Bowel Grasper, fenestrated.

W-488-30	W-489-30	W-490-30	W-491-30	W-492-30
W-488-36	W-489-36	W-490-36	W-491-36	W-492-36
W-488-43	W-489-43	W-490-43	W-491-43	W-492-43



Mangeshikar Grasping Forceps, serrated.

W-493-30	W-494-30	W-495-30	W-496-30	W-497-30
W-493-36	W-494-36	W-495-36	W-496-36	W-497-36
W-493-43	W-494-43	W-495-43	W-496-43	W-497-43



Mangeshikar Grasping Forceps, without serrated.

W-498-30	W-499-30	W-500-30	W-501-30	W-502-30
W-498-36	W-499-36	W-500-36	W-501-36	W-502-36
W-498-43	W-499-43	W-500-43	W-501-43	W-502-43

Dissecting & Grasping Forceps



Operating & Instrumental Information

Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length					
30 cm / 12" 36 cm / 14" 43 cm / 17"					



Vancaillie Adhesion Forceps, jaw fenestrated.

W-503-30 W-503-36 W-503-43	W-504-30 W-504-36 W-504-43	W-505-30 W-505-36 W-505-43	W-506-30 W-506-36 W-506-43	W-507-30 W-507-36 W-507-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Vancaillie Adhesion Forceps.

W-508-30 W-508-36 W-508-43	W-509-30 W-509-36 W-509-43	W-510-30 W-510-36 W-510-43	W-511-30 W-511-36 W-511-43	W-512-30 W-512-36 W-512-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Vancaillie Meso-ovarian Forceps.

W-513-30 W-513-36 W-513-43	W-514-30 W-514-36 W-514-43	W-515-30 W-515-36 W-515-43	W-516-30 W-516-36 W-516-43	W-517-30 W-517-36 W-517-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Vancaillie Adhesion Forceps, jaw fenestrated.

W-518-30	W-519-30	W-520-30	W-521-30	W-522-30	W-523-30
W-518-36	W-519-36	W-520-36	W-521-36	W-522-36	W-523-36
W-518-43	W-519-43	W-520-43	W-521-43	W-522-43	W-523-43



Vancaillie Adhesion Forceps.

W-524-30	W-525-30	W-526-30	W-527-30	W-528-30	W-529-30
W-524-36	W-525-36	W-526-36	W-527-36	W-528-36	W-529-36
W-524-43	W-525-43	W-526-43	W-527-43	W-528-43	W-529-43



Vancaillie Meso-ovarian Forceps.

W-530-30	W-531-30	W-532-30	W-533-30	W-534-30	W-535-30
W-530-36	W-531-36	W-532-36	W-533-36	W-534-36	W-535-36
W-530-43	W-531-43	W-532-43	W-533-43	W-534-43	W-535-43

Dissecting & Grasping Forceps



Operating & Instrumental Information

Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length					
30 cm / 12" 36 cm / 14" 43 cm / 17"					



Swolin Grasping Forceps, atraumatic bluntly polished edges and
teeth prevent from trauma and allow for gently holding the fallopian tubes or
similarly delicate structures.

W-536-30 W-536-36 W-536-43	W-537-30 W-537-36 W-537-43	W-538-30 W-538-36 W-538-43	W-539-30 W-539-36 W-539-43	W-540-30 W-540-36 W-540-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Matkowitz Grasping Forceps.

W-541-30 W-541-36 W-541-43	W-542-30 W-542-36 W-542-43	W-543-30 W-543-36 W-543-43	W-544-30 W-544-36 W-544-43	W-545-30 W-545-36 W-545-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Grasping Forceps, 2 x 4 teeth.

W-546-30 W-546-36 W-546-43	W-547-30 W-547-36 W-547-43	W-548-30 W-548-36 W-548-43	W-549-30 W-549-36 W-549-43	W-550-30 W-550-36 W-550-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Berci Multifunctional Grasping Forceps, 1 x 2 teeth,
to close gallbladder perforations or to facilitate blunt dissections.

W-551-30 W-551-36 W-551-43	W-552-30 W-552-36 W-552-43	W-553-30 W-553-36 W-553-43	W-554-30 W-554-36 W-554-43	W-555-30 W-555-36 W-555-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Swolin Grasping Forceps, atraumatic bluntly polished edges and
teeth prevent from trauma and allow for gently holding the fallopian tubes or
similarly delicate structures.

W-556-30	W-557-30	W-558-30	W-559-30	W-560-30	W-561-30
W-556-36	W-557-36	W-558-36	W-559-36	W-560-36	W-561-36
W-556-43	W-557-43	W-558-43	W-559-43	W-560-43	W-561-43



Matkowitz Grasping Forceps.

W-562-30	W-563-30	W-564-30	W-565-30	W-566-30	W-567-30
W-562-36	W-563-36	W-564-36	W-565-36	W-566-36	W-567-36
W-562-43	W-563-43	W-564-43	W-565-43	W-566-43	W-567-43



Grasping Forceps, 2 x 4 teeth.

W-568-30	W-569-30	W-570-30	W-571-30	W-572-30	W-573-30
W-568-36	W-569-36	W-570-36	W-571-36	W-572-36	W-573-36
W-568-43	W-569-43	W-570-43	W-571-43	W-572-43	W-573-43



Berci Multifunctional Grasping Forceps, 1 x 1 teeth,
to close gallbladder perforations or to facilitate blunt dissections.

W-574-30	W-575-30	W-576-30	W-577-30	W-578-30	W-579-30
W-574-36	W-575-36	W-576-36	W-577-36	W-578-36	W-579-36
W-574-43	W-575-43	W-576-43	W-577-43	W-578-43	W-579-43

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length					
30 cm / 12" 36 cm / 14" 43 cm / 17"					



Remorgida Grasping Forceps, 1 x 2 teeth.

W-580-30 W-580-36 W-580-43	W-581-30 W-581-36 W-581-43	W-582-30 W-582-36 W-582-43	W-583-30 W-583-36 W-583-43	W-584-30 W-584-36 W-584-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Remorgida Grasping Forceps, 1 x 2 teeth, insulated jaws for reduced
lateral thermal diffusion compared o standard full metal jaws.

W-585-30 W-585-36 W-585-43	W-586-30 W-586-36 W-586-43	W-587-30 W-587-36 W-587-43	W-588-30 W-588-36 W-588-43	W-589-30 W-589-36 W-589-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Shawki Cutting-Coagulation Forceps.

W-590-30 W-590-36 W-590-43			W-591-30 W-591-36 W-591-43	W-592-30 W-592-36 W-592-43
----------------------------------	--	--	----------------------------------	----------------------------------

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Remorgida Grasping Forceps, 1 x 2 teeth.

W-593-30	W-594-30	W-595-30	W-596-30	W-597-30	W-598-30
W-593-36	W-594-36	W-595-36	W-596-36	W-597-36	W-598-36
W-593-43	W-594-43	W-595-43	W-596-43	W-597-43	W-598-43



Grasping Forceps.

W-599-30	W-600-30	W-601-30	W-602-30	W-603-30	W-604-30
W-599-36	W-600-36	W-601-36	W-602-36	W-603-36	W-604-36
W-599-43	W-600-43	W-601-43	W-602-43	W-603-43	W-604-43

Dissecting & Grasping Forceps



Operating & Instrumental Information

Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length					
30 cm / 12" 36 cm / 14" 43 cm / 17"					



Manhes Micro Grasping Forceps, with multiple teeth. With a very fine uptake at the end of the jaws, the claws of which go no further than the serous membranes.

W-605-30 W-605-36 W-605-43	W-606-30 W-606-36 W-606-43	W-607-30 W-607-36 W-607-43	W-608-30 W-608-36 W-608-43	W-609-30 W-609-36 W-609-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Micro Dissecting and Grasping Forceps, atraumatic, full grasping

W-610-30 W-610-36 W-610-43	W-611-30 W-611-36 W-611-43	W-612-30 W-612-36 W-612-43	W-613-30 W-613-36 W-613-43	W-614-30 W-614-36 W-614-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Micro Dissecting and Grasping Forceps, atraumatic, short jaws, partial grasping.

W-615-30 W-615-36 W-615-43	W-616-30 W-616-36 W-616-43	W-617-30 W-617-36 W-617-43	W-618-30 W-618-36 W-618-43	W-619-30 W-619-36 W-619-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------



Manhes Micro Dissecting and Grasping Forceps, atraumatic, long jaws, partial grasping. Especially suitable for tying intra-corporeal knots.

W-620-30 W-620-36 W-620-43	W-621-30 W-621-36 W-621-43	W-622-30 W-622-36 W-622-43	W-623-30 W-623-36 W-623-43	W-624-30 W-624-36 W-624-43
----------------------------------	----------------------------------	----------------------------------	----------------------------------	----------------------------------

Dissecting & Grasping Forceps




→ Operating & Instrumental Information ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length						
30 cm / 12" 36 cm / 14" 43 cm / 17"						



Manhes Micro Grasping Forceps, with multiple teeth. With a very fine uptake at the end of the jaws, the claws of which go no further than the serous membranes.

W-625-30	W-626-30	W-627-30	W-628-30	W-629-30	W-630-30
W-625-36	W-626-36	W-627-36	W-628-36	W-629-36	W-630-36
W-625-43	W-626-43	W-627-43	W-628-43	W-629-43	W-630-43



Manhes Micro Dissecting and Grasping Forceps, atraumatic, full grasping

W-631-30	W-632-30	W-633-30	W-634-30	W-635-30	W-636-30
W-631-36	W-632-36	W-633-36	W-634-36	W-635-36	W-636-36
W-631-43	W-632-43	W-633-43	W-634-43	W-635-43	W-636-43



Manhes Micro Dissecting and Grasping Forceps, atraumatic, short jaws, partial grasping.

W-637-30	W-638-30	W-639-30	W-640-30	W-641-30	W-642-30
W-637-36	W-638-36	W-639-36	W-640-36	W-641-36	W-642-36
W-637-43	W-638-43	W-639-43	W-640-43	W-641-43	W-642-43



Manhes Micro Dissecting and Grasping Forceps, atraumatic, long jaws, partial grasping. Especially suitable for tying intra-corporeal knots.

W-643-30	W-644-30	W-645-30	W-646-30	W-647-30	W-648-30
W-643-36	W-644-36	W-645-36	W-646-36	W-647-36	W-648-36
W-643-43	W-644-43	W-645-43	W-646-43	W-647-43	W-648-43

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 36 cm,
to be use with trocars size 6 mm.

→ **Operating & Instrumental Information** ←
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Metal Handle with outer tube.

	SEMME Grasping Forceps, without connector pin for unipolar coagulation, size 5 mm, length 36 cm	
36 cm / 14"	W-649-36	
	Grasping Forceps, without connector pin for unipolar coagulation, size 5 mm, length 36 cm	
36 cm / 14"	W-650-36	
	SEMME Atraumatic Grasping Forceps and Ampulla Dilator, without connector pin for unipolar coagulation, size 5 mm, length 36 cm	
36 cm / 14"	W-651-36	

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 36 cm.

→ **Operating & Instrumental Information**
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Metal handle, without ratchet.	Metal handle, with Manhes style ratchet.	Metal handle, with hemostat style ratchet.

Insulated outer tube with Luer lock connector for cleaning, size 10mm 	
36 cm / 14"	W-652-36

Metal outer tube with Luer lock connector for cleaning, size 10mm 	
36 cm / 14"	W-653-36

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 36 cm, to be use with trocars size 11 mm.

→ **Operating & Instrumental Information**
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 10 mm.

All are Single / Double Action Jaws.

Working Length			
36 cm / 14"			
	Claw Forceps, 2 x 3 teeth, single action jaws.		
	W-654-36	W-655-36	W-656-36
	Claw Forceps, 2 x 3 teeth, short, single action jaws.		
	W-657-36	W-658-36	W-659-36
	Grasping Forceps, fenestrated, double action jaws.		
	W-660-36	W-661-36	W-662-36
	Dissection and Grasping Forceps, right angled double action jaws.		
	W-663-36	W-664-36	W-665-36

All are Double Action Jaws.

	Sawalhe Tenaculum Forceps.		
	W-666-36	W-667-36	W-668-36
	Sawalhe Tenaculum Forceps, short.		
	W-669-36	W-670-36	W-671-36
	Sawalhe Tenaculum Grasping Forceps.		
	W-672-36	W-673-36	W-674-36

Dissecting & Grasping Forceps



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 36 cm, to be use with trocars size 11 mm.

→ **Operating & Instrumental Information**
Rotating, Dismantling,
without Connector Pin For Unipolar Coagulation,
Size 10 mm.

All are Double action jaws.

Working Length			
36 cm / 14"			
 Kelly Dissecting and Grasping Forceps, long.			
	W-675-36	W-676-36	W-677-36
 Duval Grasping Forceps.			
	W-678-36	W-679-36	W-680-36
 Bowel Grasping Forceps, two rows of atraumatic teeth			
	W-681-36	W-682-36	W-683-36
 Bowel Grasping Forceps.			
	W-684-36	W-685-36	W-686-36
 Babcock Grasping Forceps.			
	W-687-36	W-688-36	W-689-36
 Babcock Grasping Forceps, rounded			
	W-690-36	W-691-36	W-692-36
 Babcock Grasping Forceps, long.			
	W-693-36	W-694-36	W-695-36
 Grasping Forceps, atraumatic for large organs, length of jaws 70 mm			
	W-696-36	W-697-36	W-698-36

Surgical Sponge Holder Dissecting & Grasping Forceps

Operating & Instrumental Information



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, 36 cm, to be use with trocars size 6 mm/11 mm..



Surgical Sponge Holder, for atraumatic dissection of tissue layers, size 10 mm, length 30 cm, 36 cm

30 cm / 12"
36 cm / 14"

W-699-36
W-700-36



Grasping Forceps, rotating, articulating, size 10mm, with connector pin for unipolar coagulation, length 36 cm.

36 cm / 14"

W-701-36

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Dismantling,
coaxial curved up, Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.



All are Single Action Jaws

Working Length	
30 cm / 12" 36 cm / 14"	
	Cuschieri Grasping Forceps , atraumatic, jaws with multiple teeth, fenestrated.
	W-702-30 W-702-36
	Cuschieri Grasping Forceps , atraumatic, jaws with multiple teeth, fenestrated, long.
	W-703-30 W-703-36
	Cuschieri Grasping Forceps , atraumatic, jaws with multiple teeth, width of jaws 4.8 mm, for atraumatic, accurate grasping.
	W-704-30 W-704-36
	Cuschieri Grasping Forceps , atraumatic, jaws with multiple teeth, width of jaws 3 mm, for atraumatic, accurate grasping.
	W-705-30 W-705-36

Dissecting & Grasping Forceps




→ Operating & Instrumental Information ←
Dismantling,
coaxial curved up, Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 9 mm.



All are Single / Double Action Jaws

Working Length	
30 cm / 12" 36 cm / 14"	
	Cuschieri Grasping Forceps , atraumatic, jaws with multiple teeth, fenestrated.
	W-706-30 W-706-36
	Cuschieri Grasping Forceps , atraumatic, jaws with multiple teeth, fenestrated, long.
	W-707-30 W-707-36
	Cuschieri Grasping Forceps , atraumatic, jaws with multiple teeth, width of jaws 4.8 mm, for atraumatic, accurate grasping.
	W-708-30 W-708-36
	Cuschieri Grasping Forceps , atraumatic, jaws with multiple teeth, width of jaws 3 mm, for atraumatic, accurate grasping.
	W-709-30 W-709-36
	Cuschieri Dissecting and Grasping Forceps.
	W-710-30 W-710-36

Dissecting & Grasping Forceps



→ Operating & Instrumental Information ←
Dismantling,
coaxial curved up, Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.



All are Single / Double Action Jaws

Working Length	
30 cm / 12" 36 cm / 14"	
	Cuschieri Grasping Forceps , atraumatic with multiple teeth, fenestrated, single action jaws.
	W-711-30 W-711-36
	Cuschieri Grasping Forceps , with multiple teeth, width of jaws 4.8 mm, for atraumatic and accurate gripping, single action jaws.
	W-712-30 W-712-36
	Cuschieri Grasping Forceps , with multiple teeth, width of jaws 3 mm, for atraumatic and accurate gripping, single action jaws.
	W-713-30 W-713-36
	Cuschieri Dissection and Grasping Forceps , double action jaws.
	W-714-30 W-714-36

Scissors, Punch Forceps & Biopsy Forceps



Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 20 cm,
to be use with trocars size 3.5 mm.

→ **Operating & Instrumental Information**
Rotating, Dismantling,
Insulated with Connector Pin For Unipolar Coagulation,
Size 5 mm.

All are Single / Double Action Jaws

Working Length	Dia.	
20 cm / 8"	3mm	

	Scissors , with serrated jaws, curved, conical, double action jaws.
W-715-20	

	Micro Hook Scissors , single action jaws.
W-716-20	

	Biopsy Punch Forceps , through cutting, single action jaws.
W-717-20	

	Blakesley Dissecting and Biopsy Forceps , single action jaws.
W-718-20	

Scissors



Multiple Puncture Approach.
 Operating instruments for Multiple
 Puncture approach length 30 cm,
 to be use with trocars size 3.5 mm.

→ **Operating & Instrumental Information**
 Rotating, Dismantling,
 without Connector Pin For Unipolar Coagulation,
 Size 5 mm.

All are Single Action Jaws

Working Length	Dia.	
30 cm / 12"	3mm	

	Koh Ultramicro Suture Scissors, sharp.
---	--

W-719-30	
----------	--

	Koh Ultramicro Suture Scissors, blunt.
---	--

W-720-30	
----------	--

Scissors, Punch Forceps & Biopsy Forceps



Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
to be use with trocars size 3.5 mm.

→ **Operating & Instrumental Information**
Rotating, Dismantling,
with Connector Pin For Unipolar Coagulation,
Size 3.5 mm.

All are Single / Double Action Jaws

Working Length	Dia.	
30 cm / 12"	3.5mm	

	Scissors , with serrated jaws, curved, conical, double action jaws.
W-721-30	

	Micro Hook Scissors , single action jaws.
W-722-30	

	Biopsy Punch Forceps , through cutting, single action jaws.
W-723-30	

	Blakesley Dissecting and Biopsy Forceps , single action jaws.
W-724-30	

Scissors, Cutting-Coagulating Forceps



Operating & Instrumental Information

Rotating, Dismantling,
with & without Connector Pin For Unipolar Coagulation,
Size 5 mm.


Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single / Double Action Jaws

Working Length			
30 cm / 12" 36 cm / 14" 43 cm / 17"			
	Manhes Scissors , serrated, single action jaws.		
	W-725-30 W-725-36 W-725-43	W-726-30 W-726-36 W-726-43	W-727-30 W-727-36 W-727-43
	Hook Scissors , single action jaws.		
	W-728-30 W-728-36 W-728-43	W-729-30 W-729-36 W-729-43	W-730-30 W-730-36 W-730-43
	Hook Scissors , jaws not crossing, single action jaws.		
	W-731-30 W-731-36 W-731-43	W-732-30 W-732-36 W-732-43	W-733-30 W-733-36 W-733-43
	Shawki Cutting-Coagulation Forceps , double action jaws.		
	W-734-30 W-734-36 W-734-43	W-735-30 W-735-36 W-735-43	

Dissecting & Grasping Forceps



Scissors

→ Operating & Instrumental Information ←
Rotating, Dismantling,
with & without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

All are Single Action Jaws

Working Length			
30 cm / 12" 36 cm / 14"			
	Micro Dissecting Scissors.		
	W-736-30 W-736-36	W-737-30 W-737-36	W-738-30 W-738-36
	Micro Dissecting Scissors, curved.		
	W-739-30 W-739-36	W-740-30 W-740-36	W-741-30 W-741-36
	Micro Hook Scissors.		
	W-742-30 W-742-36	W-743-30 W-743-36	W-744-30 W-744-36

Scissors



→ **Operating & Instrumental Information** ←
Rotating, Dismantling,
with & without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Double Action Jaws

Working Length			
30 cm / 12" 36 cm / 14" 43 cm / 17"			
	Metzenbaum Scissors , curved, length of blades 12 mm.		
	W-745-30 W-745-36 W-745-43	W-746-30 W-746-36 W-746-43	W-747-30 W-747-36 W-747-43
	Scissors , with serrated jaws, curved, spoon blades, length of blades 17 mm.		
	W-748-30 W-748-36 W-748-43	W-749-30 W-749-36 W-749-43	W-750-30 W-750-36 W-750-43
	Scissors , with serrated jaws, curved, conical.		
	W-751-30 W-751-36 W-751-43	W-752-30 W-752-36 W-752-43	W-753-30 W-753-36 W-753-43
	Scissors , straight.		
	W-754-30 W-754-36 W-754-43	W-755-30 W-755-36 W-755-43	W-755-30 W-755-36 W-755-43

Scissors



→ **Operating & Instrumental Information** ←
Rotating, Dismantling,
with & without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length			
30 cm / 12" 36 cm / 14" 43 cm / 17"			
	Manhes Micro Scissors , short jaws with rounded tips, self-sharpening.		
	W-756-30 W-756-36 W-756-43	W-757-30 W-757-36 W-757-43	W-758-30 W-758-36 W-758-43
	Manhes Micro Scissors , long jaws with rounded tips, self-sharpening.		
	W-759-30 W-759-36 W-759-43	W-760-30 W-760-36 W-760-43	W-761-30 W-761-36 W-761-43

Biopsy Forceps & Punch Forceps



→ Operating & Instrumental Information ←
Rotating, Dismantling,
with & without Connector Pin For Unipolar Coagulation,
Size 5 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 6 mm.

Single Puncture Approach.
Operating instruments for single Puncture
approach length, 43 cm, to be use with
operating laparoscopes with built-in
instrument channel.

All are Single Action Jaws

Working Length			
30 cm / 12" 36 cm / 14" 43 cm / 17"			
Blakesley Dissection and Biospy Forceps.			
	W-762-30 W-762-36 W-762-43	W-763-30 W-763-36 W-763-43	W-764-30 W-764-36 W-764-43
Manhes Biopsy Forceps.			
	W-765-30 W-765-36 W-765-43	W-766-30 W-766-36 W-766-43	W-767-30 W-767-36 W-767-43
Frangenheim Biopsy Punch Forceps, through-cutting.			
	W-768-30 W-768-36 W-768-43	W-769-30 W-769-36 W-769-43	W-770-30 W-770-36 W-770-43
Biopsy Punch Forceps, with two teeth			
	W-771-30 W-771-36 W-771-43	W-772-30 W-772-36 W-772-43	W-773-30 W-773-36 W-773-43



Handle & Outertube

→ **Operating & Instrumental Information**
Rotating, Dismantling,
for Scissors and Biopsy Forceps.



Metal handle, with Manhes style ratchet.
Size 10 mm.

Insulated outer tube with Luer lock connector for cleaning, size 10mm



36 cm / 14"

W-774-36

Metal outer tube with Luer lock connector for cleaning, size 10mm



36 cm / 14"

W-775-36

Scissors & Spoon Forceps



→ **Operating & Instrumental Information** ←
Rotating, Dismantling,
Size 10 mm.

Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 30 cm,
36 cm, to be use with trocars size 11 mm.

All are Single Action Jaws


Working Length	
30 cm / 12"	
	Scissors , heavy blades, with cylindrical grinding.
	W-776-30
	Spoon Forceps.
	W-777-30

Articulating Scissors & Biopsy Forceps



Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 30 cm, 36 cm, to be use with trocars size 6 mm/11 mm.

Operating & Instrumental Information



36 cm / 14"

W-778-36



36 cm / 14"

W-779-36

Cuschieri Scissors



→ Operating & Instrumental Information ←

Dismantling,
with and without connector pin for unipolar coagulation,
size 5 mm

Multiple Puncture Approach.

Operating instruments for Multiple Puncture
approach length 30 cm, 36 cm, to be use
with trocars with flexible sheath, size 11 mm.

All are Single Action Jaws

Working Length		
30 cm / 12" 36 cm / 14"		
	Cuschieri Scissors, serrated, bayonet-shaped.	
W-780-30 W-780-36		
	Cuschieri Scissors, serrated, coaxial curved up.	
W-781-30 W-781-36		
	Cuschieri Scissors, serrated, coaxial curved down.	
W-782-30 W-782-36		



Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach to be use
with trocars size 7 mm.

Ring Applicators

Operating & Instrumental Information



Ring Applicator.
May be fitted with 2 silastic rings.

W-783-01



Ring Applicator.
May be fitted with 2 silastic rings.


W-784-02



Ring Applicator.
May be fitted with 2 silastic rings.

W-785-03

BERCI Fascial Closure Instruments





Instruments Set for Endoscopic Axillary Lymphadenectomy

Operating & Instrumental Information

Multiple Puncture Approach.
Operating instruments for Multiple Puncture approach length 25 cm, to be use with trocars size 6 mm.

Working Length	Dia.		
25 cm / 9 3/4"	5mm		
		 Manhes Grasping Forceps , with multiple teeth, width of jaws 4.8 mm, single action jaws.	
			W-787-25
		 Metzenbaum Scissors , curved, double action jaws.	
		W-788-25	



10 mm / 36 cm

Operating & Instrumental Information

Multiple Puncture Approach.
Operating instruments with out connector pin for unipolar coagulation, for Multiple Puncture approach length 36 cm, to be use with trocars size 3.9 mm.


Working Length	Dia.		
36 cm / 14"	10mm		
		 Babcock Clamp , fenestrated, double action jaws.	
		W-789-36	



Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 36 cm,
to be use with trocars size 6 mm.

Intraoperative Cholangiography


Operating & Instrumental Information
Instruments for catheter placement,
size 5 mm

	
36 cm / 14"	W-790-36
	
36 cm / 14"	W-791-36
	
36 cm / 14"	W-792-36



Spilling Stone Recovery Forceps


Operating & Instrumental Information
for fast retrieval of spilled gallstones
during cholecystectomies.





Spilling Stone Recovery Forceps

Operating & Instrumental Information
for fast retrieval of spilled gallstones
during cholecystectomies.


Working Length	
20 cm / 8"	
	Reddick-Olsen Dissection and Grasping Forceps, S-shaped, length 20 cm.
W-794-20	
	Dissection and Grasping Forceps, S-shaped, length 20 cm.
W-795-20	
	Dissection and Grasping Forceps, S-shaped, length 20 cm.
W-796-20	
	Metzenbaum Scissors, S-shaped, length 20 cm.
W-797-20	

Single Clip Applicators for Titanium Clips



Multiple Puncture Approach.
Operating instruments for Multiple
Puncture approach length 33 cm,
to be use with trocars size 11 mm.

→ **Operating & Instrumental Information**
Dismantling, rotating,
size 10mm.



Applicator, dismantling, rotating, for pilling ligation
clips, (medium-large) with ratched to lock the jaw part holding the clip.

33 cm / 13"

W-798-33



Applicator, dismantling, rotating, for pilling ligation
clips, (medium) with ratched to lock the jaw part holding the clip.


33 cm / 13"

W-799-33

The C.C.L. Vaginal Extractor



→ **Operating & Instrumental Information**
a new extractor for
laparoscopic surgery.



Rotating Claw Forceps,
single action jaws, size 10 mm, length 36 cm.

36 cm / 14"

W-800-36



Vaginal Extractor, 11 mm
with insulated ball-sphere.


10 cm / 4"

W-801-10

The C.C.L. Vaginal Extractor



→ **Operating & Instrumental Information**
a new extractor for
laparoscopic surgery.



Semm Tissue Punch, size 10 mm, with built-in system for
multiple tissue specimen collection.

28 cm / 11"

W-802-28




Solid Tissue Removal Forceps.

18 cm / 7"

W-803-18


Gall Bladder & Endobag Extractors



Blade Speculum
for Gall Bladder Extraction
W-804-18
18cm / 7"



Speculum
for Gall Bladder Extraction
W-805-19
19cm / 7 1/2"



Koechli Endobag extractor, 3 valves,
with curved valves in two different lengths,
interchangeable, continuously adjustable.
42-806-19
19cm / 7 1/2"



W-807-19
Valve, Short, for normal
abdominal walls, 3 piece-set



W-808-19
Valve, long, for obese
abdominal walls, 3 piece-set

Gall Bladder & Endobag Extractors



Koechli Endobag Extractor
W-809-17
17cm / 6 3/4"



Fig. 1



Koechli Endobag Extractor
W-810-17
17cm 6 3/4"





Fig. 2

Gasless Laparoscopy Instruments




Forcep 1 x 2 teeth
W-811-26
26cm 10 1/4"



Metzenbaum
Scissors Curved
W-812-18
18cm 7"




Needle Holder
with ratchet
W-813-26
26cm 10 1/4"



Myoma Fixation Instrument
W-814-30
30cm 12"




Koeckerling Knot Tier,
Size 5 mm
W-815-36
36cm 14.5"



Atraumatic Clamp
W-816-26
26cm 10 1/4"



Clamp
Strong, Curved
W-817-26
26cm 10 1/4"



Atraumatic Clamp
Curved
W-818-26
26cm 10 1/4"



Tenaculum Forceps
W-819-26
26cm 10 1/4"



Claw Forceps
W-820-26
26cm 10 1/4"



WELMED
INSTRUMENTS CO.

CONTACT US

MANUFACTURERS & EXPORTERS OF SURGICAL, DENTAL, ORTHODONTIC,
DIAGNOSTIC ,PLASTIC SURGERY & ALL INSTRUMENTS.

Address: Defence Road, Near Wapda Grid Station,
Sialkot 51310 - Pakistan.

Cell : +92-333 8680640 / 340 4040061

Phone : +92 52 3554782 / 52 3558704 / 333 8680640

Email: welmedpk@gmail.com / welmedco@gmail.com
info@welmedinst.com / welmed@cyber.net.pk

Web: welmedinst.com

