



WELMED
INSTRUMENTS CO.
Quality Surgical Instruments

Manufacturers & Exporters of Surgical, Dental & Veterinary Instruments



DIAGNOSTIC

www.welmedinst.com



MINI OTOSCOPE FIBER OPTICS



WI-30-101



WI-30-102



WI-30-103



Available in colour





MINI OTOSCOPE FIBER OPTICS



WI-30-104



WI-30-105



WI-30-106



WI-30-107



Available in colour





MINI OTOSCOPE FIBER OPTICS & OPHTHALMOSCOPE SET



WI-30-108



WI-30-109



WI-30-110



Available in colour





METAL BODY OTOSCOES



WI-30-111



WI-30-112



WI-30-113




WI-30-114





METAL BODY OTOSCOPES



WI-30-115




WI-30-116



WI-30-117



WI-30-118





METAL BODY & PLASTIC OPHTHALMOSCOPES



WI-30-119



Also Available in Set



WI-30-120



Also Available in Set



OPHTHALMOSCOPES



WI-30-121



WI-30-122



WI-30-123



WI-30-124

Available in colour





TONGUE DEPRESSOR



WI-30-125



PEN TORCH



WI-30-126



WI-30-127

Available in colour





ENT SETS



WI-30-128



BULBS / LAMPS



Set



DERMATOSCOPE SET



Features:

All metal construction for durability
Bayonet locking system for easy engagement
and disengagement from the handle
High resolution lenses system with magnification

- 1- Zipper case
- 2- Bulb for dermatoscope



WI-30-129



MACINTOSH CONVENTIONAL



Fig,4



Fig,3



Fig,2



Fig,1



Handle

Advantages

Handles

Brass
Universal Fitting with
non magnetic pin

Blades

Stainless Steel
Non magnetic
with magnetic pin



WI-30-135



MILLER CONVENTIONAL



Fig,4



Fig,3



Fig,2



Fig,1



WI-30-131

Miller English Set of 4
Blade Size: 1,2,3,4 1 Medium Handle.
1 small & 1 Large spare Lamp Packed in
Strong Plastic Casing with Foam Inserting.



Large Lamp



Small




MILLER MEGALIGHT

Miller Fix tube HBT (Fiber Optic)
Cold Light Guide

FEATURES

303/304 Surgical stainless steel
Smooth edges, high quality surface finish
Provides with 3mm Fiber Optic offers no Black Marks
fully Autoclavable and easier to clean satin finish.





MILLER FIBER OPTIC



Fig,4



Fig,3



Fig,2



Fig,1



Handle

Advantages

Economical optic Fiber blades
Fiber optic illumination assembly is detachable and is easy to clean, maintain and re-polished. Thousands of fiber strings are firmly secured in stain less steel pipe cover with saves from damage or breakage.



WI-30-133



FIBER OPTIC LARYNGOSCOPES

MACINTOSH FIBER OPTIC BLADES



Fig, 4



Fig, 3



Fig, 2



Fig, 1




WI-30-134




MACINTOSH FIBER OPTIC



Fig,4



Fig,3



Fig,2



Fig,1



Handle




Advantages

Economical optic blades.
Fiber optic illumination assembly is detachable and is easy to clean, maintain and re-polished. Thousands of fiber strings are firmly secured in stain less steel pipe cover with saves from damage or breakage.

WI-30-136



SINGLE USE / CONVENTIONAL




WI-30-137

Advantages

Low Cost: Magnetic Stainless steel
Blades with plastic connector
and Fiber plastic pipe.



WI-30-138



Advantages

Very Low Cost: Magnetic Stainless steel
Blades with Screw type fitting
and Fiber plastic pipe.



HANDLES

FIBER OPTIC



CONVENTIONAL



Brass Handles
Highy quality
Long life
Universal fitting with non
magnetic pin



WELMED
INSTRUMENTS CO.

Quality Surgical Instruments

Defence Road, Near Wapda Grid Colony
Sialkot, 51310 - Pakistan.

Tel : +92-52-3554782 - 3558704

Fax: +92-52-3550296

Cell: +92-333-8680640

E-mail: welmed@cyber.net.pk

: wmm@brain.net.pk

: info@welmedinst.com

Web : www.welmedinst.com

ISO 13485 / 9001
CE Mark