



WELMED

INSTRUMENTS CO.

Serving Humanity Since 1908


MANUFACTURER OF ORTHODONTIC INSTRUMENTS

ORTHODONTICS



www.welmedinst.com

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-186
Castroviejo



DI-17-187
Bone Caliper Double Scale
15.5 cm



DI-17-190
Boley



DI-17-188
Separating Plier for Separators



DI-17-189
Ligature Tying Plier



DI-17-191
Beerendonk



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-101 TC
Wire Cutter upto
022" x .028" Ligature Wires
(0.22 mm x 0.51 mm)



DI-17-102 TC
Hard Wire Cutter
022" x .028"
(0.56 mm x 0.71 mm)



DI-17-103 TC
Side Nipper for Lock Pins and
Ligatures for Soft Thin Wire
Ø .016" (0.4 mm)



DI-17-104 TC
Wire Cutter Angled 45°
upto .015" (0.38 mm)



DI-17-105 TC
Micro Light Wire Cutter upto
Max 0.15"



DI-17-106 TC
Micro Distal & Cutter
Safety Hold Hard Wire
.022" x .025"
(0.55 mm x 0.64 mm)



DI-17-107 TC
Wire Cutter upto
.014" (0.36 mm)

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-108 TC

Distal End Cutter Safety Hold

Re : .022" x .028" (0.55 mm x 0.71 mm)


Ro : .020" x .016" (0.51 mm x 0.41 mm)

DI-17-109 TC

Micro Distal End Cutter - Safety Hold

Re : upto .022" x .028" (0.55 mm x 0.71 mm)

Ro : .020" (0.50 mm)



DI-17-110 TC

Micro Distal End Cutter
 Safety Hold

upto 0.21 mm x 0.70 mm
 14 cm

DI-17-111 TC

Distal End Cutter

Safety Holding (Round Shape)

Re : 0.56 mm x 0.75 mm

Ro : 0.51 mm

DI-17-112

Posterior Band Remover

DI-17-113

Direct Bonding Bracket
 Remover with Pad

DI-17-114 TC

Detailing Step Plier Special

Double Side Offset Head

Forms Bayonet Bands.

Forms Right Hand or Left Hand

Step-up or Step-Down

Bends In 1 mm Increments.

Remover with Pad



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-176
Beebee
(Cur)

DI-17-175
Beebee
(Dir)




DI-17-177
Universal



DI-17-178
Mosquito Forceps
(Dir)

DI-17-179
Mosquito Forceps
(Cur)

DI-17-180 A .018"
DI-17-180 B .022"



DI-17-180 A
Aluminum Height Gauge



DI-17-181
Mershon Band
Pusher



DI-17-182
Bend Setter with
Tin Insert 15 cm



DI-17-183
Bend Pusher
Plastic Handle



DI-17-184
Mathieu Without TC
14 cm



DI-17-185
Mathieu
14 cm




Mathieu With TC
14 cm



DI-17-184T
THIN



DI-17-184K
KOCHER



DI-17-184B
HEAVY



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri




DI-17-115
Direct Bonding Attachment
Bracket Remover (Angled) with Stop Screw



DI-17-116
Direct Bonding Bracket
Bracket Remover Narrow, with Stop Screw



DI-17-117
TC and without TC
Tweed Angle Arch Banding
Blade 14 mm



DI-17-118
TC and without TC
Tweed Angle Arch Banding
Blade 10 mm



DI-17-119
Bend Seating Plier Will Firmly
Grip The Bracket's Tie Wing Aiding
In The Seating of The Band. Designed
with Serrated Tips For Maximum Grip
That Prevent Slippage During use.



DI-17-120
TC and without TC
Tweed Angle Arch Banding
Blade 6.5 mm




DI-17-121
TC and without TC
Arch Forming Plier



DI-17-122
Legature Arch Bending Plier

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-160 **DI-17-161** **DI-17-162** **DI-17-163** **DI-17-164**



DI-17-165
 Distal Bender Simple to use
 Bends tails of springs Works on all Stainless
 Steel and most NITI and Beta III Wires



DI-17-166
 NITI Ice Instruments Turns NITI wire soft and
 flexible Immediately Designed for Intraoral use
 Durable Stainless Steel design



DI-17-167
 Positive Positioner Precise Placement
 of any direct bond attachment
 Light weight design Hard Anodized aluminum



DI-17-169
 Bracket Positioning Gauge



DI-17-168
 Taylor Tweezer 15 cm



DI-17-170 Twirl-On Director Applicator



DI-17-171 Legatur Director



DI-17-172 Legatur Director



DI-17-173 Band Pusher and Schure Scaler



DI-17-174 Band Pusher and Schure Scaler

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-123
TC and without TC
Utility Plier Slim Tip



DI-17-124
TC and without TC
Bend Seating Plier (Dir)



DI-17-125
TC and without TC
How Plier with a 45° Angle
of Tips (Cur)



DI-17-126
TC and without TC
Weingart Well Grip Plier
with Serrated



DI-17-127
Tweed Plier Loop Forming
with Stop Screw
Re : 0.56 mm x 0.71 mm



DI-17-128
Tweed Plier Loop Forming
with Stop Screw
Re : 0.56 mm x 0.71 mm



DI-17-129
Angled Bell Hook Crimping Plier



DI-17-130
Bell Hook Crimping Plier



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri

DI-17-151
Universal

DI-17-152
New Universal

DI-17-153
Nance 15 cm

DI-17-154
Schwarz

DI-17-155
Waldschs 15 cm

DI-17-156
Heavy Wire Cutter
for Hard Wires upto $\varnothing$ 1.5 mm

DI-17-157
Niti Tip Back Plier

DI-17-158
Aderer 14 cm

DI-17-159
Versatile Light Wire Bending Pier.
Grooves Assist with Precise Bending
and Closing of Loops. Serrated Surface
Grip Wire upto .020" (51 mm)

Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-141
 Palatal Arch Plier
 Allow Easy Placement or Removal



DI-17-142
 Flat Nose Plier with Smooth
 Breaks for Bending



DI-17-143
 Adams Plier for Wire Bending




DI-17-144
 Round Nose Plier
 for Wire Bending



DI-17-145
 Double Rounded Jaw
 Wire Bending Plier



DI-17-146
 Mckellops
 16 cm



DI-17-147
 Mckellops
 14 cm



DI-17-148
 Reynolds



DI-17-149
 Bird Beak Plier Bending Round,
 Square or Rectangular




DI-17-150
 Bird Beak Plier Bending Round,
 Square or Rectangular



Orthodontic Instruments
ортодонтических инструментов
Instrumentos ortodônticos
Ortodonti Aletleri



DI-17-131
ADERER
Three Jaw Wire Bending Plier
(Thin Large Points)



DI-17-132
ADERER
Three Jaw Wire Bending Plier
(Medium Points)



DI-17-133
YOUNG
Loop Bending Plier




DI-17-134
JOHNSON
Contouring Plier



DI-17-135
Nance Closing Loop Plier



DI-17-136
ADERER
Three Jaw Wire Bending Plier



DI-17-137
Light Wire Pliers with Cutter Beatz



DI-17-138
Clasp Bending Plier



DI-17-139
Arch Forming Plier



DI-17-140
Arch Forming Plier (Plain)



WELMED

INSTRUMENTS CO.

Serving Humanity Since 1908

Defence Road, Near Wapda Grid Colony
Sialkot, 51310 - Pakistan.

Tel : +92-52-3554782 - 3558704

Fax: +92-52-3550296

Cell: +92-333-8680640 

E-mail: welmedpk@gmail.com

: wmm@brain.net.pk

: welmed@cyber.net.pk

Website: info@welmedinst.com

ISO 13485 / 9001
Mark
FDA USA (certified)